



COVID-19 related ARDS – Non-invasive respiratory support

DL-45340-2021

The development of clinical data on COVID-19 has been very rapid, resulting in a large amount of data being generated in a very short time. However, hard evidence still appears to be scarce. This article has been written to the best of our knowledge based on selected literature and opinions of clinical experts. It does not represent a summary of all available literature and therefore does not claim to be exhaustive. As COVID-19 is a very complex disease, you should always refer to the original literature mentioned in this article, other relevant literature and the circumstances of the individual case when deciding on the right ventilation strategy for your patients. It is also strongly advised to follow your national/local guidelines and standards.

There seem to be high hopes that non-invasive respiratory support would prevent the need for endotracheal intubation (ETI) and invasive ventilation. However, there is also the fear of missing the right point for ETI in the individual patient. There is not much evidence specifically for COVID-19 patients, therefore recommendations are not available on all aspects on non-invasive respiratory support and

the question of when to decide for ETI. In the following we will take a practical approach: We compiled the recommendations of 4 different guidelines and complemented this with information from recent publications. This is not to be seen as a review reflecting all currently available data, but more as a practical approach to get at least some guidance.

Reviewed Guidelines:

- German S3 Guideline – Recommendations for the therapy of hospitalised patients with COVID-19, Version 4.1, February 2021 [referred to as **GS3**]³⁴
- Surviving Sepsis Guidelines on the Management of Adults with Coronavirus Disease 2019 (COVID-19) in the ICU: First Update; March 2021 [referred to as **SSC**]³⁵
- ERS Guideline for the Management of hospitalized adults with coronavirus disease 2019 (COVID-19): A European Respiratory Society living guideline, January 2021 [referred to as **ERS**]³⁶
- Australian guideline for clinical care of people with COVID-19, National COVID-19 Clinical Evidence Taskforce [referred to as **NCCET**]³⁷

Recommendations for non-invasive respiratory support

Conventional Oxygenation:

The SSC Guideline is the only reviewed guideline to specifically comment on conventional oxygen therapy: The guideline suggests to start supplemental oxygen if peripheral SaO₂ is

<92% (strong), but recommends it at an SaO₂ <90% (strong). The guideline recommends to maintain SaO₂ no higher than 96% in patients with acute hypoxemic respiratory failure (AHRF).



HFNC AND CPAP/NIV:

Treatment with High Flow Nasal Cannula (HFNC) or CPAP/NIV is either suggested or (conditionally) recommended by all guidelines, in part with specific additional information.

There is a uniform recommendation not to delay intubation.

Close monitoring important. All guidelines call for close monitoring of patients for improvement or further deterioration to avoid the delay of intubation, as this may be associated with poor outcomes.

GS3:

- Recommendation to carry out a trial of HighFlow or CPAP/ NIV for patients with hypoxemic resp. failure (P/F=100-300) under permanent monitoring and intubation readiness. Few direct comparisons of HFNC, CPAP and NIV are currently available. Therefore, evidence for non-COVID Patients needs to be evaluated.

SSC:

- Suggests using HFNC over NIPPV (weak); if HFNC is not available and absence of urgent indication for ETI, the guideline suggests NIPPV trial with close monitoring and short-interval assessment for worsening.
- If HFNC is not available and absence of urgent indication for ETI, the guideline suggests NIPPV trial with close monitoring and short-interval assessment for worsening.

ERS:

- Conditionally recommends that hospitalized patients with COVID-19 should receive HFNC or CPAP. The authors urge to not delay intubation if patients fail to respond to this measure.

NCCET:

- Conditional recommendation for HFNC to be considered in patients unable to maintain SaO₂ of ≥92% despite conventional oxygen delivery at >6l/min or an FiO₂ = 0,4.
- Consider using NIV therapy for patients with hypoxemia associated with COVID-19.
- Generally, the guideline recommends strict attention to staff safety (unconditional use of PPE) and the use of single rooms or negative pressure rooms wherever possible.
- Contact, droplet and airborne precautions are to be in place.
- The guidelines recommends against the usage of HFNC and NIV in COVID-19 patients in shared wards, emergency department cubicles or during interhospital patient transfer/ retrieval.

i

HFNC OR CPAP/NIV AND THE QUESTION ON WHAT WORKS BETTER FOR NIV: HELMET OR MASK

HFNC OR CPAP/NIV is a matter of uncertainty. There is one recent larger study in which 2 approaches were tested on 100 COVID-19 patients³⁸:

1. Helmet for 48h followed by HFNC
2. HFNC only.

RESULT: There was no difference with respect to primary outcome (respiratory support-free days). However, significantly fewer patients with helmet + highflow were intubated. A large NHS study, the RECOVERY trial, is ongoing, randomizing COVID-19 patients with acute respiratory failure to receive HFNC, CPAP, NIV or conventional oxygen therapy.³⁹ Primary outcome parameters are intubation rate and 30-day mortality. At the time of writing, already more than 1,000 patients were recruited, and the results are urgently awaited.³⁸

CPAP/NIV – helmet or mask. Not much evidence is available with regards to COVID-19, data from non-COVID trials therefore need to be taken into account. One example is an RCT from 2016 comparing NIV with helmet and mask.⁴⁰ This trial demonstrated that the helmet may be superior to the mask: More than 60% of patients receiving NIV by mask required secondary intubation, quite a high NIV failure rate; only 18% of the helmet group required intubation. This led to a significantly reduced 90-day mortality. As reason for the good result favoring the helmet, the authors discussed the applicability of a higher PEEP in helmets.³⁸

RELEVANT IN DAILY CLINICAL PRACTICE? But is the information of which non-invasive way of support is chosen really that clinically relevant? Maybe not. The choice depends on the acceptance of the individual patient. Not every method is equally suited for all patients. Furthermore, in daily clinical practice the mentioned methods are frequently used in parallel: HFNC during meals, helmet during proning, and mask-CPAP in a rehab chair. Future research needs to deliver evidence on any significant advantage of one over the other and for or against combinations.³⁸

(AWAKE) PRONE POSITIONING is not explicitly recommended in all guidelines. The Australian guideline recommends considering prone positioning for at least 3 hours per day if tolerated. The ERS conditionally recommends delivering non-invasive respiratory support with or without adjunct strategies, such as prone positioning, if no immediate indications for invasive mechanical ventilation are present. Furthermore, the ERS states that prone positioning of non-intubated patients with acute hypoxemic respiratory failure due to COVID-19 pneumonia has been recently tested across different settings, wards, or in ICUs as an adjunct to conventional oxygen therapies. Large heterogeneity across these experiences can be recognised. They differ in terms of patient selection, type of oxygen therapy support used, setting, timing and duration of the intervention, and therefore provide variable results. Despite

this heterogeneity, reports document a significant improvement in oxygenation and respiratory rate upon prone positioning, and the majority were able to tolerate the procedure.

In a recent publication, Leasa et al. declare that after initiation of HFNO or CPAP, patients may be encouraged to assume the prone position, particularly if the $\text{PaO}_2/\text{FiO}_2$ ratio is below 200.⁴¹



REMEMBER: Awake prone positioning requires the maintenance of close monitoring and the presence of well-skilled health care staff.



IMPORTANT PRECAUTIONS

Precautions. All reviewed guidelines call for the strict use of appropriate PPE to protect clinical staff, while some provide recommendations for the environment, in which non-invasive respiratory support is provided.

Consider available precautions to limit aerosol spread at patient side, such as surgical masks over HFNC interface. In case of NIV, independent of patient interface (mask or helmet), check if your ventilator offers an Anti-Air Shower functionality. In case of hose disconnection from the patient interface, the ventilator detects the disconnection and reduces the device-given high flow (e.g. for leakage compensation) to a minimum level. This may help to reduce a possible spread of droplets and aerosols in the room exposing clinical staff to a higher risk of infection. Accordingly, reconnection is detected, and ventilation will continue automatically with the previous settings.



In our article on ventilating patients with COVID-19-associated ARDS, we reviewed relevant literature and four current guidelines to provide a practical overview. For references and details, please visit our website: www.draeger.com/covid-ventilation



REFERENCES

1. Williams GW et al, *Anesthesiology* 2021;134:270; Pfortmueller CA et al, *Best Practice & Research Clinical Anaesthesiology* 2020; <https://doi.org/10.1016/j.bpa.2020.12.011>
2. Gattinoni L et al, COVID-19 Does Not Lead to a “Typical” Acute Respiratory Distress Syndrome; *Am J Respir Crit Care Med.* 2020 May 15;201(10):1299-1300. doi: 10.1164/rccm.202003-0817LE.
3. Grasselli G et al., Pathophysiology of COVID-19-associated acute respiratory distress syndrome: a multicentre prospective observational study, *Lancet Respir Med.* 2020 Dec;8(12):1201-1208. doi: 10.1016/S2213-2600(20)30370-2. Epub 2020 Aug 27.
4. Chiumello D et al. Physiological and quantitative CT-scan characterization of COVID-19 and typical ARDS: a matched cohort study. *Intensive care medicine.* 2020;1-10. doi:10.1007/s00134-020-06281-2
5. Ferrando C. et al., Clinical features, ventilatory management, and outcome of ARDS caused by COVID-19 are similar to other causes of ARDS, *Intensive Care Med.* 2020 Dec;46(12):2200-2211. doi: 10.1007/s00134-020-06192-2. Epub 2020 Jul 29.
6. Cummings M. J. et al., Epidemiology, clinical course, and outcomes of critically ill adults with COVID-19 in New York City: a prospective cohort study, *Lancet.* 2020 Jun 6;395(10239):1763-1770. doi: 10.1016/S0140-6736(20)31189-2. Epub 2020 May 19.
7. Botta M. et al., Ventilation management and clinical outcomes in invasively ventilated patients with COVID-19 (PROVENT-COVID): a national, multicentre, observational cohort study, *Lancet Respir Med.* 2021 Feb;9(2):139-148. doi: 10.1016/S2213-2600(20)30459-8. Epub 2020 Oct 23.
8. Gattinoni L. et al., COVID-19 pneumonia: different respiratory treatments for different phenotypes?, *Intensive Care Med.* 2020 Jun;46(6):1099-1102. doi: 10.1007/s00134-020-06033-2. Epub 2020 Apr 14.
9. Marini JJ et al., Integrating the evidence: confronting the COVID-19 elephant, *Intensive Care Med.* 2020 Oct;46(10):1904-1907. doi: 10.1007/s00134-020-06195-z. Epub 2020 Jul 25
10. Gattinoni L et al., The Respiratory Drive: An Overlooked Tile of COVID-19 Pathophysiology, *Am J Respir Crit Care Med.* 2020 Oct 15;202(8):1079-1080. doi: 10.1164/rccm.202008-3142ED
11. Younes Zaid et al, Platelet reactivity to thrombin differs between patients with COVID-19 and those with ARDS unrelated to COVID-19, Published online 2021 Jan 27. , doi: 10.1182/bloodadvances.2020003513
12. Chalmers JD et al, *Eur Resp J*, Management of hospitalised adults with coronavirus disease 2019 (COVID-19): a European Respiratory Society living guideline, *Apr 15;57(4):2100048.* doi: 10.1183/13993003.00048-2021. Print 2021 Apr.
13. Tonelli R et al., Spontaneous Breathing and Evolving Phenotypes of Lung Damage in Patients with COVID-19: Review of Current Evidence and Forecast of a New Scenario, *J Clin Med.* 2021 Mar; 10(5): 975. Published online 2021 Mar 2. doi: 10.3390/jcm10050975
14. Osuchowski MF et al., The COVID-19 puzzle: deciphering pathophysiology and phenotypes of a new disease entity, *Lancet Respir Med.* 2021 May 6;S2213-2600(21)00218-6. doi: 10.1016/S2213-2600(21)00218-6. Online ahead of print
15. Esnault P et al, High Respiratory Drive and Excessive Respiratory Efforts Predict Relapse of Respiratory Failure in Critically Ill Patients with COVID-19, *Am J Respir Crit Care Med.* 2020 Oct 15; 202(8): 1173–1178. Published online 2020 Oct 15. doi: 10.1164/rccm.202005-1582LE
16. Takeshi Y et al, *CCM* 2013;41:536-545
17. Tonelli R et al., Early Inspiratory Effort Assessment by Esophageal Manometry Predicts Noninvasive Ventilation Outcome in De Novo Respiratory Failure A Pilot Study, *Am J Respir Crit Care Med.* 2020 Aug 15;202(4):558-567. doi: 10.1164/rccm.201912-2512OC
18. Gattinoni L et al., The Respiratory Drive: An Overlooked Tile of COVID-19 Pathophysiology, *Am J Respir Crit Care Med.* 2020 Oct 15;202(8):1079-1080. doi: 10.1164/rccm.202008-3142ED
19. Apigo et al. *Crit Care.* 2020 Jul 31;24(1):477
20. Carteaux G et al., Failure of Noninvasive Ventilation for De Novo Acute Hypoxemic Respiratory Failure: Role of Tidal Volume, *Crit Care Med.* 2016 Feb;44(2):282-90. doi: 10.1097/CCM.0000000000001379
21. Tobin MJ et al., *Am J Respir Crit Care Med.* 2020 Aug 1;202(3):356-360. doi: 10.1164/rccm.202006-2157CP.)
22. Vaporidi K et al., *Am J Respir Crit Care Med.* 2020 Jan 1;201(1):20-32. doi: 10.1164/rccm.201903-0596SO
23. Spinelli E et al., Respiratory drive in the acute respiratory distress syndrome: pathophysiology, monitoring, and therapeutic interventions, *Intensive Care Med.* 2020; 46(4): 606–618. Published online 2020 Feb 3. doi: 10.1007/s00134-020-05942-6
24. Bertoni M et al., A novel non-invasive method to detect excessively high respiratory effort and dynamic transpulmonary driving pressure during mechanical ventilation, *Crit Care.* 2019 Nov 6;23(1):346. doi: 10.1186/s13054-019-2617-0
25. Tomasino S et al, *Semin Cardiothorac Vasc Anesth* 2020;24:287-292
26. Moraes CCA et al. Bedside electrical impedance tomography unveils respiratory “chimera” in COVID-19. *Am J Respir Crit Care Med.* 2021;203(1):120-1
27. Kotani T et al., Roles of Electrical Impedance Tomography in Determining a Lung Protective Strategy for Acute Respiratory Distress Syndrome in the Era of Coronavirus Disease 2019, *JAMA J.* 2021;4(2):81-85. DOI: 10.31662/jmaj.2021-0014
28. Gattinoni L. et al., COVID-19 and ARDS: the baby lung size matters, *Intensive Care Med.* 2021; 47(1): 133–134. Published online 2020 Dec 4. doi: 10.1007/s00134-020-06324-8
29. Writing Group for the Alveolar Recruitment for Acute Respiratory Distress Syndrome Trial (ART) Investigators et al., ‘Effect of Lung Recruitment and Titrated Positive End-Expiratory Pressure (PEEP) vs Low PEEP on Mortality in Patients With Acute Respiratory Distress Syndrome: A Randomized Clinical Trial’, *JAMA*, 27 September 2017, <https://doi.org/10.1001/jama.2017.14171>
30. Zhao Z et al., ‘Positive End-Expiratory Pressure Titration with Electrical Impedance Tomography and Pressure–Volume Curve in Severe Acute Respiratory Distress Syndrome’, *Annals of Intensive Care* 9, no. 1 (December 2019): 7, <https://doi.org/10.1186/s13613-019-0484-0>

31. Zarantonello F et al., 'Prone Position and Lung Ventilation/Perfusion Matching in Acute Respiratory Failure Due to COVID-19', *American Journal of Respiratory and Critical Care Medicine*, 21 May 2020, rccm.202003-07751M, <https://doi.org/10.1164/rccm.202003-07751M>
32. Guérin C et al., 'A Prospective International Observational Prevalence Study on Prone Positioning of ARDS Patients: The APRONET (ARDS Prone Position Network) Study', *Intensive Care Medicine* 44, no. 1 (2018): 22-37, <https://doi.org/10.1007/s00134-017-4996-5>.
33. Wiggermann N et al., 'Prone Positioning of Patients With COVID-19: A Review of Equipment and Methods', *Human Factors*, 16 August 2020, 0018720820950532, <https://doi.org/10.1177/0018720820950532>
34. German S3 Guideline – Recommendations for the therapy of hospitalised patients with COVID-19, Version 4.1, February 2021
35. Surviving Sepsis Guidelines on the Management of Adults with Coronavirus Disease 2019 (COVID-19) in the ICU: First Update; March 2021; doi: 10.1097/CCM.0000000000004899.
36. ERS Guideline for the Management of hospitalized adults with coronavirus disease 2019 (COVID-19): A European Respiratory Society living guideline, January 2021, DOI: 10.1183/13993003.00048-2021
37. Australian guideline for clinical care of people with COVID-19, National COVID-19 Clinical Evidence Taskforce (Link)
38. Grieco DL et al., Effect of Helmet Noninvasive Ventilation vs High-Flow Nasal Oxygen on Days Free of Respiratory Support in Patients With COVID-19 and Moderate to Severe Hypoxemic Respiratory Failure: The HENIVOT Randomized Clinical Trial, *JAMA*. 2021 May 4;325(17):1731-1743. doi: 10.1001/jama.2021.4682
39. Perkins GD et al., RECOVERY- Respiratory Support: Respiratory Strategies for patients with suspected or proven COVID-19 respiratory failure; Continuous Positive Airway Pressure, High-flow Nasal Oxygen, and standard care: A structured summary of a study protocol for a randomised controlled trial, *Trials*. 2020 Jul 29;21(1):687. doi: 10.1186/s13063-020-04617-3
40. Patel BK et al., Effect of Noninvasive Ventilation Delivered by Helmet vs Face Mask on the Rate of Endotracheal Intubation in Patients With Acute Respiratory Distress Syndrome: A Randomized Clinical Trial, *JAMA*. 2016 Jun 14;315(22):2435-41. doi: 10.1001/jama.2016.6338
41. Leasa D et al., Knowledge translation tools to guide care of non-intubated patients with acute respiratory illness during the COVID-19 Pandemic, *Crit Care*. 2021 Jan 8;25(1):22. doi: 10.1186/s13054-020-03415-2
42. Papoutsis E et al, *Crit Care* 2021;25:121
43. Raoof S et al., High-Flow, Noninvasive Ventilation and Awake (nonintubated) prone in Patients with Coronavirus Disease 2019 with Respiratory Failure, *Chest*. 2020 Nov;158(5):1992-2002. doi: 10.1016/j.chest.2020.07.013. Epub 2020 Jul 15.]
44. Carteaux G et al., Failure of Noninvasive Ventilation for De Novo Acute Hypoxemic Respiratory Failure: Role of Tidal Volume, *Crit Care Med*. 2016 Feb;44(2):282-90. doi: 10.1097/CCM.0000000000001379
45. Pisano A et al., Indications for Tracheal Intubation in Patients With Coronavirus Disease 2019 (COVID-19), *J Cardiothorac Vasc Anesth*. 2021 May; 35(5): 1276-1280
46. Roca O et al., An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High-Flow Therapy, *Am J Respir Crit Care Med*. 2019 Jun 1;199(11):1368-1376. doi: 10.1164/rccm.201803-0589OC
47. Ricard JD et al., Use of nasal high flow oxygen during acute respiratory failure, *Intensive Care Med*. 2020 Sep 8 : 1-10. doi: 10.1007/s00134-020-06228-7
48. Zucman N et al., Prediction of outcome of nasal high flow use during COVID-19-related acute hypoxemic respiratory failure, *Intensive Care Med*. 2020 Jul 15 : 1-3. doi: 10.1007/s00134-020-06177-1
49. Duan J et al., Assessment of heart rate, acidosis, consciousness, oxygenation, and respiratory rate to predict noninvasive ventilation failure in hypoxemic patients, *Intensive Care Med*. 2017 Feb;43(2):192-199. doi: 10.1007/s00134-016-4601-3. Epub 2016 Nov 3.]
50. Brault C et al, *AJRCCM* 2020, doi 10.1164/rccm.202005-2025LE
51. Telias I et al., 'Airway Occlusion Pressure As an Estimate of Respiratory Drive and Inspiratory Effort during Assisted Ventilation', *American Journal of Respiratory and Critical Care Medicine* 201, no. 9 (1 May 2020): 1086-98, <https://doi.org/10.1164/rccm.201907-1425OC>.
52. Positionspapier zur praktischen Umsetzung der apparativen Differenzialtherapie der akuten respiratorischen Insuffizienz bei COVID-19: Deutsche Gesellschaft für Pneumologie und Beatmungsmedizin e.V. (DGP), https://www.researchgate.net/publication/340865640_Positionspapier_zur_praktischen_Umsetzung_der_apparativen_Differenzialtherapie_der_akuten_respiratorischen_Insuffizienz_bei_COVID-19_Deutsche_Gesellschaft_fur_Pneumologie_und_Beatmungsmedizin_eV_DGP
53. Rose L et al., 'Automated versus Non-Automated Weaning for Reducing the Duration of Mechanical Ventilation for Critically Ill Adults and Children', *Cochrane Database Syst Rev* 6 (2014).
54. Lellouche F et al., 'A Multicenter Randomized Trial of Computer-Driven Protocolized Weaning from Mechanical Ventilation', *American Journal of Respiratory and Critical Care Medicine* 174, no. 8 (15 October 2006): 894-900, <https://doi.org/10.1164/rccm.200511-1780OC>.
55. Taenaka H. et al., Individualized ventilatory management in patients with COVID-19-associated acute respiratory distress syndrome, *Respir Med Case Rep*. 2021; 33: 101433. Published online 2021 May 31. doi: 10.1016/j.rmcr.2021.101433
56. Zhao Z. et al., The use of electrical impedance tomography for individualized ventilation strategy in COVID-19: a case report, *BMC Pulm Med*. 2021; 21: 38. Published online 2021 Jan 22. doi: 10.1186/s12890-021-01411-y

IMPRINT

GERMANY
 Drägerwerk AG & Co. KGaA
 Moislinger Allee 53–55
 23542 Lübeck

www.draeger.com