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www.draeger.com

REFERENCIAS:

- 1 Organización Mundial de la Salud, Oficina regional europea, <http://www.euro.who.int/en/health-topics/disease-prevention/nutrition/a-healthy-lifestyle/body-mass-index-bmi>; consulta del 3 de noviembre de 2017.
- 2 Rossaint R., Werner CH.: *Die Anästhesiologie: Allgemeine und spezielle Anästhesiologie, Schmerztherapie und Intensivmedizin*, 3. Auflage, Springer Verlag.
- 3 Salome CM, Munoz PA, Berend N, et al. Effect of obesity on breathlessness and airway responsiveness to methacholine in non-asthmatic subjects. *Int J Obes (Lond)* 2008;32:502-9.
- 4 Eriksen J, Andersen J, Rasmussen JP, et al. Effects of ventilation with large tidal volumes or positive end-expiratory pressure on cardiorespiratory function in anesthetized obese patients. *Acta Anaesthesiol Scand* 1978;22:241-8.
- 5 Hackney JD, Crane MG, Collier CC, et al. Syndrome of extreme obesity and hypoventilation: studies of etiology. *Ann Intern Med* 1959;51:541-52.
- 6 Pelosi P, Croci M, Ravagnan I, et al. The effects of body mass on lung volumes, respiratory mechanics, and gas exchange during general anesthesia. *Anesth Analg* 1998;87:654-60.
- 7 Hodgson LE, Murphy PB, Hart N. Respiratory management of the obese patient undergoing surgery. *J Thorac Dis*. Mayo de 2015;7(5):943-52.
- 8 Ferretti A, Giampiccolo P, Cavalli A, et al. Expiratory flow limitation and orthopnea in massively obese subjects. *Chest* 2001;119:1401-8.
- 9 Steier J, Jolley CJ, Seymour J, et al. Neural respiratory drive in obesity. *Thorax*. Agosto de 2009;64(8):719-25.
- 10 Pankow W, Podszus T, Gutheil T, et al. Expiratory flow limitation and intrinsic positive end-expiratory pressure in obesity. *J Appl Physiol (1985)* 1998;85:1236-43.
- 11 Smetana GW. Preoperative pulmonary evaluation. *N Engl J Med* 1999;340:937-44.
- 12 Reinius H, Jonsson L, Gustafsson S, et al. Prevention of atelectasis in morbidly obese patients during general anesthesia and paralysis: a computerized tomography study. *Anesthesiology* 2009;111:979-87.
- 13 Eichenberger A, Proietti S, Wicky S, et al. Morbid obesity and postoperative pulmonary atelectasias: an underestimated problem. *Anesth Analg*. Diciembre de 2002;95(6):1788-92.
- 14 Salome CM.: Physiology of obesity and effects on lung function; *J Appl Physiol* 108: 206–211, 2010.
- 15 Rubinstein I, Zamel N, DuBarry L, et al. Airflow limitation in morbidly obese, nonsmoking men. *Ann Intern Med* 1990;112:828-32.
- 16 Tanoubi I, Drolet P, Donati F. Optimizing preoxygenation in adults. *Can J Anaesth*. Junio de 2009;56(6):449-66.
- 17 Nestler, C., Wrigge H. et al.: Individualized positive end-expiratory pressure in obese patients during general anesthesia: a randomized controlled clinical trial using electrical impedance tomography; *Br J Anaesth*. 16 de octubre de 2017. doi: 10.1093/bja/aex192. [Publicación electrónica antes de la impresión].
- 18 Schumann, R.: Pulmonary Physiology of the Morbidly Obese and the Effects of Anesthesia. *International Anesthesiology Clinics*. 51 (3): 41-51, Verano 2013 DOI: 10.1097/AIA.0b013e3182981252.
- 19 Sharp JT, Henry JP, Sweany SK, et al. The total work of breathing in normal and obese men. *J Clin Invest* 1964;43:728-39.
- 20 Zerah F, Harf A, Perlemuter L, et al. Effects of obesity on respiratory resistance. *Chest* 1993;103:1470-6.
- 21 Kress JP, Pohlman AS, Alverdy J, et al. The impact of morbid obesity on oxygen cost of breathing (VO(2)RESP) at rest. *Am J Respir Crit Care Med* 1999;160:883-6.
- 22 Mandal S, Hart N. Respiratory complications of obesity. *Clin Med* 2012;12:75-8.
- 23 Young T, Palta M, Dempsey J, et al. The occurrence of sleep-disordered breathing among middle-aged adults. *N Engl J Med* 1993;328:1230-5.
- 24 Young T, Peppard PE, Taheri S. Excess weight and sleep-disordered breathing. *J Appl Physiol (1985)* 2005;99:1592-9.
- 25 Resta O, Foschino-Barbaro MP, Legari G, et al. Sleep-related breathing disorders, loud snoring and excessive daytime sleepiness in obese subjects. *Int J Obes Relat Metab Disord* 2001;25:669-75.
- 26 Lam KK et al.: Obstructive sleep apnea, pain, and opioids: is the riddle solved? *Curr Opin Anaesthesiol*. Febrero de 2016; 29(1): 134–140. Publicado online el 28 de nov. 2015. doi:10.1097/ACO.0000000000000265 PMID: PMC4927322.
- 27 Serpa Neto A, Cardoso SO, Manetta JA, Pereira VG, Esposito DC, Pasqualucci Mde O, et al. Association between use of lung-protective ventilation with lower tidal volumes and clinical outcomes among patients without acute respiratory distress syndrome: a meta-analysis. *JAMA*. 2012;308(16):1651–9.
- 28 Hemmes SN, Serpa Neto A, Schultz MJ. Intraoperative ventilatory strategies to prevent postoperative pulmonary complications: a meta-analysis. *Curr Opin Anaesthesiol*. 2013;26(2):126–33.
- 29 Murphy C et al.: Airway management and oxygenation in obese patients; *Can J Anesth/J Can Anesth (2013)* 60:929–945.

REFERENCIAS:

- 30 Ladhka K. et al.: Intraoperative protective mechanical ventilation and risk of postoperative respiratory complications: hospital based registry study, *BMJ* 2015;351:h3646.
- 31 Futier E et al.: Noninvasive Ventilation and Alveolar Recruitment Maneuver Improve Respiratory Function during and after Intubation of Morbidly Obese Patients; *Anesthesiology* 2011; 114:1354 – 63.
- 32 Ball L, Pelosi P: Intraoperative mechanical ventilation in patients with non-injured lungs: time to talk about tailored protective ventilation?; *Ann Transl Med.* Enero de 2016;4(1):17.
- 33 Futier E, Marret E, Jaber S; Perioperative Positive Pressure Ventilation: An Integrated Approach to Improve Pulmonary Care, *Anesthesiology.* Agosto de 2014;121(2):400-8.
- 34 Andreas Güldner y Marcelo Gama de Abreu: Intraoperative protective ventilation reduces postoperative pulmonary complications, *Anästhesiol Intensivmed Notfallmed Schmerzther* 2015; 50(9): 524-528.
- 35 Entrevista con el Prof. Dr. Hermann Wrigge, University Hospital Leipzig, Klinik und Poliklinik für Anästhesiologie und Intensivtherapie, 03.11.2016.
- 36 Rothen HU, Neumann P, Berglund J, Valtysson J, Magnusson A, Hedenstierna G.; Dynamics of re-expansion of atelectasis during general anaesthesia.; *Br J Anaesth.* 1999 Apr;82(4):551-6.
- 37 Tusman G, Groisman I, Fiolo FE, Scandurra A, Arca JM, Krumrick G, Bohm SH, Sipmann FS; Noninvasive monitoring of lung recruitment maneuvers in morbidly obese patients: the role of pulse oximetry and volumetric capnography. *Anesth Analg.* Enero de 2014;118(1):137-44.
- 38 Pelosi P, Gama de Abreu M, Rocco PR; New and conventional strategies for lung recruitment in acute respiratory distress syndrome. *Crit Care.* 2010;14(2):210.
- 39 Suarez-Sipmann, Böhm SH, Tusman G, Pesch T, Thamm O, Reissmann H, Reske A, Magnusson A, Hedenstierna G.; Use of dynamic compliance for open lung positive end-expiratory pressure titration in an experimental study.; *Crit. Care Med.* Enero de 2007; 35(1):214-21.
- 40 Hodgson C et al. A randomised controlled trial of an open lung strategy with staircase recruitment, titrated PEEP and targeted low airway pressures in patients with acute respiratory distress syndrome, *Crit Care.* 2011; 15(3): R133. Publicación online del 2 de junio de 2011.
- 41 Santos RS, Moraes L, Samary CS, Santos CL, Ramos MB, Vasconcellos AP, Horta LF, Morales MM, Capelozzi VL, Garcia CS, Marini JJ, Gama de Abreu M, Pelosi P, Silva PL, Rocco PR; Fast Versus Slow Recruitment Maneuver at Different Degrees of Acute Lung Inflammation Induced by Experimental Sepsis, *Anesth Analg.* Abril de 2016;122(4):1089-100. doi: 10.1213/ANE.0000000000001173.
- 42 Futier E et al.: Positive end-expiratory pressure improves end-expiratory lung volume but not oxygenation after induction of anaesthesia; *Eur J Anaesthesiol* 2010; 27:508-513.
- 43 Barbosa FT et al.: Positive end-expiratory pressure (PEEP) during anaesthesia for prevention of mortality and postoperative pulmonary complications (artículo de revisión); 2014, *The Cochrane Collaboration.* Publicado por John Wiley & Sons, Ltd.
- 44 Brismar B, Hedenstierna G, Lundquist H, et al. Pulmonary densities during anesthesia with muscular relaxation—a proposal of atelectasis. *Anesthesiology* 1985;62:422-8.
- 45 Neumann P, Rothen HU, Berglund JE, et al. Positive end-expiratory pressure prevents atelectasis during general anaesthesia even in the presence of a high inspired oxygen concentration. *Acta Anaesthesiol Scand* 1999;43:295-301.
- 46 Tokics L, Hedenstierna G, Strandberg A, et al. Lung collapse and gas exchange during general anesthesia: effects of spontaneous breathing, muscle paralysis, and positive end-expiratory pressure. *Anesthesiology* 1987;66:157-67.
- 47 Clarke JP, Schuitemaker MN, Sleight JW. The effect of intraoperative ventilation strategies on perioperative atelectasis. *Anaesth Intensive Care* 1998;26:262-6.
- 48 Maracajá-Neto LF, Verçosa N, Roncally AC, et al. Beneficial effects of high positive end-expiratory pressure in lung respiratory mechanics during laparoscopic surgery. *Acta Anaesthesiol Scand* 2009;53:210-7.
- 49 Meininger D, Byhahn C, Mierdl S., et al. Positive end-expiratory pressure improves arterial oxygenation during prolonged pneumoperitoneum. *Acta Anaesthesiol Scand* 2005;49:778-83.
- 50 Futier E, Constantin JM, Paugam-Burtz C, et al. A trial of intraoperative low-tidal-volume ventilation in abdominal surgery. *N Engl J Med.* 2013;369:428-37.
- 51 Talab HF, Zabani IA, Abdelrahman HS, Bukhari WL, Mamoun I, Ashour MA, Sadeq BB, El Sayed SI; Intraoperative ventilatory strategies for prevention of pulmonary atelectasis in obese patients undergoing laparoscopic bariatric surgery. *Anesth Analg.* Noviembre de 2009;109(5):1511-6.
- 52 Bluth T et al.: Protective intraoperative ventilation with higher versus lower levels of PEEP in obese patients (PROBESE): study protocol for a randomized controlled trial; *Trials.* 28 de abril de 2017; 18(1):202. doi: 10.1186/s13063-017-1929-0.