

Case Study – Long-Term Ventilation With Variable PSV

Case: A 73-year old male patient, while working in his garden, fell about 12 feet from a ladder. He was found by his wife who called an ambulance. When the paramedics arrived, he was observed to have an oxygen saturation SpO_2 of 65 %. The patient was then intubated and ventilated.




DIAGNOSIS AND INITIAL TREATMENT

Due to a suspected right-sided tension pneumothorax, thoracic suction was applied and the patient brought to the emergency department of the University Hospital Dresden, where a radiological examination was performed. The examination revealed a right-side hemopneumothorax, bilateral rib fractures, a sternal fracture, and fractures of thoracic vertebra T4 and T11. After surgical treatment and internal stabilization of the 11th thoracic vertebra, the patient was moved to the surgical intensive care unit where volume-controlled ventilation was applied. He was placed under the care of an intensivist. As a result of the pleural effusions, atelectasis, and tracheobronchitis, the patient started showing signs of respiratory insufficiency with pronounced hypoxemia ($PaO_2/FiO_2=173$). As expected, proper oxygenation and ventilation would be difficult to maintain because the multiple fractures, unstable thorax, and resultant paradoxical breathing.

METHOD

The spontaneous breathing of the patient was first supported with conventional pressure support ventilation (PSV) and a constant support pressure of 10 cmH_2O , flow trigger = 3.0 l/min, PEEP = 12 cmH_2O and $FiO_2 = 0.4$. This made it possible to achieve an average tidal volume of 514 ml and a respiratory frequency of 21/min, which in turn resulted in a minute volume of 10.8 l/min, a PaO_2/FiO_2 of 250 and a $PaCO_2$ of 53.1 mmHg. Given the reduced variability in tidal volume, which was about 10% and thus significantly lower than the variability of the tidal volume in spontaneously breathing patients (approximately 26%)¹, as well as the difficulty of finding the right support pressure for the patient, which at the same time provided some relief and allowed an adequate gas exchange, the ventilation mode was switched to variable pressure support.

Window for activating Variable Pressure Support in the Evita Infinity V500 ventilator. In addition to setting the support pressure (ΔP_{supp}), as with conventional pressure support ventilation, the variability of the pressure support (Pressure Var.) can be set in a range from 0-100%. Pressure Var. (in percentage) refers to the range in which ΔP_{supp} can vary. Since ΔP_{supp} values are subject to normal distribution, the coefficient of variation is approx. 17% (50%/3) at a Pressure Var = 50%.



D-62688-2012



The chest x-ray revealed a left-side pleural effusion with adjacent partial atelectasis, right basal partial atelectasis, a right thorax suction drainage with the tip projecting toward the right midfield and right thoracic subcutaneous emphysema.



Six weeks after the previous admittance and assisted spontaneous breathing with pressure support ventilation, the lung conditions improved and the patient was transferred to a rehabilitation clinic.

Mechanical ventilation was supported with a Dräger Evita Infinity V500 SW 2.n* ventilator. The ventilator is equipped with the new Variable Pressure Support option. Variable pressure support provides the ability to apply distributed pressure support within defined limits stochastically, i.e., randomly, as described by Gama de Abreu and coworkers². In case of variable pressure support, the desired degree of pressure support variability is set in addition to the usual pressure support parameters (e.g., ramp, cycle off criterion, trigger and pressure support). This results in a breathing pattern where the level of pressure support varies from breath to breath.

The support pressure variability was initially set to 90% (range adjustable from 0-100%). With a value, for example, of 100%, the pressure support would amount to a minimum of 0 cmH₂O and a maximum of twice the pressure support value, with the extreme values occurring only very rarely. The specified pressure support variability of 90% resulted in a coefficient of variation of 30% in the support pressure. The reason for this is that, with a normal distribution, these two variables (variability of the support pressure and coefficient of variation) have a ratio of about 3:1. Since all the other ventilator settings remained unchanged, this setting resulted in a tidal volume of approximately 530 ml (average).

RESULT

Compliance was improved by approximately 20% during the initial first hour, while the gas exchange parameters with a PaO₂/FiO₂ of 248 and a PaCO₂ of 50 mmHg remained almost constant. Over a period of approximately seven hours, the pressure support variability was reduced to 60% and maintained at that level.

Over a period of ten days, the patient was treated successfully with a variable pressure support of 8 cmH₂O and variability in the range of 40-60% while continuing to suffer from persistent thoracic instability. The patient was tracheotomized before his transfer to a rehabilitation hospital. Conventional pressure support ventilation was used then only intermittently to reduce the work of breathing, since the severely fatigued respiratory muscles would recover over a lengthy period of time.

* in countries the product name ist Infinity Acute Care System Workstation Critical Care

SUMMARY

The critical finding of this case is that while using variable pressure support - different, randomly generated pressure support values could be applied during a long-term weaning period. This was made possible in spite of the variability of the patient's spontaneous breathing being affected by the unstable thorax. The use of variable pressure support avoided having to make frequent changes to the ventilator settings to ensure the patient felt comfortable, the gas exchange was adequate, and compliance of the respiratory system improved as well.

Key conclusions:

- Variable pressure support leads to increased variability of the tidal volume which is similar to spontaneous respiration variability.
- Using variable pressure support required fewer ventilator adjustments and thus relieves staff from constant need to adjust parameters.
- With variable pressure support, it was possible to cover a wide range of ventilation pressures by optimizing the gas exchange and increasing patient comfort levels.

Expectations for the future use of variable pressure support

- Variable pressure support makes it possible to reduce pressure support compared to traditional pressure support ventilation, thus allowing for faster weaning.
- Since different pressure support values are generated, the always changing needs of the patient could be better met in terms of respiratory support.
- Variable pressure support could help to reduce the number of necessary assisted ventilation adjustments.
- Variable pressure support could replace traditional pressure support ventilation.

Brief description of the general research situation

According to experimental studies, variable pressure support leads to better oxygenation, lower airway pressures, and fewer ventilator-associated lung injuries compared to traditional pressure support ventilation.^{2,3,4}

Additionally, the breathing effort with variable pressure support is lower than with traditional pressure support ventilation and proportional assist ventilation (proportional pressure support).⁵

Early clinical experiences suggest that variable pressure support can cover a wide range of pressure support to meet the needs of the patient.

LITERATUR

- | | | | |
|----|--|----|--|
| 1. | Tobin MJ, Mador MJ, Guenther SM, Lodato RF, Sackner MA: Variability of resting respiratory drive and timing in healthy subjects. <i>Journal of Applied Physiology</i> 1988; 65: 309-317 | 4. | Carvalho AR, Spieth PM, Güldner A, Cuevas M, Carvalho NC, Beda A, Spieth S, Stroczyński C, Wiedemann B, Koch T, Pelosi P, de Abreu MG: Distribution of regional lung aeration and perfusion during conventional and noisy pressure support ventilation in experimental lung injury. <i>J Appl Physiol</i> 2011; 110: 1083-92 |
| 2. | Gama de Abreu M, Spieth P, Pelosi P, Carvalho AR, Walter C, Schreiber-Ferstl A, Aikele P, Neykova B, Hübler M, Koch T: Noisy pressure support ventilation: A pilot study on a new assisted ventilation mode in experimental lung injury. <i>Crit Care Med</i> 2008; 36: 818-827 | 5. | Spieth P, Güldner A, Beda A, Carvalho NC, Nowack T, Krause A, Rentzsch I, Suchantke S, Thal S, Engelhard K, Kasper M, Koch E, Pelosi P, Gama de Abreu M: Comparative effects of proportional assist and variable pressure support ventilation on lung function and damage in experimental lung injury. <i>Crit Care Med</i> ; in press |
| 3. | Spieth PM, Carvalho AR, Güldner A, Kasper M, Schubert R, Carvalho NC, Beda A, Dassow C, Uhlig S, Koch T, Pelosi P, Gama de Abreu M: Pressure support improves oxygenation and lung protection compared to pressure-controlled ventilation and is further improved by random variation of pressure support. <i>Crit Care Med</i> 2011; 39: 746-55 | | |



Prof. Dr. med. habil. Marcelo Gama de Abreu, DEWA

Managing Chief Physician at the Department of Anesthesiology and Intensive Care Medicine (Director: Prof. Dr. med. Thea Koch), University Hospital Carl Gustav Carus, TU Dresden

Academic Career

- Electronic Engineering Degree (Federal University of Rio de Janeiro – UFRJ)
- Medical Studies (University of the State of Rio de Janeiro – UERJ)
- Master of Science in Biomedical Engineering (Federal University of Rio de Janeiro – UFRJ)
- PhD at the University of Heidelberg
- Post-doctoral work and supernumerary professor at the Technical University of Dresden

Professional Career

- Chief physician at the Department of Anesthesiology and Intensive Care, University Hospital Carl Gustav Carus (Medical Director: Prof. Dr. med. T. Koch)
- Executive Senior Physician for Clinical and Experimental Research at the Department of Anesthesiology and Intensive Care, University Hospital Carl Gustav Carus (Medical Director: Prof. Dr. med. T. Koch)
- Managing Senior Physician at the Department of Anesthesiology and Intensive Care Medicine at the University Hospital Carl Gustav Carus (Medical Director: Prof. Dr. med. T. Koch)
- Member of management of the Department of Anesthesiology and Intensive Care, University Hospital Carl Gustav Carus, Dresden, Germany
- Member of the management of the University Hospital Carl Gustav Carus (advisory body to the medical board of the University Hospital Carl Gustav Carus)
- Elected member of the Research Committee of the European Society of Anesthesiology, Brussels, Belgium
- Chairman of the Subcommittee Respiration (SC5) of the European Society of Anesthesiology, Brussels, Belgium

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany

As of August 2015:

Dräger Medical GmbH changes to Drägerwerk AG & Co. KGaA.

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com