

Dräger SmartPilot® View Software – visualisation of anaesthetic levels

The software visualises the complex synergies of anaesthetic drugs and predicts their effects based on pharmacodynamic modelling for both the current status and the prospective course of general anaesthesia. SmartPilot View turns abstract device data into comprehensive visual information to support the decision on more precise and patient-optimised titration of anaesthetic drugs.

Improved drug dosing

Decision-making support

Improved workflow efficiency

D-46302-2021

Comprehensive visualisation

Optimised OR and PACU utilisation

Benefits

Improved drug dosing

SmartPilot View is a software which calculates and predicts the combined effects of commonly used analgesic and hypnotic drugs based on population models. It can support the clinician in intraoperatively administering drugs in a patient-individualized and optimized manner what may contribute to:

- better hemodynamic stability
- better controlled patient condition, both during induction and during maintenance
- avoiding over-dosage and a well-managed recovery (faster awakening) to reduce risk of postoperative delirium

Decision-making support

SmartPilot View supports the decision-making by turning complex pharmacokinetic and pharmacodynamic calculations into an easy-to-interpret visual information:

- displays the current anaesthesia level and hypnotic-analgesic balance (orange dot)
- displays 20 minute prediction for the anaesthesia level (white arrow)
- communicates with selected infusion pumps and based on the current pump rates and the vaporizer settings the "What if" function shows what effect a rate change or bolus would have
- build-in pre-prediction ("What-if" function) for IV-pump rates and anaesthesia device settings

Improved workflow efficiency

SmartPilot View supports our workflow by:

- easing the documentation of frequently performed user interventions in the logfiles over the entire surgical procedure via the event markers, e.g. LOC (Loss of Consciousness), Intub. (Intubation), Cut, Move, Extub. (Extubation) and Other. The event markers can be moved and deleted on the time scale at any time.
- automatically receiving drug dosing information from Dräger anaesthesia machines and connected compatible syringe pumps to reduce the manual data entry to a minimum

Comprehensive visualisation

SmartPilot View allows the anaesthesiologist instant access to the status of anaesthesia and the administered drugs during the whole surgical procedure. Following parameters support the intuitive and fast compilation of patient status:

- Isobole lines in 2D diagram representing areas of probability how patient will react to a defined pain stimulus, and indicating level of anaesthesia e.g. TOL 90, TOL 50, TOSS
- history function on the 2D diagram to track the anaesthesia level from the beginning of the case
- calculated Noxious Stimulus Response Index (NSRI) displayed as a numeric value and as trend graph
- display of the past (40 min), present and preview (20 min) of calculated effect-site concentrations of administered hypnotics, opioids, and muscle relaxants

Benefits

- Bispectralindex (BIS) signal quality index SQI in %
- continuity of case from induction room into main operation room

Optimised OR and PACU utilisation

SmartPilot View can contribute to optimising OR and PACU utilisation times:

- bargraph indicator to estimate the wake-up time
- patient-optimised drug dosing may help to reduce risk of overdosing and, hence, prolonged PACU stay

Specific features in combination with Dräger Anaesthesia Workstations

The SmartPilot View runs on medical-grade touch panel PCs* and can be combined with Dräger anaesthesia devices as e.g. Primus® family, Perseus® A500, Zeus® Infinity® Empowered and Atlan® family. These combinations can offer further workflow enhancements:

- data transfer of measured ventilation and gas concentration values
- "What if" function to show the calculated effect of a volatile agent setting before it is even confirmed (only in combination with Zeus IE)
- automatic transfer of Vapor® 3000 and D-Vapor® 3000 hand dial setting as well as measured ventilation and gas concentration values (only in combination with Perseus A500 comprising VaporView functionality)

* Exception: integrated SmartPilot View option of Zeus IE runs on the main screen of the Zeus IE anaesthesia machine.

System Components



D-36084-2021

Medical Grade PC

The panel PC is a state-of-the-art medical grade computer that supports clinical software applications and networked information at the point of care to provide relevant information for clinical decision making. With its quiet and hygienic fanless design it is ideal for use in direct patient vicinity.



D-30739-2017

Infinity® Acute Care System

Transform your clinical workflow with Infinity® Acute Care System. Its multiparameter monitor integrates with its networked medical-grade workstation, giving you real-time vital signs, access to clinical hospital systems and data management applications for a comprehensive range of patient information and powerful analysis tools at the point-of-care.



MT-1105-2007

Infinity® Omega Solution

This two-screen solution provides access to images, lab results and other clinical data, while displaying real-time vital signs data at the point of care on a 20" touchscreen. Bedside/transport monitor continuously monitors patient during transport and at the bedside.

Technical Data

SmartPilot® View

Software version	3.0n
------------------	------

Classification

Classification in accordance with Directive 93/42/EEC Annex IX	Class IIb
--	-----------

General

SmartPilot® View displays pharmacokinetic, pharmacodynamic (PK/PD) modelling information for use with total intravenous (TIVA), balanced anaesthesia and volatile anaesthesia.

Display of haemodynamic measurements and processed EEG (BIS® Aspect).

Drug and dosage data entry for intravenous drugs: manually or automatically from the syringe pumps. For inhaled anaesthetics: automatically from anaesthesia device.

Event markers allow the patient's individual status to be documented with regard to the observed anaesthesia level. The markers may be deleted or repositioned if needed.

Support of USB stick for case recording.

Editable drug database.

PK/PD models

Interaction between volatile and intravenous hypnotics and opioids visualised through two-dimensional graph with 15 minute prediction.

PK models for each single drug visualised in a time-based graph with 40 minute trend and 20 minute prediction.

Manual entry of drug bolus and rates at current time or retrospectively followed by recalculation of the case, incl. prediction.

PK models

Propofol	Schnider
Remifentanyl	Minto
Alfentanil	Scott
Fentanyl	Scott
Sufentanil	Gepts
Volatiles	Bailey
Pancuronium	Rupp
Rocuronium	Wierda
Vecuronium	Cronnelly
Mivacurium	Laurin

PD models

Sequential model	Bouillon
------------------	----------

Supported drugs

Hypnotics:	
Volatile drugs	Desflurane, Isoflurane, Sevoflurane
Intravenous drugs	Propofol
Opioids (intravenous)	Fentanyl, Remifentanyl, Alfentanil, Sufentanil
Muscle relaxants*	Pancuronium, Rocuronium, Mivacurium, Vecuronium

*The impact of muscle relaxants on the drug interaction is not taken into account.

Patient demographic data

Height	150–200 cm (59.1 to 78.7 in)
Weight	40–140 kg (89 to 308 lbs)
Age	18–90 years
Body mass index (BMI)	<= 35 ASA Level < IV

The limitations of the patient demographic data are due to the calculation of lean body mass using the formula from James [James WPT. Research on obesity. London. Her Majesty's Stationery Office. (ISBN 0-11-4500347). 1976].

Technical Data

System compatibility

Anaesthesia devices	Zeus®, Zeus® IE, Perseus® A500, Primus®, Primus® IE, Atlan® family
PC-platform	Dräger C700 for IT with or without Infinity® Explorer or other medical-grade PC (system requirements): – W7 32/64bit – 1-GHz-CPU (x86-x64) – 1 GB RAM (at least 500 MB available RAM at startup of SPV) – 1 GB available Hard Disk Drive space – Screen resolution: 1,280 x 768
Monitoring	Dräger Infinity® Delta, Dräger Infinity® Delta XL, Dräger Zeus® / Zeus® Infinity® Empowered with integrated haemodynamic monitoring, IACS patient monitor, Philips monitor (MP50, MP70, MX800)

Display of patient measurements and processed EEG (BIS® Aspect) provided by the basic device:

Heart rate	1/min
Blood pressure (ART M or NIBP)	mmHg or kPa
etCO ₂	mmHg or kPa
BIS Index and SQI	

Infusion pumps

Dräger	IVenus**
B. Braun	Perfusor® Space with SpaceCom**
Fresenius	Orchestra® Modul DPS with Base A
Fresenius	Orchestra® Modul DPS with Base Primea*
Fresenius	Injectomat® Agilia with Link + Agilia
Alaris®	GH, PK*, TIVA with Gateway
Terumo	Terumo Terufusion® TE-332, TE-371, TE-372
Terumo	Terumo Terufusion® TE-SS800 with TE-RS800
Syramed	µSP6000 with UniQue Concept Box

Range of supported devices (e.g. Patient Monitors, Anaesthesia devices and infusion pumps) is expanding: please contact your Dräger representative for the current list of compatible devices.

*TCl pumps are supported

**up to four infusion pumps

Perfusor® Space with SpaceCom is the registered trade mark of Company B.Braun.

Orchestra® Modul DPS with Base A, Orchestra® Modul DPS with Base Primea and Injectomat® Agilia with Link Agilia are registered trade mark of Company Fresenius.

Terumo Terufusion® TE-332, TE-371, TE- 372 and Terumo Terufusion® TE-SS880 with TE-RS800 are registered trade mark of Company Terumo.

Notes

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06460-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel: +65 6872 9288
Fax: +65 6259 0398
asia.pacific@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

