Alarm Monitoring System GMS-Gateway/ GMS-Gateway VoluCount™

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**Instructions for Use
Installation Instructions**

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Dräger Medical GmbH

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use/ Installation Instructions

Any use of the device requires full understanding and strict observation of these instructions. The device is only to be used for purposes specified here.

Maintenance

The device must be inspected and serviced regularly by trained service personnel at six-monthly intervals (and a record kept). Repair and general overhaul of the device may only be carried out by trained service personnel. We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Only authentic Dräger spare parts may be used for maintenance. Observe chapter "Maintenance Intervals".

DANGER

This device is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Liability for proper function or damage

The liability for the proper function of the device is irrevocably transferred to the owner or operator to the extent that the device is serviced or repaired by personnel not employed or authorized by DrägerService or if the device is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

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Intended Use

GMS Gateway

Interface of the Dräger Alarm Monitoring System to external systems using a LonWorks[®] interface for transmission of alarm and operational signals from the central gas supply system to a visualisation system (e.g. to a building management system).

NOTE

The GMS gateway is not intended for use in the patient environment or in rooms used for medical purposes. It is not approved for use within the area control unit.

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GMS Gateway VoluCount™

Interface of the Dräger Alarm Monitoring System to external systems using a LonWorks[®] interface for transmission of the gas consumption values to a visualisation system (e.g. to a building management system).


NOTE

The GMS gateway VoluCount™ is not intended for use in the patient environment or in rooms used for medical purposes. It is not approved for use within the area control unit.

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- 1 Service LED, LonWorks® side
- 2 Service LED, GMS side
- 3 Status DUO LED, LonWorks® side
- 4 Status DUO LED, GMS side
- 5 24 V LED

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
Installation

The GMS gateway and the GMS gateway VoluCount™ are constructed as housings for DIN rail installation. The unit is installed in a switch box or in a separately available plastic housing (G 50 748) by:

- Snapping onto the top hat rail (35 mm).
- Fastening of the connector terminals as labeled.

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
NOTE

The shielding of the G 50 750 GM network system cable is not connected to the GMS gateway resp. GMS gateway VoluCount™. If the equipment is supplied via system cable, ensure that the shielding is connected at the other end of the cable.

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Testing and Commissioning

The GMS gateway and the GMS gateway VoluCount™ are components of the Dräger Alarm Monitoring System.

NOTE

The system can only go into operation after an acceptance test has been carried out by trained and qualified personnel.

Note national regulations!

Upon completion of the installation or maintenance work:

- A test and acceptance program must be carried out by specialists.

These tests ensure:

- that the safety requirements designed to protect patients and staff have been met and
- that the performance characteristics of the central gas supply system and the Dräger Gas Management are met.

- Test results are to be documented in written form.

Following these acceptance tests, the operational device is handed over to the user together with the documentation, and the operating staff is instructed on how to operate the system.

- The handover procedure is to be protocolled.

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Variante 1 (G 51 297):

"GMS Gateway A V1_01.XIF"
(digitales Interface=Alarm)
ID 94:52:A9:46:00:09:04:02

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LonWorks® Interface

GMS Gateway

The GMS gateway is constructed using the following objects:

- Node object (object 0)
- Data collectors object (object 1)
- Area control units object (object 2)

Device Templates

For each interface version a specialised device template is available in the form of an XIF file.

Version 1 (G 51 297):

"GMS Gateway A V1_01.XIF"
(digital interface=alarm)
ID 94:52:A9:46:00:09:04:02

Version 2 (G 51 296):

"GMS Gateway AP V2_00.XIF"
(mixed digital/analog interface = alarm + pressure)
ID 94:52:A9:46:00:09:04:03

GMS Gateway VoluCount™

The GMS gateway VoluCount™ is constructed using the following objects:

- Node object (object 0)
- Volucount object (object 1):
A manufacturer-specific object of the Dräger Medical GmbH, as there is no standard LonMark® object available for the required function.

Device Templates

For the GMS gateway VoluCount™ a device template is available in the form of an XIF file.

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
Node Object

The node object serves to monitor all objects in the node regarding self-tests, area inspections, error control, file transfer, etc.

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Prescribed Network Variables

A request for the status and checks of individual node objects is available through nviRequest in the format SNVT_obj_request. These include:

- Activation and deactivation of single or all node objects
- Status request of individual node objects
- Deletion of the status bits of all node objects
- Request for self-testing of all node objects
- Display of error and alarm messages of the objects

The requests "RQ_NORMAL", "RQ_REPORT_MASK", and "RQ_UPDATE_STATUS" are supported for all objects.

The results of the requested tests and settings with nviRequest are reported with nvoStatus in format SNVT_obj_status as bit information.

The node object supports the status "in alarm" and "comm_failure".

If the interprocessor communication of the assembly group is interrupted, both error messages are sent after the corresponding request.

Optional Network Variables

The nvoFileDirectory in format SNVT_address contains the start address of the configuration area.

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Network Variables in the Manufacturer-specific Area

The output variable nvoGmsState in the format SNVT_state is located in the manufacturer-specific area.

The bits are interpreted as follows:

Bit	Value 0	Value 1
0	Intercommunication of the controllers is OK	Communication not possible / defective
1	All devices on the GMS side can be polled without error	At least one device is unreachable
2-15	Not used	

Configuration Properties of the Node Object

The node object contains two configuration properties:

- Software versions and neuron ID of the GMS side. This is a configuration property of 14 bytes in size, with the "read only" attribute, i.e., it is only readable if the device is online. This is a UCPT type 1.

This configuration property is structured as follows:


- Software version of LonWorks® side with four bytes in the format MAJOR, MINOR, REVISION_HIGH, REVISION_LOW
 - Neuron ID of GMS side with six bytes
 - Software version of GMS side with four bytes in the format MAJOR, MINOR, REVISION_HIGH, REVISION_LOW
- A MaxSendTime as a heartbeat for the transmission of the nvoGmsState variables, in the range of 0 to 6,553.4 seconds. The application periodically checks for any change to this configuration parameter. Default value: 60 seconds

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GMS Gateway – Data Collectors Object

This object represents the data from up to 8 GMS data collectors. Each data collector has up to 64 bits of status information, so that the corresponding network variable has the format SNVT_state_64.

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Network Variables for the Data Collectors Object

There are only eight output variables, structured as an array. The data point from input 1 of a data collector is displayed in nvoDataCollect[X].bit0.

Any change in the value of a variable is transmitted immediately.

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Configuration Properties of the Data Collectors Object

There are different configuration properties for this object which are responsible for the transmission of output variables:

- **MaxSendTime for the Data Collectors Object**
This configuration parameter from SCPT type 49 enables the periodic transmission of all output variables, thus saving the parametrisation of a number of individual heartbeats. When a value is unequal to 0, the variable with the next higher index is transmitted, after the configured time runs out. After index 7 has been sent, the next transmission is the variable with index 0.
Default value of the parameter: 0 seconds

NOTE

This transmission does not take into account whether a variable is actually bound. For this reason, the transmission cycle time is always the same. If the time is parameterised, for instance, with 1 second, all variables are then transmitted every 8 seconds. This is only "visible" with bound variables.

- **MaxSendTime for each Output Variable of the Data Collectors Object**
If MaxSendTime is not to be used for the overall object, there is an individual heartbeat parameter available, for each output variable, also of SCPT type 49, with a value configurable between 0 and 6,553.4 seconds.
Default value for these parameters: 0 seconds

NOTE

To make use of these individual times, the MaxSendTime for the data collectors object must be 0 seconds.

Alarms Supported in the Node Object

When testing this object, the status "in alarm" and "comm failure" are supported. Here, "in alarm" is active if at any position in the eight output variables of the object, the value 1 is set.

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Area Control Units Object

To represent the data for the area control units object, there are two versions available:

- Version 1: digital display of "Alarm"
- Version 2: digital and analogue display of "Alarm" and "Pressure"

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
Area Control Units Object Version 1 (alarm, digital display)

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The largest area control unit can accommodate up to 5 gases, so that a maximum of 10 digital error messages, in the form of high and low alarms are given. For this reason, the variables of format SNVT_state have 16 individual bits.

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Network Variables for the Area Control Units Object

There are only 50 output variables, structured as an array. The low alarm for medium 1 of an area control unit is shown in `nvoAcuStates[X].bit0`, and the high alarm for medium 1 in `nvoAcuStates[X].bit1`. Thus, for an area control unit with 5 gases, max. `nvoAcuStates[X].bit9` can be assigned to this variable corresponding to the high alarm for medium 5. Any change in the values of a variable is transmitted immediately.

Configuration Properties of the Area Control Units Object

There are different configuration properties for this object as well, which are responsible for the transmission of output variables:

- **MaxSendTime for the Area Control Units Object**
This configuration parameter of SCPT type 49 enables periodic transmission of all output variables of the object. When a value is unequal to 0, the variable with the next higher index is transmitted, after the time runs out. After index 49 has been sent, the next transmission is the variable with index 0.
Default value of the parameter: 0 seconds

NOTE

This transmission does not take into account whether a variable is actually bound. For this reason, the transmission cycle time is always the same. If the time is parameterised, for instance, with 1 second, then the transmission of all variables is done every 50 seconds. This is only "visible" with bound variables.

- **MaxSendTime for each Output Variable of the Area Control Units Object**
If `MaxSendTime` is not to be used for the overall object, there is an individual heartbeat parameter available, for each output variable, also of SCPT type 49, with a value configurable between 0 and 6,553.4 seconds.
Default value for these parameters: 0 seconds

NOTE

To make use of these individual times, the `MaxSendTime` for the area control units object must be 0 seconds.

Alarms Supported for the Node Object

When testing this object, the status "in alarm" and "comm failure" are supported. Here, "in alarm" is active if at any position in the 50 output variables of the object the value 1 is set.

"Comm failure" is active if a node can no longer be reached in GMS.

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
Area Control Units Object Version 2 (alarm and pressure, digital and analogue display)

This object is used to display the digital data of up to 50 area control units.

The largest area control unit can accommodate up to 5 gases, so that a maximum of 10 digital error messages, in the form of high and low alarms along with 5 analogue pressure values are given. For this reason, the variables of format UNVT_contain are in a user-defined format.

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```
typedef struct
{
    unsigned char vk5_state[2];
    SNVT_press gas_pressure[5];
}S_UNVT_CONTAIN;
```

Network Variables for the Area Control Units Object

There are only 50 output variables, structured as an array. The structure of the container variables UNVT_contain correspond to the following:

```
typedef struct
{
    unsigned char vk5_state[2];
    SNVT_press gas_pressure[5];
}S_UNVT_CONTAIN;
```

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The bits of vk5_state are defined as follows:

Vk5_state[0]	Wertigkeit/ Value	Information	Vk5_state[1]	Wertigkeit/ Value	Information
Bit 7	2 ⁷ =128	Gas 1 Low	Bit 7	2 ⁷ =128	Gas 5 Low
Bit 6	2 ⁶ =64	Gas 1 High	Bit 6	2 ⁶ =64	Gas 5 High
Bit 5	2 ⁵ =32	Gas 2 Low	Bit 5	2 ⁵ =32	Frei/ Free
Bit 4	2 ⁴ =16	Gas 2 High	Bit 4	2 ⁴ =16	Frei/ Free
Bit 3	2 ³ =8	Gas 3 Low	Bit 3	2 ³ =8	Frei/ Free
Bit 2	2 ² =4	Gas 3 High	Bit 2	2 ² =4	Frei/ Free
Bit 1	2 ¹ =2	Gas 4 Low	Bit 1	2 ¹ =2	Frei/ Free
Bit 0	2 ⁰ =1	Gas 4 High	Bit 0	2 ⁰ =1	Frei/ Free

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Any change to the value of the digital status bits is transmitted immediately.

Analogue pressure values are recorded in an array of type SNVT_press of 5 bytes in the container variables. The manner in which changes to pressure values are processed in the gateway is determined by the configuration properties defined below.

An unused pressure value is shown as \$7fff.

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Configuration Properties of the Area Control Units Object

There are different configuration properties for this object as well, which are responsible for the transmission of output variables:

- **MaxSendTime for the Area Control Units Object**
This configuration parameter of SCPT type 49 enables periodic transmission of all output variables of the object. When a value is unequal to 0, the variable with the next higher index is transmitted after the time runs out. After index 49 has been sent, the next transmission is the variable with index 0.
Default value of the parameter: 0 seconds

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- **MinDeltaLevel for the Area Control Units Object**
This configuration parameter of SCPT type 88 enables transmission when an analogue value changes only after deviation by at least MinDeltaLevel, in % of the current value. When a value is unequal to 0, only the container value is sent if the difference between the current analogue value and the last analogue value sent is larger than MinDeltaLevel. MinDeltaLevel has no effect on transmission through the MaxSendTime for the area control units object. MinDeltaLevel can be set between 0 % and 100 %.
Default value of the parameter: 0 %
- **MaxSendTime for each Output Variable of the Area Control Units Object**
If MaxSendTime is not to be used for the overall object, there is an individual heartbeat parameter available, for each output variable, also of SCPT type 49, with a value configurable between 0 and 6,553.4 seconds.
Default value for these parameters: 0 seconds

**NOTE**

To make use of these individual times, the MaxSendTime for the area control units object must be 0 seconds.

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- **MinSendTime for each Output Variable of the Area Control Units Object**
There is a MinSendTime available for each output variable. It is not possible to update the output variables any quicker than the value configured here. The parameters MinDeltaLevel and MinSendTime thus ensure that a changing pressure value does not lead to a higher rate of transmission of the output variables, and thus a higher load on the bus.
MinSendTime can be configured at a value between 0 and 6,553.4 seconds.
Default value for these parameters: 0 seconds

**NOTE**

To make use of these individual times, the MaxSendTime for the area control units object must be 0 seconds.

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
GMS Gateway VoluCount™ – Volucount Object

This object is used to represent the data of up to 25 Volucounts.

The object is structured as follows:

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```
typedef struct
{
    unsigned char vc_state[1];
    unsigned long   gas_flow[5];
    SNVT_count_f gas_count[5];
}UNVT_CONTAIN_VC;
```

Network variables of the Volucount Object

There are only 25 output variables, structured as an array.

The container variable UNVT_contain_vc has the following structure:

```
typedef struct
{
    unsigned char vc_state[1];
    unsigned long   gas_flow[5];
    SNVT_count_f gas_count[5];
}UNVT_CONTAIN_VC;
```

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The bits of the nviResetStates are defined as follows:

nviResetState	Reset of	nviResetState	Reset of
Bit 1	Channel 1	Bit 9	-
Bit 2	Channel 2	Bit 10	-
Bit 3	Channel 3	Bit 11	-
Bit 4	Channel 4	Bit 12	-
Bit 5	Channel 5	Bit 13	-
Bit 6	-	Bit 14	-
Bit 7	-	Bit 15	-
Bit 8	-	Bit 16	-

If the value of the digital status is changed, it is transmitted immediately.

The gas flow and gas count values are transmitted in an array of size 5. The configuration properties, which are described in the following, determine how the changes of gas flow values are processed in the gateway. During transmission, gas count values are firmly linked with the status and gas flow values.

Configuration Properties of the Volucount Object

For this object, different configuration properties exist, which are responsible for the transmission of the output variables.

– MaxSendTime for the Volucount Object

This configuration parameter from SCPT type 49 enables the periodic transmission of all output variables of the object. When a value is unequal to 0, the variable with the next higher index is transmitted, after the configured time runs out. After index 49 has been sent, the next transmission is the variable with index 0.
Default value of the parameter: 0 seconds

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NOTE

This transmission does not take into account whether a variable is actually bound. For this reason, the transmission cycle time is always the same. If the time is parameterised, for instance, with 1 second, all variables are then transmitted every 25 seconds. This is only "visible" with bound variables.

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– **MinDeltaLevel for the Volucount Object**

This configuration parameter of SCPT type 88 enables transmission when an analogue value changes only after deviation by at least MinDeltaLevel, in % of the current value. When a value is unequal to 0, the container value is only sent if the difference between the current analogue value and the last analogue value sent is larger than MinDeltaLevel. MinDeltaLevel has no effect on transmission through the MaxSendTime for the Volucount object. MinDeltaLevel can be set between 0 % and 100 %. Default value of the parameter: 0 %

– **MaxSendTime for each Output Variable of the Volucount Object**

If MaxSendTime is not to be used for the overall object, there is an individual heartbeat parameter available, for each output variable, also of SCPT type 49, with a value configurable between 0 and 6,553.4 seconds.

Default value for these parameters: 0 seconds

NOTE

To make use of these individual times, the MaxSendTime for the Volucount object must be 0 seconds.

– **MinSendTime for each output variable of the Volucount Object**

If MaxSendTime is not to be used for the overall object, a parameter for an individual minimum send interval for the transmission of changed analogue values of SCPT type 50 is available for each output variable.

MinSendTime can be used to parameterise a value between 0 and 6,553.4 seconds.

Default value for these parameters: 0 seconds

NOTE

To make use of these individual times, the MaxSendTime for the Volucount object must be 0 seconds.

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When testing this object, the status "in alarm" and "comm failure" are supported.

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

Operation

General

The GMS gateway resp. the GMS gateway VoluCount™ is the data interface between the Dräger Alarm System and a LonWorks®-capable external system, e.g. a building management system.

The device has 5 lights for status display, used to show the status of the data network.


- 1 Green 24 VDC light: the device is ready for operation.

Service LEDs

- 2 Both controller boards include a service LED which is controlled by the Echelon® firmware. The normal state is ON for 0.5 sec after power up, then permanently OFF.

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

Status DUO LED on the GMS Side

- 1 The status DUO LED signals the following operational states:
- Steady green: Normal operating state. No errors, target state
 - Flashing green at 1Hz: At least one node in the GMS cannot be reached through communications
 - Flashing red at 1Hz: No communication possible with the LonWorks[®] controller
 - Flashing yellow five times at 1Hz: Reaction to a LonWorks[®] Wink

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Status DUO LED on the LonWorks[®] Side

- 2 The status DUO LED signals the following operational states:
- Steady green: Normal operating state. No errors, target state
 - Flashing red: No communication possible with the GMS controller
 - Flashing yellow five times at 1Hz: Reaction to "Identify GMS"

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DANGER

Do not disinfect!

Do not use cleaning agents which release chlorine or oxygen!

Do not use solvents such as benzine or ether! This may damage the surface.

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Maintenance Intervals

Daily

Check operation of the system by visual inspection.

Regularly, at least once every six months

Functional test and visual inspection of the entire system carried out by experts.

Contracting a service agreement with DrägerService is recommended.

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Fault – Cause – Remedy

Fault	Cause	Remedy
Green 24 VDC light does not light up.	Power supply interrupted.	Have power supply repaired by experts.
Status LED on GMS side permanently lit red or yellow.	Hardware defective	Have repaired by experts.
Status LED on LonWorks [®] side permanently lit red or yellow.	Hardware defective	Have repaired by experts.
Status LED on LonWorks [®] side flashing green.	Hardware defective	Have repaired by experts.

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Technical Data

Ambient conditions

Operating temperature	+15 to +40 °C
Relative humidity	max. 95 %

Operating values

Power supply	24 VAC
Current consumption	60 mA

LON¹ data

transmission

Transmission rate	78 kBaud FT10
Line length	depends on line type
Electrical properties	compliant with Echelon- FTT 10A transceiver speci- fications

Software version information

GMS gateway

Generation of application

SW version on GMS side: GMG_0101 (G 51 299)	01 00 00 01
SW version on LonWorks [®] side: GMS Gateway A V1_01 (G 51 297)	01:00:00:01
GMS Gateway AP V2_00 (G 51296)	02:00:00:00

GMS gateway VoluCount™

SW version on GMS side: VCG_0100 (G 41 893)	01:00:00:00
SW version on LonWorks [®] side: VC gateway V1_00 (G 41 892)	01:00:00:00

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Order List

Description	Order number
Shield terminal	G 41 784
Power supply, 30 watt/top hat rail	G 50 728
GM network system cable	G 50 750

¹ LON: Local Operation Network Fa. Echelon

¹ LON: Local Operation Network, manufactured by Echelon

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Index

Area Control Units Object	14
Version 1	14
Alarmer supported for the node object	15
Configuration properties	15
Network variables	15
Version 2	16
Alarmer supported for the node object	18
Configuration properties	17
Network variables	16
Data Collectors Object	12
Alarmer supported in the node object	13
Configuration properties	13
Network variables	12
Fault – Cause – Remedy	25
For Your Safety and that of Your Patients	3
Installation	6
Intended Use	4
LonWorks® Interface	8
Maintenance Intervals	24
Visual check	24
Visual check and display test	24
Node Object	9
Configuration properties	11
Network variables, manufacturer-specific area	11
Network variables, optional	10
Network variables, prescribed	10
Operation	22
General	22
Service LEDs	22
Status DUO LED on the GMS side	23
Status DUO LED on the LonWorks® side	23
Order List	26
Purify	24
Technical Data	26
Testing and Commissioning	7
Volucount Object	19
Alarmer supported for the node object	21
Configuration properties	20
Network variables	19
What's What	5

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Directive 93/42/EEC
concerning Medical Devices

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© Dräger Medical GmbH
Drägerwerk AG & Co. KGaA