

Instructions for use

Infinity Acute Care System



WARNING

To properly use accessories, read and comply with all relevant instructions for use.

Monitoring accessories

Typographical conventions

- 1 Consecutive numbers indicate steps of action, with the numbering restarting with "1" for each new sequence of actions.
- Bullet points indicate individual actions or different options for action.
- Dashes indicate the listing of data, options, or objects.
- (A) Letters in parentheses refer to elements in the related illustration.
- A Letters in illustrations denote elements referred to in the text.

Trademarks

- Infinity®
- Acute Care System
- Hemo4®
- Hemo2®
- MPod®
- MCable™
- MonoLead®

are trademarks of Dräger.

- rainbow®
 - Masimo®
 - Masimo® SET® (Signal Extraction Technology)
- are registered trademarks of Masimo Corporation.

- Nellcor™
- OxiMax™
- OxiBand™
- Durasensor™

are trademarks of Nellcor Puritan Bennett, LLC.

Definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the equipment or other property.

NOTE

A NOTE provides additional information intended to avoid inconveniences during operation.

This page intentionally left blank.

Contents

For your safety and that of your patients . . .	7
List of accessories	8
List of accessories	11
Infinity Acute Care System accessories overview	13
General IACS accessories	16
IACS / M540 options	18
ECG monitoring – reusable leadsets	19
ECG monitoring – single-patient use leadsets .	20
ECG monitoring – MonoLead	21
ECG monitoring – neonatal	23
Invasive blood pressure (IBP) with Hemo2 Pod or Hemo4 pod	24
IBP Monitoring with MPod – QuadHemo or MCable – Dual Hemo	25
Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors	26
Non-invasive blood pressure (NIBP) monitoring – neonatal	28
Pulse oximetry monitoring – Masimo rainbow SET (rainbow and M-LNCS)	29
Pulse oximetry monitoring – Masimo rainbow SET (LNCS)	30
Pulse oximetry monitoring – Masimo SET (LNCS)	31
Pulse oximetry monitoring – Masimo SET (M- LNCS)	33
Pulse oximetry monitoring – Nellcor OxiMax . .	34
Cardiac (C.O.) output monitoring	35
CO ₂ monitoring	36
Temperature monitoring	37
MCable – Analog/Sync	40

This page intentionally left blank.

For your safety and that of your patients

Strictly follow these instructions for use	8
General precaution	8
Training	8
Sterile accessories	8
Accessories	9
Connected devices	9
Installing accessories	9

Strictly follow these instructions for use Training

WARNING

To maintain patient safety, adhere to all **WARNINGS** and **CAUTIONS** on equipment labels or appearing in the instructions for use for the Infinity Acute Care System (IACS) applications.

Training for users is available via the Dräger organization responsible (see www.draeger.com).

Sterile accessories

CAUTION

Read all cleaning instructions (for example, originating from the disinfectant manufacturer and the hospital) carefully *before* cleaning the device. Refer to the chapter entitled "Cleaning" for device-specific cleaning instructions in the Infinity Acute Care System M540 instructions for use.

CAUTION

Risk of medical device failure and of patient injury.

Do not use sterile-packaged accessories if the packaging has been opened, is damaged or if there are other signs of non-sterility. Single-use accessories must not be reused, reprocessed, or resterilized.

General precaution

CAUTION

Even accessories designed to be reused have a limited service life. Handling and reprocessing can increase wear and markedly shorten service life (e.g., disinfectant residues can attack the material more intensely during autoclaving). If signs of wear become visible, such as cracks, deformation, discoloration, peeling, etc., affected accessories must be replaced.

List of accessories

Some articles are not available worldwide because they are not approved in all countries of the world.

Technical documentation is available on request.

Accessories

WARNING

Risk due to incompatible accessories

Dräger has tested only the compatibility of accessories listed in the current list of accessories. If other, incompatible accessories are used, there is a risk of patient injury due to medical device failure.

Dräger recommends that the medical device is only used together with accessories listed in the current list of accessories.

Connected devices

WARNING

Risk of electric shock and of device malfunction

Any connected devices or device combinations not complying with the requirements mentioned in these instructions for use can compromise the correct functioning of the medical device and lead to an electric shock. Before operating the medical device, strictly comply with the instructions for use of all connected devices or device combinations.

Installing accessories

CAUTION

Risk of device failure

Install accessories to the basic device in accordance with the instructions for use of the basic device. Make sure that there is a safe connection to the basic device.

Strictly observe instructions for use.

CAUTION

Risk of incorrect use

The accessories are not available individually. Only one copy of the instructions for use is included in the bulk package and must therefore be kept in a location accessible for users.

This page intentionally left blank.

List of accessories


Infinity Acute Care System accessories overview	13	Non-invasive blood pressure (NIBP) monitoring – neonatal	28
General IACS accessories	16	NIBP hose	28
System components	16	Single-patient use neonatal NIBP cuffs	28
Cables and mouse	16	Pulse oximetry monitoring – Masimo rainbow SET (rainbow and M–LNCS)	29
Keyboards	16	Masimo rainbow SET cables and sensors	29
Miscellaneous hardware	16	Masimo M–LNCS sensors	29
Mounts for workstation	17	Pulse oximetry monitoring – Masimo rainbow SET (LNCS)	30
Mounts for power supply	17	Masimo SpO ₂ cables	30
Mounts for the Cockpit and M500/M540	17	Masimo single-patient use LNCS sensors	30
IACS / M540 options	18	Masimo SpO ₂ accessories (not shown)	30
ECG monitoring – reusable leadsets	19	Pulse oximetry monitoring – Masimo SET (LNCS)	31
ECG monitoring – single-patient use leadsets	20	Masimo SpO ₂ cables	31
ECG monitoring – MonoLead	21	Masimo disposable sensors	31
ECG accessories	21	Masimo reusable sensors	32
MonoLead	21	Masimo SpO ₂ accessories (not shown)	32
ECG monitoring – neonatal	23	Pulse oximetry monitoring – Masimo SET (M–LNCS)	33
Neonatal leadsets and accessories	23	Masimo SET cables and sensors	33
Invasive blood pressure (IBP) with Hemo2 Pod or Hemo4 pod	24	Masimo M–LNCS sensors	33
IBP pods and cables	24	Pulse oximetry monitoring – Nellcor OxiMax	34
Transducer adapter cables	24	Nellcor SpO ₂ cables	34
IBP Monitoring with MPod – QuadHemo or MCable – Dual Hemo	25	Nellcor single-patient use sensors	34
IBP cables	25	Nellcor reusable sensor	34
Transducer adapter cables	25	Cardiac (C.O.) output monitoring	35
Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors	26	C.O. cables, catheters, and thermistors	35
NIBP hose	26	Cardiac output kits (not shown)	35
Reusable adult/pediatric NIBP cuffs	26	CO₂ monitoring	36
Single-patient use adult/pediatric NIBP cuffs	27	CO ₂ cables	36
		Disposable airway adapters	36
		Reusable airway adapters	36




























List of accessories

Additional accessories (not shown)	36
Temperature monitoring	37
Temperature cables	37
Skin probes	37
Temperature cables	38
General purpose probes	38
MCable – Analog/Sync	40
Analog/sync cables	40

Infinity Acute Care System accessories overview

The following diagram is an overview of the accessories available with the Infinity Acute Care System.



CO ₂		CO ₂ & Analog/Sync		Analog & Sync		Temp			
 <p>+</p>		 <p>+</p>		 <p>+</p>		 <p>+</p>			
 <p>Pediatric Adult</p>		 <p>CO₂ Analog/Sync</p>				 <p>+</p>			
ECG – reusable		ECG – MonoLead		ECG – single-use		ECG – neo		NIBP	
								<p>Adult/Child Neonatal</p> 	
		 <p>3</p>				<p>+</p> 		<p>+</p>	
		 <p>5</p>							
		 <p>4 6</p>							

087

This chapter lists the accessories available with the Infinity Acute Care System (IACS).

The chapter is divided into the following sections:

- "General IACS accessories" on page 16
- "IACS / M540 options" on page 18
- "ECG monitoring – reusable leadsets" on page 19
- "ECG monitoring – single-patient use leadsets" on page 20
- "ECG monitoring – MonoLead" on page 21
- "ECG monitoring – neonatal" on page 23
- "Invasive blood pressure (IBP) with Hemo2 Pod or Hemo4 pod" on page 24
- "IBP Monitoring with MPod – QuadHemo or MCable – Dual Hemo" on page 25
- "Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors" on page 26
- "Non-invasive blood pressure (NIBP) monitoring – neonatal" on page 28
- "Pulse oximetry monitoring – Masimo rainbow SET (LNCS)" on page 30
- "Pulse oximetry monitoring – Masimo rainbow SET (rainbow and M-LNCS)" on page 29
- "Pulse oximetry monitoring – Nellcor OxiMax" on page 34
- "Cardiac (C.O.) output monitoring" on page 35
- "CO₂ monitoring" on page 36
- "Temperature monitoring" on page 37
- "MCable – Analog/Sync" on page 40

On the following pages, each parameter section includes a diagram depicting parameter-specific accessories, along with one or more tables displaying accessory names and part numbers. The letters displayed in the diagrams correspond to the letters in the tables. Because some accessories are too numerous to list in the diagram, they do not have a letter associated with them.

Part numbers are subject to change without prior notice. Contact your Dräger representative for the most current part numbers.

Failure to use Dräger-approved accessories may compromise the proper functioning of the IACS.

General IACS accessories

The following table lists the system components, power cables, system cables, and keyboards for the IACS.

System components

The following parts are components in the IACS. They are not available for purchase individually. Contact your Dräger representative for additional information.

Part	Part No.
Infinity M500	MS20407
Infinity M540	MS20401
Infinity P2500	MS22277
Infinity PS50	8416800

Cables and mouse

Part	Part No.
Mouse	MS23352
Mouse pad	MS24519
System cable, 1.5 m (4.9 ft)	MS20135
System cable, 2.5 m (4.9 ft)	MK08707
System cable, 3.5 m (14.7 ft)	MK08708
System cable, 4.5 m (14.7 ft)	MS20614
IACS Y-cable for PS50	MS20345
System cable, 2.2 m (7.2 ft) <i>(Note that this cable is only to be used in combination with Zeus IE)</i>	MK04408
Ethernet patch cable, 4.9 m (16.1 ft)	4726399
Ethernet cables, 2.4 m (7.9 ft)	4726381
Ethernet cables, 1.2 m (3.9 ft)	4726373
Nurse call cable	8417370
Export protocol cable	MS22948

Keyboards

Part	Part No.
Standard keyboard (Russian)	MS24520
Washable keyboard (Czech)	MS22286
Washable keyboard (Denmark)	MS20790
Washable keyboard (English)	MS22288
Washable keyboard (French)	MS22295
Washable keyboard (German)	MS22287
Washable keyboard (Greek)	MS23692
Washable keyboard (Hungarian)	MS22282
Washable keyboard (Italian)	MS22283
Washable keyboard (Norwegian)	MS22284
Washable keyboard (Portuguese)	MS22291
Washable keyboard (Spanish)	MS22294
Washable keyboard (Swedish/Finnish)	MS22285
Washable keyboard (Turkish)	MS22292

Miscellaneous hardware

Part	Part No.
Transport dock and clamp	MS28144
Cable hook	MS26297
SpO ₂ MCable mount	MS26266
Rainbow SET MCable mount	MS28576
Rolling stand for M540	RS-0025-05C

Mounts for workstation

Part	Part No.
Mount for workstation	MS17139
Mount bracket rail GCX	MS15934

Mounts for power supply

Part	Part No.
Wall mount adapter CPS/IDS power supply	4720061
Rail mount for CPS/IDS power supply	4720095
Wall bracket CPS/IDS power supply	4720129
Power supply rollstand bracket CPS/IDS	4720103

Mounts for the Cockpit and M500/M540

Part	Part No.
Dräger lift arm, front rail, VESA, +300	G96548
Dräger double arm, front rail, VESA	G96715

NOTE

All hardware from GCX can only be purchased directly from GCX corporation or one of its distributors.


IACS / M540 options

The following table lists the available IACS options-- Check with your local sales representative for options that might be available for free during certain ongoing promotions.

The menu selections on the M540 or the C700/C500 will not be available until the corresponding options have been enabled.

Option	M540 standalone	C700	C500	Part No.
M540 wireless	Yes	No	No	MS16266
Full arrhythmia	Yes	Yes	Yes	MS22225
Multi IBP	Yes	Yes	Yes	MS20504
12-lead ECG	Yes	Yes	Yes	MS20508
12 channels	No	No	Yes	MS20505
16 channels	No	Yes	Yes	MS20506
Web enabled layouts	No	Yes	Yes	MS20511
Web enabled tabs	No	Yes	Yes	MS20512
View editor	No	Yes	Yes	MS20515
Physiological calculations	No	Yes	Yes	MS20516

ECG monitoring – reusable leadsets



ECG accessories

	Part	Part No.
A	ECG extension cable 3/5/6, 2 m (6.5 ft) ¹⁾	MS16256
B	ECG "4-lead" spacer, 10 pieces	MS22832

NOTE


¹⁾ The extension cable works with MonoLeads

Traditional leadsets

	Part	Part No.
C	ECG grabber 3-lead, IEC1 (Eu- rope), 1 m (3.3 ft)	MS16159
	ECG grabber 3-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MS16231
D	ECG grabber 5-lead, IEC1 (Eu- rope), 1 m (3.3 ft)	MS16158
	ECG grabber 5-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MS16546
E	ECG grabber 6-lead, IEC1 (Eu- rope), 1 m (3.3 ft)	MS16157
	ECG grabber 6-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MS16547
F	ECG grabber 4-lead, IEC1 (Eu- rope), 1 m (3.3 ft) ¹⁾	MS22365
	ECG grabber 4-lead, IEC2 (AHA/US) 1 m (3.3 ft) ¹⁾	MS22364

074

ECG monitoring – single-patient use leadsets



ECG accessories

	Part	Part No.
A	ECG extension cable 3/5/6, 2 m (6.5 ft) ¹⁾	MS16256
B	ECG "4-lead" spacer, 10 pieces	MS22832

NOTE


¹⁾ The extension cable works with MonoLeads, traditional leads, and single-patient use leads. The extension cable is not compatible with ECG 4-lead grabber sets.

Traditional leadsets

	Part	Part No.
C	ECG grabber 3-lead, IEC1 (Eu- rope), 1 m (3.3 ft)	MP00875
	ECG grabber 3-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP00877
D	ECG grabber 5-lead, IEC1 (Eu- rope), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP00879
	ECG grabber 5-lead, IEC2 (AHA/US), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP00881
E	ECG grabber 6-lead, IEC1 (Eu- rope), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP03122
	ECG grabber 6-lead, IEC2 (AHA/US), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP03123
F	ECG grabber 4-lead, IEC1 (Eu- rope), 1 m (3.3 ft)	MP03120
	ECG grabber 4-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03121

483

ECG monitoring – MonoLead



ECG accessories

	Part	Part No.
A	ECG extension cable 3/5/6, 2 m (6.5 ft) ¹⁾	MS16256
B	ECG "4-lead" spacer, 10 pieces	MS22832

NOTE

¹⁾ The extension cable works with certain MonoLead sets, traditional leads, and single-patient use leads (see below for more information). The extension cable is not compatible with ECG MonoLead 4-lead.


MonoLead

	Part	Part No.
C	ECG MonoLead 3-lead, IEC1 (Europe), 2 m (6.5 ft)	MS14555
	ECG MonoLead 3-lead, IEC2 (AHA/US), 2 m (6.5 ft)	MS14556
	ECG MonoLead 3-lead, IEC1 (Europe), 4.1 m (13.5 ft) ¹⁾	MS28561
	ECG MonoLead 3-lead, IEC2 (AHA/US), 4.1 m (13.5 ft) ¹⁾	MS28557
D	ECG MonoLead 5-lead, IEC1 (Europe), 2.5 m (8.2 ft)	MS14559
	ECG MonoLead 5-lead, IEC2 (AHA/US) 2.5 m (8.2 ft)	MS14560
	ECG MonoLead 5-lead, IEC1 (Europe), 5.3 m (17.4 ft) ¹⁾	MS28559
	ECG MonoLead 5-lead, IEC2 (AHA/US) 5.3 m (17.4 ft) ¹⁾	MS28558
¹⁾ Note: The indicated accessories must NOT be used in combination with the extension cable MS16256		

List of accessories

E	ECG MonoLead 6-lead, IEC1 (Europe), 3 m (10 ft)	MS14683
	ECG MonoLead 6-lead, IEC2 (AHA/US), 3 m (10 ft)	MS14682
	ECG MonoLead 6-lead, IEC1 (Europe), 5.7 m (18.7 ft) ¹⁾	MS17185
	ECG MonoLead 6-lead, IEC2 (AHA/US), 5.7 m (18.7 ft) ¹⁾	MS28560
F	ECG MonoLead 4-lead, IEC1 (Europe), 1.7 m (5.6 ft)	MS14693
	ECG MonoLead 4-lead, IEC2 (AHA/US), 1.7 m (5.6 ft)	MS14694
	ECG MonoLead 4-lead, IEC1 (Europe), 3.7 m (12.1 ft)	MS20476
	ECG MonoLead 4-lead, IEC2 (AHA/US), 3.7 m (12.1 ft)	MS20475
¹⁾ Note: The indicated accessories must NOT be used in combination with the extension cable MS16256		

ECG monitoring – neonatal




Neonatal leadsets and accessories

WARNING

To prevent burns during electrosurgery (ESU), do not use the ECG neonatal 3-lead adapter cable with the ECG extension cable (MS16256).

	Part	Part No.
A	ECG neonatal 3-lead adapter cable, 2 m (6.5 ft), single-pin	MS25951
B	ECG electrode KittyKat 300 pieces	5195024
C	ECG 4-lead spacer, 10 pieces	MS22832

Invasive blood pressure (IBP) with Hemo2 Pod or Hemo4 pod



IBP pods and cables

	Part	Part No.
A	M540 Hemo Connect cable, 1.5 m (4.9 ft) ¹⁾	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft) ¹⁾	MS23898
B	Hemo2 pod (<i>not shown</i> – no longer available)	4319435
	Hemo4 pod (no longer available)	4315961

NOTE


¹⁾ The M540 Hemo Connect cable is used with the Hemo2 pod/Hemo4 pod.

Transducer adapter cables

	Part	Part No.
C	Hemo pod adapter, 10-pin	3375958
D	Hemo pod adapter, SensoNor, 7-pin	MS23598
E	Hemo pod adapter, Baxter-Edwards	5196980
F	Hemo pod adapter, Abbott-Medex	5196998
G	Hemo pod adapter, Becton/Dickinson	3375941

0183

IBP Monitoring with MPod – QuadHemo or MCable – Dual Hemo



IBP cables

	Part	Part No.
A	MCable – Dual Hemo	MS20783
B	M540 - Hemo Connect cable, 1.5 m (4.9 ft) ¹⁾	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft) ¹⁾	MS23898
C	MPod – QuadHemo	MS20725


NOTE

¹⁾ The M540 Hemo Connect cable is used with the MPod – QuadHemo (not necessary for the MCable – Dual Hemo).

Transducer adapter cables

	Part	Part No.
D	Abbott-Medex (Transtar)	MS22535
E	Baxter-Edwards	MS22147
F	Becton/Dickinson	MS22148
G	Utah Medical	MS22534
H	Dräger 7-pin	MS22533
I	Dräger 10-pin	MS22532

Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors



NIBP hose

The following single-patient use and reusable cuffs and hose use Dräger connectors.

	Part	Part No.
A	Connection tube, adult/pediatric 3.7 m (12.15 ft)	MP00953

Reusable adult/pediatric NIBP cuffs

NOTE

All cuffs are latex-free.

	Part	Part No.
B	Infant arm cuff, XXS 8 to 13 cm (3 to 5 in)	MP00911
	Child arm cuff, XS 12 to 19 cm (5 to 7 in)	MP00912
	Small adult arm cuff, S 17 to 29 cm (6 to 10 in)	MP00913
	Adult arm cuff, M 23 to 33 cm (9 to 13 in)	MP00915
	Adult arm cuff, M+ 23 to 43 cm (9 to 17 in)	MP00916
	Adult arm cuff, M++ 23 to 53 cm (9 to 21 in)	MP00917
	Large adult arm cuff, L 31 to 40 cm (12 to 16 in)	MP00918
	Large adult arm cuff, L+ 31 to 55 cm (12 to 22 in)	MP00919
	Large adult thigh cuff, XL 38 to 50 cm (15 to 20 in)	MP00921


082

Single-patient use adult/pediatric NIBP cuffs

See the illustration on page 26.

	Part	Part No.
B	Infant arm cuff, XXS, 8 to 13 cm (3 to 5 in) 10 pieces	MP00924
	Child arm cuff, XS, 12 to 19 cm (5 to 7 in) 10 pieces	MP00925
	Small adult arm cuff, S, 17 to 29 cm (6 to 11 in) 10 pieces	MP00926
	Adult arm cuff, M, 23 to 33 cm (9 to 13 in) 10 pieces	MP00928
	Adult arm cuff, M+, 23 to 43 cm (9 to 17 in) 10 pieces	MP00929
	Adult arm cuff, M++, 23 to 53 cm (9 to 21 in) 10 pieces	MP00930
	Adult arm cuff, L, 31 to 40 cm (12 to 16 in) 10 pieces	MP00931
	Adult arm cuff, L+, 31 to 55 cm (12 to 22 in) 10 pieces	MP00932
	Adult thigh cuff, XL, 38 to 50 cm (15 to 20 in) 10 pieces	MP00934

Non-invasive blood pressure (NIBP) monitoring – neonatal



NIBP hose

The following single-patient use and reusable cuffs and hoses use universal connectors.


	Part	Part No.
A	Neonatal NIBP hose, 2.4 m (7.87 ft)	2870298

Single-patient use neonatal NIBP cuffs

	Part (arm circumference)	Part No.
B	Neonatal cuff #1 3.1 to 5.7 cm (1.2 to 2.2 in)	2870181
	Neonatal cuff #2 4.3 to 8.0 cm (1.7 to 3.1 in)	2870199
	Neonatal cuff #3 5.8 to 10.9 cm (2.3 to 4.3 in)	2870207
	Neonatal cuff #4 7.1 to 13.1 cm (2.8 to 5.1 in)	2870215
	Neonatal cuff #5 8.3 to 15.0 cm (3.2 to 6 in)	2870173

0180

Pulse oximetry monitoring – Masimo rainbow SET (rainbow and M–LNCS)



Masimo rainbow SET cables and sensors


All Masimo rainbow SET sensors and cables, red LNC cables, M–LNCS sensors can only be purchased directly from Masimo or one of its distributors.

	Part	Part No.
A	Infinity MCable – Masimo rainbow SET	MS27003

Masimo M–LNCS sensors

	Part
B	rainbow direct connect sensor
C	rainbow intermediate cable
D	Single-patient use and reusable rainbow or M–LNCS sensors

Pulse oximetry monitoring – Masimo rainbow SET (LNCS)



Masimo SpO2 cables

	Part	Part No.
A	Infinity MCable – Masimo rainbow SET	MS27003
B	Red LNC intermediate cable, LNCS (all Masimo rainbow SET sensors and cables, red LNC cables, M–LNCS sensors can only be purchased directly from Masimo or one of its distributors).	

The following tables list the compatible sensors available for purchase from Dräger.


Masimo single-patient use LNCS sensors

	Part	Part No.
C	Adtx, adult, 20 pieces	MP00790
	Inf/Inf-3, infant, 20 pieces	MP00791
	Neo/Neo-3, neonatal, 20 pieces	MP00792
	NeoPt/NeoPt-3, preterm, 20 pieces	MP00794
	Pdtx, pediatric, 20 pieces	MP00793
D	DCI adult digit	MP00796
	DCIP pediatric digit	MP00795
	TC-I adult ear	MP00788
	YI adult/pediatric/infant/neonatal multi-site	MP00789
	TF-I adult forehead	MP00799

Masimo SpO2 accessories (not shown)

Part	Part No.
Adhesive SpO2 square, 12 pieces	7496875
Clothing clip, 5 pieces	7496883
Petite wrap-adult/pedi/neo	7496909
Posey wrap, 12 pieces	7496933
Sensor tape replacement, 100 pieces	7496925

Pulse oximetry monitoring – Masimo SET (LNCS)



Masimo SpO2 cables

	Part	Part No.
A	Infinity MCable – Masimo SET	MS20667
B	Masimo intermediate cable, LNCS 1.2 m (3.9 ft)	MP00800
	Masimo intermediate cable, LNCS 3 m (10 ft)	MP00797

The following tables list the compatible sensors available for purchase from Dräger.

Masimo disposable sensors

	Part	Part No.
C	Atdx, adult, 20 pieces	MP00790
	Inf/Inf-3, infant, 20 pieces	MP00791
	Neo/Neo-3, neonatal, 20 pieces	MP00792
	NeoPt/NeoPt-3, preterm 20 pieces	MP00794
	Pdtx, pediatric, 20 pieces	MP00793

Masimo reusable sensors


See illustration on page 31.

	Part	Part No.
D	DCI adult digit	MP00796
	DCIP pediatric digit	MP00795
	TC-I adult ear	MP00788
	YI adult/pediatric/ infant/neonatal multi-site	MP00789
	TF-I adult forehead	MP00799

Masimo SpO₂ accessories (not shown)

Part	Part No.
Adhesive SpO ₂ square, 12 pieces	7496875
Clothing clip, 5 pieces	7496883
Petite wrap-adult/pedi/neo	7496909
Posey wrap, 12 pieces	7496933
Sensor tape replacement, 100 pieces	7496925

Pulse oximetry monitoring – Masimo SET (M–LNCS)



Masimo SET cables and sensors

All Masimo SET M–LNCS sensors and M–LNC cables can only be purchased directly from Masimo or one of its distributors.


	Part	Part No.
A	Infinity MCable – Masimo SET	MS20667
B	M–LNC intermediate cable	

Masimo M–LNCS sensors

	Part
C, D	Single-patient use and reusable M–LNCS sensors

497

Pulse oximetry monitoring – Nellcor OxiMax



Nellcor SpO2 cables

	Part	Part No.
A	Infinity MCable – Nellcor OxiMax	MS20668
B	Nellcor intermediate cable, 1.2 m (3.9 ft)	MP00748
	Nellcor intermediate cable, 3 m (9.8 ft)	MS20979

Nellcor single-patient use sensors


The following tables list the compatible sensors available for purchase from Dräger.

	Part	Part No.
C	OxiMax MaxA, 46 cm (18 in), 24 pieces	MX50065
	OxiMax MaxAL, 91 cm (36 in), 24 pieces	MX50071
	OxiMax MaxI, infant, 24 pieces	MX50067
	OxiMax MaxN, neo/adult, 24 pieces	MX50068
	OxiMax MaxP, pediatric, 24 pieces	MX50066
	Nellcor, OxiMax Max-R	MX50069
	Nellcor DURA-Y SpO2 sensor, D-YS	MX01004
	Nellcor SpO2 sensor, OxiMax Max-Fast, 24 pieces	MX50070
	OxiBand Oxi-P/I, ped/inf, 100 pieces	2170221
	Nellcor Oxiband wrap A/N, 100 pieces	8201012

Nellcor reusable sensor

	Part	Part No.
D	Durasensor DS100A, adult, finger	7262764

Cardiac (C.O.) output monitoring



C.O. cables, catheters, and thermistors

	Part	Part No.
A	M540 Hemo Connect cable, 1.5 m (4.9 ft)	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft)	MS23898
B	MPOD – QuadHemo	MS20725
	Hemo2 pod (not shown)	4319435
	Hemo4 pod (not shown)	4315961
C	Intermediate cable, 1 m (3.3 ft)	3368458
D	Intermediate cable to catheter	8419160
E	Thermistor cable, Ohmeda	8420077
F	Thermistor T-piece, 25 pieces	5741975
G	Thermistor cable, Baxter-Edwards	8539983


Cardiac output kits (not shown)

	Part	Part No.
	Cardiac output kit, Baxter ¹⁾	7265155
	Cardiac output kit, Ohmeda ¹⁾	7265163

NOTE

¹⁾ Kit includes intermediate cable, catheter cable, and thermistor cable.

CO2 monitoring



CO2 cables

	Part	Part No.
A	MCable - Mainstream CO2	6871950
B	Y-cable CO2 and analog/sync ¹⁾	MS22259

NOTE

¹⁾ The Y-cable allows simultaneous use of the MCable - Mainstream CO2 and MCable – Analog/Sync.

Disposable airway adapters

	Part	Part No.
C	Airway adapter main, adult	MP01062
	Airway adapter main, pediatric	MP01063

Reusable airway adapters


	Part	Part No.
C	Airway adapter main, adult	6870279
	Airway adapter main, pediatric	6870280

Additional accessories (not shown)

Part	Part No.
Calibration set	8412710
CO2 test gas cylinder	6850435
Parking holder for CO2 sensor	8416352
Test filter	6870281

073

Temperature monitoring



Temperature cables

	Part	Part No.
A	Adapter cable 7-pin connector for temperature probe with 1/4-inch connector	5198333
	Dual temperature Y cable, 7-pin	5592154

Skin probes

WARNING

Temperature skin probes are not for use in electro-surgery.

Adult reusable


	Part	Part No.
B	Temperature skin probe, 1.5 m (4.9 ft)	4329822
	Temperature skin probe, 3 m (9.8 ft)	5204669

Adult single-patient use

	Part	Part No.
C	Temperature skin probe, adult, 1.6 m (5.2 ft)	MP00995
	Temperature skin probe, adult, 3 m (9.8 ft)	MP00996

Pediatric/neonatal single-patient use

	Part	Part No.
C	Temperature skin probe, pediatric, 1.6 m (5.2 ft)	MP00997
	Temperature skin probe, pediatric, 3 m (9.8 ft)	MP00998



Temperature cables

	Part	Part No.
A	Adapter cable 7-pin connector for temperature probe with 1/4-inch connector	5198333
	Dual temperature Y cable, 7-pin	5592154

General purpose probes

WARNING
Temperature probe covers contain latex.

NOTE
Cover internally placed reusable temperature sensors with temperature probe sheaths.

Adult reusable general purpose probes

	Part	Part No.
B	Temperature probe, adult 1.5 m (4.9 ft)	4329889
	Temperature probe, adult 3 m (9.8 ft)	5204644
C	Core probe covers, 20 cm, 10 pieces (contains latex)	7014616

Pediatric reusable general purpose probes

	Part	Part No.
B	Temperature probe, pediatric, 1.5 m (4.9 ft)	4329848
	Temperature probe, pediatric, 3 m (9.8 ft)	5204651

Adult single-patient use general purpose probes

D	Temperature probe, adult, 1.6 m (5.2 ft)	MP00991
	Temperature probe, adult, 3 m (9.8 ft)	MP00992

Pediatric single-patient use general purpose probes

D	Temperature probe, pediatric, 1.6 m (5.2 ft)	MP00993
	Temperature probe, pediatric, 3 m (9.8 ft)	MP00994

NOTE

1) Temperature core probes can be used during electrosurgery and in esophageal and rectal applications.

NOTE

Temperature probes with 7-pin connectors are compatible with Hemo4 pod, Hemo2 pod, MPod – QuadHemo (see page 24).


MCable – Analog/Sync

Analog/sync cables

	Part	Part No.
A	Y-cable CO ₂ and analog/sync ¹⁾	MS22259
B	MCable – Analog/Sync	MS20662
C	Analog output cable, 5 m (16.4 ft)	4314618
	QRS sync output cable, 3 m (9.8 ft)	4314667

NOTE

¹⁾ The y-cable allows simultaneous use of the MCable - Mainstream CO₂ and MCable – Analog/Sync.



071

This page intentionally left blank.

These instructions for use only apply to
Infinity Acute Care System – Monitoring Accessories SW VG2

with the Serial No.:

If no Serial No. has been filled in by Dräger,
these Instructions for Use are provided for general information only and are not intended for use with any specific machine or unit.

This document is provided for customer information only, and will not be updated or exchanged without customer request.



Directive 93/42/ECC
concerning medical products

 **Manufacturer**
Dräger Medical GmbH
 Moislinger Allee 53 – 55
D-23558 Lübeck
Germany
 +49 451 8 82-0
FAX +49 451 8 82-20 80
 <http://www.draeger.com>

Distributed in USA by
Draeger Medical, Inc.
 3135 Quarry Road
Telford, PA 18969-1042
U.S.A.
 (215) 721-5400
(800) 4DRAGER
(800 437-2437)
FAX (215) 723-5935
 <http://www.draeger.com>

MS29258 – RI03 – en
© Dräger Medical GmbH
Edition: 2 – 2013-02
(Edition: 1 – 2012-06)



Dräger reserves the right to make modifications
to the equipment without prior notice