Gas Management Monitor 3G - Basic/ 6G - Basic Monitor 3G - AMS/ 6G - AMS

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WARNING

To properly use this medical device, read and comply with these Instructions for Use.

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Safety Information Definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the medical device or other property.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation.

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Content

For Your Safety and that of Your Patients	3
Intended Use	4
Acceptance and Handover to the User	5
Operation	6
Normal operation	7
Alarm	7
Function test	8
Communication error	9
Remote display	9
Combination possibilities	11
Cleaning and Disinfecting	12
Maintenance Intervals	13
Fault - Cause - Remedy	14
Technical Data	15
System Overview	16

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For Your Safety and that of Your Patients

Strictly follow these Instructions for Use

WARNING

Any use of the medical device requires full understanding and strict observation of all sections of these Instructions for Use. The medical device must only be used for the purpose specified under "Intended Use" on page 4.

Strictly observe all WARNING and CAUTION statements throughout these Instructions for Use and all statements on medical device labels.

Failure to observe these safety information statements constitutes a use of the medical device that is inconsistent with its intended use.

Maintenance

WARNING

The medical device must be inspected and serviced regularly by professionals who possess the required qualifications due to their training and their experience. Repair of the medical device must also be performed by trained personnel with additional product-specific DrägerService training.

Dräger recommends that a service contract is obtained with DrägerService and that all repairs are performed by DrägerService. Dräger further recommends that only authentic Dräger repair parts are used for maintenance. If the above are not complied with, the correct functioning of the medical device may be compromised. See chapter "Maintenance Intervals".

Not for use in areas of explosion hazard

WARNING

This medical device is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

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Intended Use

These Instructions for Use apply to the use of the Monitor 3G - Basic, Monitor 6G - Basic, Monitor 3G - AMS and Monitor 6G - AMS emergency alarm panels.

Monitor 3G - Basic and Monitor 6G - Basic are used as clinical emergency alarm panels in accordance with DIN EN ISO 7396-1 for monitoring the pressure of medical gases at the area stop valve. In particular, Monitor 3G - Basic and Monitor 6G - Basic are used for retrofits of area control units without independent electronic control units.

Monitor 3G - Basic processes signals for up to 3 gases, Monitor 6G - Basic for up to 6 gases.

The pressure is monitored by means of external switch gauges or appropriate isolated contacts.

The alarms generated can be transmitted to the Dräger Alarm Management System via a network connection. In addition, Monitor 3G - Basic and Monitor 6G - Basic feature a relay driver output with one multi-messaging signal contact per gas type.

Monitor 3G - AMS and Monitor 6G - AMS are used as additional informative alarm panels for the clinical emergency signals from the area control units.

In this case, the pressure is monitored exclusively in the area control unit.

The alarms are transmitted via the network connection. Both direct connection to an individual area control unit and integration into the Dräger Alarm Management System are possible.

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Acceptance and Handover to the User

The Monitor 3G - Basic, Monitor 6G - Basic, Monitor 3G - AMS and Monitor 6G - AMS emergency alarm panels are part of the Dräger Gas Management.

NOTE

The system may only be operated after acceptance by specialized staff.

Observe national regulations!

In the European Economic Area (EEA), for example, DIN EN ISO 7396-1 applies.

After completion of the installation or after maintenance work, the specialized expert will perform a test and acceptance program.

This test determines:

- whether the safety requirements designed to protect patients and staff have been met

and

- whether the performance characteristics of the central gas supply system and the Dräger Gas Management are met.

The test results are to be documented in written form.

Following these acceptance tests, the operational pipeline system is handed over to the operating organization together with the documentation and the operating personnel are instructed on how to operate the system.

A formal record must be kept of the handover procedure.

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
Operation

Emergency alarm panels and subsidiary alarm panels monitor the operating status of central gas supply systems for medical gases and vacuum, evaluate the data and indicate the operating status of the monitored system visually and/or acoustically.

In the "AMS" enhancement, the operating status is transmitted via the Dräger Alarm Management System and can be simultaneously displayed on operation alarm panels.

Emergency alarm panels are available in 2 versions:

- Monitor 3G - Basic (G 50726)
- Monitor 6G - Basic (G 50727)

Subsidiary alarm panels are available in 2 versions:

- Monitor 3G - AMS (G 50736)
- Monitor 6G - AMS (G 50737)

The Monitor 3G - Basic and Monitor 6G - Basic emergency alarm panels are installed in a separate location from the area control unit and feature 3 or 6 display units respectively.

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

The display units show the status of the monitored gas. They consist of 4 elements:

- Status indicator lights
- Label panel
- Audio alarm
- Key for muting the audio alarm and for the function test.

The Monitor 3G - AMS subsidiary alarm panel is an emergency alarm panel that is installed in a separate location from the corresponding area control unit. It can indicate alarms for up to 3 gases.

The Monitor 6G - AMS subsidiary alarm panel is an emergency alarm panel that is installed in a separate location from the corresponding area control unit. It can indicate alarms for up to 6 gases.

Normal operation

If the pipeline pressure measured in the area control unit is within a tolerance of

- +/-20 % of the nominal pressure for O₂, N₂O, CO₂, AIR and AIR 800 or
- higher than -0.34 bar relative for VAC,

1 the green light will be on.

Alarm

If any pipeline pressure is outside the tolerance, an alarm is activated.

- If the pressure is greater than
- +20 % of the nominal pressure for O₂, N₂O, CO₂, AIR and AIR 800 or
 - lower than -0.4 bar relative for VAC,
- 1 the green light turns off,
 2 the red light flashes.
 • An audio alarm is activated.

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

If the pressure is less than

- 20 % of the nominal pressure for O₂, N₂O, CO₂, AIR and AIR 800:

- 1 the green light turns off,
 - 3 the red light flashes.
- An audio alarm is activated.

The audio alarm can be muted for approx. 10 to 15 minutes:

- 4 Press the »●« key,
- the audio alarm is muted.
 - The red light changes from flashing to steady light.

If the reason for the alarm activation still remains after approx. 10 to 15 minutes:

- the audio alarm will reactivate and
- the red light will start flashing again.

If the signal line from the switch gauge to the emergency alarm panel is interrupted, an alarm will also be indicated.

When the cause of the alarm has been eliminated:

- 1 the red light turns off,
 - 2 the green light lights up again.
- The audio alarm stops.

Function test

A function test can be carried out on an emergency alarm panel when no alarms are active.



- Press and hold the »●« keys once each.

Monitor 3G - Basic and Monitor 3G - AMS:

- 3 All red lights are on in addition to the already activated green lights.

Monitor 6G - Basic and Monitor 6G - AMS:

- 3 All red lights on one side are on in addition to the already activated green lights.
- The audio alarm is activated.

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
Any relays connected for potential-free transmission are not activated during a function test. Thus, a Building Management System, which is possibly connected, will not be affected by the test.

Communication error

In the "AMS" enhancement, the communication is permanently monitored via the Dräger Alarm Management System. If the emergency alarm panel or the subsidiary alarm panel detects a communication error with the connected operating alarm panels, that operating status is indicated on an emergency alarm panel:

- The green light flashes.

Alarms will continue to be indicated by red lights on the Monitor 3G - Basic and the Monitor 6G - Basic.

Remote display

If the system is installed with the Dräger Alarm Management System or connected relays, the status of the emergency alarm panel is transmitted as a message to downstream units.

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Combination possibilities

Combinations of	Monitor 3G - Basic	Monitor 6G - Basic	Monitor 3G - AMS	Monitor 6G - AMS
ACU 3 gases with no display electronics	Direct connection of switch gauge from one ACU	Direct connection of switch gauge from 2 ACUs	Not possible	Not possible
ACU 3 gases with display electronics	Connection via relay module (G 18989) for one multi-messaging signal per gas type from one ACU	Connection via one relay module each (G 18989) for one multi-messaging signal per gas type from 2 ACUs	Connection via network line from Monitor G from one ACU	Connection via network line from Monitor G from 2 ACUs
ACU 5 gases equipped with 1 to 3 gases	Connection via relay module (G 41742) for high and low alarm per gas type from one ACU	Connection via one relay module each (G 41742) for high and low alarm per gas type from 2 ACUs	Connection via Alarm System Interface (G 41743) from one ACU	Connection via Alarm System Interface (G 41743) from 2 ACUs
ACU 5 gases equipped with 4 to 5 gases	Not all gases can be displayed	Connection via two relay modules (G 41742) for high and low alarm per gas type from one ACU	Not all gases can be displayed	Connection via Alarm System Interface (G 41743) from one ACU
Integration into the Alarm Management System	Transmission of signals to the Alarm Management System via optional transceiver (G 50719)	Transmission of signals to the Alarm Management System via optional transceivers (2 x G 50719)	Transceiver for Alarm Management System already installed	Transceiver for Alarm Management System already installed
Housing	AP set 130 (G 50744) UP set 130 (G 50740)	AP set 250 (G 50745) UP set 250 (G 50741)	AP set 130 (G 50744) UP set 130 (G 50740)	AP set 250 (G 50745) UP set 250 (G 50741)

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Cleaning and Disinfecting

CAUTION

Do not use solvents such as benzine or ether!

- Wipe surface with a disposable cloth.

NOTE

Wipe disinfection only!

Follow the directions for use issued by the disinfectant manufacturer.

Surface disinfectants are to be used for disinfection. For material compatibility use disinfectants based on

- aldehydes,
- alcohols,
- quaternary ammonium compounds.

Damage to materials may occur if disinfectants are used which are based on

- phenols,
- halogen-releasing compounds,
- strong organic acids,
- oxygen-releasing compounds.

Users in the Federal Republic of Germany are advised to use the disinfectants found in the current DGHM list (DGHM: German Society for Hygiene and Microbiology).

The DGHM list also gives the active ingredients of each disinfectant.

The following surface disinfectants are recommended for material compatibility:

Surface disinfectant	Manufacturer
Terralin	Schülke & Mayr GmbH
Sagromed Spray	Schülke & Mayr GmbH
Incidin Perfect	Ecolab GmbH & Co. OHG
Sokrena	Bode Chemie GmbH

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Maintenance Intervals

Daily

Perform a visual inspection and function test of the entire system to ensure that it is ready for operation.

Regularly, at least every six months

Function test and visual inspection of complete system by trained service personnel.

It is recommended to conclude a service agreement with DrägerService.

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Fault - Cause - Remedy

Fault	Cause	Remedy
None of the lamps light up in an emergency	Power failure	Have the power supply repaired by experts.
Gauge pressure correct but emergency alarm panel shows error message (pressure too high or pressure too low)	Sensor contact defective	Repair by experts.
Green light flashing	Communication problems between emergency alarm panel and area control unit	Repair by experts.

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24 V DC, 24 V AC, or 12 V DC

Power consumption

Max. power during operation
for Monitor 3G - Basic and
Monitor 3G - AMS 5 W
for Monitor 6G - Basic and
Monitor 6G - AMS 10 W

Relays for potential-free transmission

Relay drivers for the external connection of relays are available.	Relay driver data: 1 driver per gas type, voltage = operating voltage, rectified, max. driver current 50 mA, not short circuit-proof
--------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

Data transfer

AMS

Transfer rate:	78.13 kbps
Conductor length:	depends on type of conductor
Electrical specifications:	acc. to Echelon-FTT 10A Transceiver Specification

Noise emission
acc. to DIN EN ISO 7396-1 > 57 dB

UMDNS code 18-046

Universal Medical Device
Nomenclature System

Visual display units acc. to DIN EN ISO 7396-1

Systemübericht

Benennung	Sachnummer
Monitor 3G - Basic	G 50 726
Monitor 6G - Basic	G 50 727
Monitor 3G - AMS	G 50 736
Monitor 6G - AMS	G 50 737
Zubehör:	
HW- / UP-Set 130	G 50 740
HW- / UP-Set 250	G 50 741
AP-Set 130	G 50 744
AP-Set 250	G 50 745
Rangierverteiler / Platine	G 50 732
Systemkabel GM-Netz	G 50 750
Rangierblock 3-pol. / Hutschiene	G 50 723
Transceiver FTT10A	G 50 719
Relaismodul	G 18 989

System Overview

Description	Order no.
Monitor 3G - Basic	G 50 726
Monitor 6G - Basic	G 50 727
Monitor 3G - AMS	G 50 736
Monitor 6G - AMS	G 50 737
Accessories:	
HW / UP set 130	G 50 740
HW / UP set 250	G 50 741
AP set 130	G 50 744
AP set 250	G 50 745
Junction block / board	G 50 732
System cable for GM network	G 50 750
Junction block, 3-pole / top hat rail	G 50 723
Transceiver FTT 10A	G 50 719
Relay module	G 18 989

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Directive 93/42/EEC
concerning Medical Devices

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Edition: 9 – 2015-01

(Edition: 1 – 1999-03)

Dräger reserves the right to
make modifications to the
device without prior notice.

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© Dräger Medical GmbH
Dräger Medical GmbH
Drägerwerk AG & Co. KGaA