Medical Air Compressor

de	Ergänzung, Seite 2	ja	補足説明書, ページ 103
en	Supplement, page 8	ko	부록, 페이지 108
fr	Complément, page 14		
es	Anexo, página 20		
nl	Aanvullingsdocument, pagina 26		
it	Foglio informativo, pagina 32		
pt	Documento suplementar, página 38		
sv	Tilläggsinformation, sida 44		
tr	Ek bilgi sayfası, sayfa 50		
pl	Suplement, strona 56		
cs	Dodatek, strana 62		
ro	Fișă de completare, pagina 68		
sk	Dodatok, strana 74		
hu	Kiegészítő dokumentum, oldal 80		
ru	Приложение, стр. 86		
bg	Допълнение, Страница 92		
zh	增补, 第 98 页		

WARNING

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1) Dauer des Einschaltstroms: 100-300 ms

Medical Air Compressor

Supplement to the instructions for use

Medical Air Compressor with part number: 9039713, Edition: 2 – 2010-09 and higher

- Keep this supplement with the instructions for use of the medical device.

The supplement updates the information of the instructions for use in the following chapters:

Definition of target groups

For this product, users and experts are defined as target groups.

These target groups must have received instruction in the use of the product and must have the necessary training and knowledge to use, install, reprocess, maintain, or repair the product. The target groups must be able to understand the language of this document.

The product must be used, installed, reprocessed, maintained, or repaired exclusively by defined target groups.

Users

Users are persons who use the product in accordance with its intended use.

Experts

Experts are persons who perform repair or complex maintenance work on the product.

Experts must have the necessary knowledge and experience with complex maintenance work on the product.

For your safety and that of your patients

General safety information

used with this medical device.

The following WARNING and CAUTION statements apply to general operation of the medical device. WARNING and CAUTION statements specific to subsystems or particular features of the medical device appear in the respective sections of these instructions for use or in the instructions for use of another product being

Accessories

WARNING

Risk due to incompatible accessories

Dräger has tested only the compatibility of accessories listed in the current list of accessories or in separate declarations by Dräger. If other incompatible accessories are used, there is a risk of patient injury due to medical device failure.

Dräger recommends that the medical device is only used together with accessories listed in the current list of accessories.

Connection to other devices

Device combinations approved by Dräger meet the requirements of the following standards:

- IEC 60601-1 (3rd edition)
Medical electrical equipment
Part 1: General requirements for safety and essential performance
 - IEC 60601-1-2
Medical electrical equipment
Part 1-2: General requirements for safety and essential performance
Collateral standard: Electromagnetic compatibility; Requirements and tests
 - IEC 60601-1-8
Medical electrical equipment
Part 1-8: General requirements for safety and essential performance
Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
- IEC 60601-1 (2nd edition)
Medical electrical equipment
Part 1: General requirements for safety
 - IEC 60601-1-1
Medical electrical equipment
Part 1-1: General requirements for safety
Collateral standard: Safety requirements for medical electrical systems
 - IEC 60601-1-2
Medical electrical equipment
Part 1-2: General requirements for safety
Collateral standard: Electromagnetic compatibility; Requirements and tests
 - IEC 60601-1-4
Medical electrical equipment
Part 1-4: General requirements for safety
Collateral standard: Programmable electrical medical systems
 - IEC 60601-1-8
Medical electrical equipment
Part 1-8: General requirements for safety
Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems

If Dräger devices are connected to other Dräger devices or third-party devices and the resulting combination is not approved by Dräger, the correct functioning of the devices may be compromised. The operating organization must make sure that the device combination conforms to the applicable requirements of the relevant standards.

Strictly observe assembly instructions and instructions for use for each connected device.

Product-specific safety information

WARNING

Risk of malfunction

Unallowed modifications to the medical device lead to malfunctions.

This medical device may not be changed without permission from Dräger.

WARNING

Risk of infection

The personnel responsible may become infected with disease-causing germs.

Before each use, maintenance activity, return, or disposal, reprocess the device and device components in accordance with the instructions for use. Comply with the hospital's hygiene regulations!

WARNING

Risk of electric shock

There are conducting components under the housing cover.

Do not remove the housing cover.

Application

Environment of use

WARNING

Risk of explosion and fire

This medical device is neither approved nor certified for use in areas where oxygen concentrations greater than 25 Vol%, combustible or explosive gas mixtures are likely to occur.

Area of use for medical air compressor

The medical device may only be used in intensive care units and in anesthetic recovery rooms.

Overview

Symbols

Symbol	Explanation
	Storage temperature
	Atmospheric pressure

Symbol	Explanation
	Relative humidity
	Fuse
 nom. 45 kg (100 lbs)	Nominal weight (for information, see chapter "Technical data" on page 13)

Assembly and preparation

Before initial operation

WARNING

Risk of electric shock and device malfunction

If the device is connected to a power socket with incorrect mains voltage or without a protective ground, the user can be injured and the device damaged.

Only connect the power cable to power sockets with a protective ground, see "Technical data".

WARNING

Risk of electric shock

The connection of devices to auxiliary power sockets can lead to an increased leakage current. If the protective ground of one of these devices fails, the leakage current may rise above the permissible values. Only connect with the approval of the respective device manufacturer. Have the leakage current checked by service personnel.

If the permissible value is exceeded, use a mains power socket on a wall instead of the auxiliary power socket of the device.

Operation

Mains power supply

The power supply unit is permanently connected with the medical air compressor.

Switching off the device

The mains plug of the medical device must be accessible so that the device can be disconnected from the mains power supply at any time.

Disconnecting the medical device from the mains power supply:

- Pull out the mains plug.

Technical data

Ambient conditions

Height above sea level

Max. 3000 m (9843 ft)

Operating data

GMDN Code: Global Medical Device Nomenclature – worldwide nomenclature system for medical devices

31253

Use of latex

This device is latex-free.

Pthalates

The biocompatibility of applied parts has been tested in accordance with ISO 10993 for intact skin and breathing gases.

Weight

Nominal weight (device ready for operation, including trolley base 8411520)

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These instructions for use apply only to

Medical Air Compressor

with Serial No.:

If no Serial No. has been filled in by Dräger, these instructions for use are provided for general information only and are not intended for use with any specific device.

These instructions for use are provided for customer information only and will only be updated or exchanged upon customer request.

Directive 93/42/EEC concerning medical devices

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Edition: 2 – 2015-02

(Edition: 1 – 2011-11)

Dräger reserves the right to make modifications to the medical device without prior notice.

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© Dräger Medical GmbH
U.S.A.
Dräger Medical GmbH
Drägerwerk AG & Co. KGaA