Gas Management Monitor L/ Monitor LL/ Monitor LLT

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Instructions for Use

Inhaltsverzeichnis**Contents**

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Dräger Medical GmbH

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the product requires full understanding and strict observation of these instructions. The product is only to be used for the purposes specified here.

Maintenance

The product must be inspected and serviced by experts at regular 6-monthly intervals (and a record kept). Repairs to the product may only be carried out by experts. We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Only authentic Dräger spare parts may be used for maintenance. Observe chapter "Maintenance Intervals".

Not for use in areas of explosion hazard

The product is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Liability for proper function or damage

Liability for the proper function of the product is irrevocably transferred to the owner or operator to the extent that the product is serviced or repaired by personnel not employed or authorized by DrägerService or if the product is used in a manner not conforming to its intended use. Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Drägerwerk Aktiengesellschaft are likewise not modified by the recommendations given above.

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Intended Use

Monitor L, Monitor LL, Monitor LLT

Operating alarm panels for the display of information from connected data collectors and emergency alarm panels for use in medical gas, compressed air and vacuum pipeline systems.

For display of the operating state of a monitored pipeline system in a Dräger data network.

Monitor L : 22 LED displays

Monitor LL : 46 LED displays

Monitor LLT : 48 LED displays and
text message display

The different types of operating alarm panel are selected according to the specific requirements of each project.

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Acceptance and Handover to the User

The Monitor L, Monitor LL and Monitor LLT emergency alarm panels are an integral part of the Dräger Gas Management medical gas, compressed air and vacuum pipeline system.

The unit may not be operated until it has been tested and commissioned by trained and qualified personnel.

Observe the national regulations.
For example EN 737-3 in Europe.

A test and acceptance program must be carried out by specialists upon completion of the installation or maintenance work.

These tests determine:

- whether the safety requirements designed to protect patients and staff have been met

and

- whether the performance characteristics of the medical gas pipeline system and Dräger Gas Management are met.

Written records of the tests must be kept.

Following these acceptance tests, the operational pipeline system is handed over to the user together with the documentation and the operating personnel are instructed on how to operate the system.

A formal record must be kept of the handover procedure.

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Operation

Operating alarm panels show states of the monitored medical gas, compressed air and vacuum pipeline system visually and/or acoustically. In the "Network" enhancement the state is also transferred via the Dräger data network at the same time and can be represented on operating alarm panels.

Operating alarm panels are available in 3 versions:

- Monitor L
- Monitor LL
- Monitor LLT

Monitor L and Monitor LL both indicate the states of connected medical gas pipeline equipment by means of lamps. Monitor LLT can also show the states of connected equipment via a text display in addition to the lamps.


Up to 1,000 pieces of information with instructions for personnel can be shown on this text display.

Monitor L, Monitor LL and Monitor LLT comprise:

- State display lamps
- Label panel
- Audio alarm
- Keys for acknowledging or muting the audio alarm and for the function test
- Additional text display with 2 keys for "scrolling" through the text data (Monitor LLT only)

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Lamp panel


Monitor L and Monitor LL comprise a single lamp panel only. On this lamp panel are displayed operating messages, caution messages with and without audio alarm and warning messages with audio alarm from the connected emergency alarm panels of Monitor G, Monitor 3G, Monitor 6G and the connected data collectors. Each message is explained by the label next to the lamp. Whether a message is an operating message, caution message or warning message is determined at configuration of the Dräger Gas Management System.

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Text panel


In addition to the lamp panel, the Monitor LLT also has a text panel with another 3 lamps. Information on equipment monitored or on current messages can be shown on the Monitor LLT text panel. Up to 1,000 text messages can be displayed in list form on a Monitor LLT. Messages are listed by time of occurrence. The 3 lamps in the text panel show the overall status of the Monitor LLT. The exact function is described below.

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
Status lamps on the LLT

The Monitor LLT has 3 additional status lamps as well as the lamp panel. These status lamps show the overall status.

- 1 If only operating messages are displayed, the green lamp is lit.
- 2 If there is a caution message, the yellow lamp is lit.
- 3 If there is at least one warning message, the red lamp is lit.

The current highest message state can be read from the lamps. The caution and warning lamps flash until the corresponding message is acknowledged by the user.

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Allocation of messages and priorities

The operating, caution and warning messages shown are determined at installation of the Dräger Gas Management System.

Normal operation

In normal operation the green lamps are lit to indicate operation of the equipment connected. All red and yellow lamps are off. No audio alarm sounds.

Operating message

In normal operation of equipment connected via a data collector, a green lamp provides an operating message.

- green lamp lit
- audio alarm off



The message with the highest priority is always indicated acoustically, i.e. the audio alarm sounds while the warning message appears at the same time.

Caution message with no audio alarm

The audio alarm does not sound when a caution message is displayed.

- yellow lamp flashes
- audio alarm off

The message with the highest priority is always indicated acoustically, i.e. the audio alarm sounds while the warning message appears at the same time.

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

Caution message with audio alarm

A caution with audio alarm is activated automatically if a caution message with audio alarm from an emergency alarm panel or data gatherer is displayed.

- yellow lamp flashes
- audio alarm sounds

Warning message with audio alarm

A warning with audio alarm is activated automatically if a warning message with audio alarm from an emergency alarm panel or data gatherer is displayed.

- red lamp flashes
- audio alarm sounds

Automatic silencing of alarm

If a caution with audio alarm is displayed on the operating alarm panel and the cause of that message is remedied, the audio alarm is switched off automatically. The lamp corresponding to that message continues to flash until acknowledged by the operator.

Acknowledgment


Caution messages and warning messages need to be acknowledged:

- Press »« key
- All flashing lamps change to steady light
- The audio alarm stops

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Text messages

Further information in addition to the messages already shown on the lamp panel is given in the text panel on Monitor LLT. Information of lower priority can therefore also be displayed as text messages in the text panel. Text messages can for example show the operation of equipment.


The »↓« and »↑« keys can be used to "scroll" between the different messages. The most recent message is displayed as the default. The »↓« key scrolls to the previous message and the »↑« key scrolls to the next message.

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
Function test

A function test can be run on Monitors L, LL and LLT when no warnings or cautions are active.

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Communications fault

If the operating alarm panel detects a communications fault with a unit, all the lamps allocated to that unit flash green irrespective of their normal colour.

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Cleaning and Disinfecting

Do not use solvents such as benzene or ether.

- Wipe surface with a disposable cloth

Wipe disinfection only

- Follow the directions for use issued by the disinfectant manufacturer.

Use surface disinfectants. For material compatibility use disinfectants based on

- aldehydes,
- alcohols,
- quaternary ammonium compounds.

Damage to materials may occur if disinfectants are used which are based on

- phenols,
- halogen-releasing compounds,
- strong organic acids,
- oxygen-releasing compounds.

Users in the Federal Republic of Germany are advised to use the disinfectants on the current DGHM list (DGHM: German Society for Hygiene and Microbiology).

The DGHM list gives the composition of each disinfectant.

For countries where the DGHM list is not available, we recommend the types of disinfectant given above.

Examples of suitable disinfectants:

- Terralin (Schülke & Mayr, Norderstedt)
- Sagromed Spray (Schülke & Mayr, Norderstedt)
- Incidin Perfect (Henkel KG aA)
- Sokrena (Bode Chemie GmbH&C)

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Maintenance Intervals

Daily

- Check operation of complete system by visual inspection and function test.

Regularly, at least every six months

- Function test and visual inspection of complete system by experts. It is recommended to conclude a service agreement with DrägerService.

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Fault	Cause	Remedy
None of the emergency alarm panel lamps is lit	Power failure	Repair by experts
Display lamps are flashing green	Communications fault in the associated emergency alarm panel or data gatherer	Repair by experts
Function test reveals unexpected result	Audio alarm faulty Lamp faulty Text display faulty	Repair by experts

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Data transfer

LON¹⁾

Transfer rate:	78.13 kbps
Conductor length:	depends on type of conductor
Electrical specifications:	to Echelon-FTT 10A Transceiver Specification

Sound level to EN 737-3	> 57 dB
----------------------------	---------

Code UMDNS 18-046

Universal Medical Device
Nomenclature System

Visual display units	to DIN EN 737-3
----------------------	-----------------

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System Overview

The following components are available for the Gas Management data network:

Designation/Description	Order No.
Monitor CC	G 50 718
Monitor G comprising	
Basic monitor	G 50 720
Front panel Set 1/Monitor G	G 50 701
Front panel Set 2/Monitor G	G 50 702
Front panel Set 3/Monitor G	G 50 703
Monitor G/Network comprising	
Basic monitor/network	G 50 700
Front panel Set 1/Monitor G	G 50 701
Front panel Set 2/Monitor G	G 50 702
Front panel Set 3/Monitor G	G 50 703
Monitor 3G	G 50 726
Monitor 3G/network	G 50 706
Monitor 6G	G 50 727
Monitor 6G/network	G 50 707
Monitor L/network	G 50 712
Monitor LL/network	G 50 711
Monitor LLT/network	G 50 710
Basic data collector/network	G 50 716
Data gatherer upgrade 16D	G 50 717

The components of the system are interconnected via a LON network.
Up to 64 network units may be interconnected in one segment. It is possible to create several segments and connect these up by routers. Some of the products in the Gas Management system network comprise more than one network unit.

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Directive 93/42/EEC
concerning Medical Devices

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90 29 586 - GA 6918.010 de/en
© Dräger Medical GmbH
Edition: 7 - 2015-01
Subject to alteration

As of 2015-08:
Dräger Medical GmbH
changes to
Drägerwerk AG & Co. KGaA