Gas Management Monitor L/Monitor LL/Monitor LLT

Montageanweisung

Installation Instructions

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Contents

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Dräger Medical GmbH

For Your Safety

Strictly follow the installation instructions

Any use of the apparatus during installation requires full understanding and strict observation of these instructions. For operation, the appropriate instructions for use should be followed.

Follow the installation regulations

Installation should only be carried out by VDE qualified experts.

Installation must be in accordance with the relevant DIN/VDE standards or other relevant standards. In Germany, for example, VDE 0100 and VDE 0160 apply.

Follow the static regulations

This is an electronic product. The general regulations on static-sensitive components should be followed when working on the electronics.

Switch off power supply

The product must be voltage-free before and during work.

Not for installation in areas of explosion hazard

Do not install the product in areas of explosion hazard. This product is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Liability for proper function or damage

Liability for the proper function of the product is irrevocably transferred to the owner or operator to the extent that the product is incorrectly installed by personnel not employed or authorized by Dräger. Dräger cannot be held liable for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of Dräger's terms of sale and delivery are likewise not modified by the recommendations given above.

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Intended use

These installation instructions apply to installation of operating alarm panels Monitor L, Monitor LL and Monitor LLT in the Dräger Gas Management system network and the data gatherer including upgrade module and necessary accessories.

Operating alarm panels for the display of information from connected data gatherers and emergency alarm panels for use in medical gas, compressed air and vacuum pipeline systems.

For display of the operating state of a monitored pipeline system in a Dräger data network.

Monitor L : 22 LED displays
Monitor LL : 46 LED displays
Monitor LLT : 48 LED displays and
text message display

The different types of operating alarm panel are selected according to the specific requirements of each project.

Monitor L, Monitor LL and Monitor LLT are components of the Dräger Gas Management system. They are designed for use in an integrated network to display the status of the system monitored. They are not intended for use unless connected to emergency alarm panels or data gatherers.

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
Operating alarm panels

Preliminary installation

Installation of Monitor L, Monitor LL and Monitor LLT operating alarm panels are all described in this chapter. The corresponding installation set should be used in each case.

Hollow walls


- Select installation boxes according to monitor:
 - Monitor L: HW/UP set 250 (G 50 741)
 - Monitor LL: HW/UP set 310 (G 50 742)
 - Monitor LLT: HW/UP set 400 (G 50 743)
- Cutout for appropriate HW/UP set in hollow walls (plasterboard or panel walls with plasterboard).

Fit HW/UP set to hollow walls:

- Select cable bushing,
- Fit cable chafe guard,
- Feed cable through hole in installation box.
- Insert installation box into prepared cutout and align horizontally (arrow indicates top).
- Press installation box firmly into place,
- Bend out retaining clips according to wall thickness.

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- Bend out retaining clips,
- Prevent plaster entering holes at side e.g. by covering with adhesive tape.
- Select cable bushing.
- Fit cable chafe guard.
- Feed cable through hole in installation box.
- Fit plaster guide to installation box.
- Align installation box horizontally and fit flush with plaster, allowing front edge of installation box to project from the raw wall by the thickness of the plaster (about 1 cm).
- Arrow indicates top.

Surface installation

The AP set is used for surface installation.

- Select AP set according to Monitor:
 - Monitor L: AP set 250 (G 50 745)
 - Monitor LL: AP set 310 (G 50 746)
 - Monitor LLT: AP set 400 (G 50 747)
- See back of surface housing for hole position.
- Insert dowel.
- Fit screw centre top and tighten.
- Select cable bushing on back of housing,
- Break out pre-pressed hole.
- Feed cable through hole into housing.
- Close cable entries on building side e.g. using PG screw fittings:
 - Sets 250 and 310: 5xPG11
 - Set 400: 6xPG11
 - for all sets: 1xPG9 and 1xPG13.5 or direct connection to cable duct.
- Hang surface housing in place.
- Tighten screws bottom left and right.
- Leave seal supplied in housing for final assembly.

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
Preparing connecting cable

For operating voltage: either 24V AC or 24V DC

Max. power at 24V DC in operating mode:

- 8,5 W for Monitor L
- 16,5 W for Monitor LL
- 20,5 W for Monitor LLT

Jumpering terminal

The jumpering terminal in the housing base of AP/UP set 400 is not required.

Network cable

GM network system cable (G 50 750); this cable comprises 2x2 pairs communications leads and 2 supply leads for the 24V supply to the emergency and operating alarm panels including shielding.

Alternative cable types should be agreed with the design department for each project.

GM network system cable:

- Pair 1, shielded: La/Lb
- Pair 2, shielded: Reserve
- Pair 3, unshielded: 24V supply

- Remove sheath about 5 cm, separate shield from leads.
- Remove shield except for shielded conductor, cover with shrinkdown plastic tubing, lead through shielded conductor.
- Strip leads.

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
00537002

Connecting network cable and power supply

The power supply to Monitor L, LL and LLT is via external power supplies and uses 24V conductors.

When using 24V DC: Ensure correct polarity.

Monitor L

The cable is connected directly to the Monitor L.

- Strip leads.

Connections:

| Conductors | Monitor L terminal |
|----------------------------|--------------------------------------|
| Network shielded conductor | 1 |
| Network La | 2 |
| Network Lb | 3 |
| DC supply: | |
| Conductor 0V (GND) | 4 |
| Conductor +24V | 5 |
| AC supply: | |
| Conductors 24V ~ | 4 and 5 (independent of polarity) |

- Lay cable in housing so that there is no contact between electronics and shielded conductor.

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

Fitting the unit

UP/HW installation

- Insert unit into UP/HW box.
- On Monitor L and LL, the »E« button is at the bottom right.
- On Monitor LLT, the text display is at the top of the module.
- Fasten using the four M3 screws supplied.
- Cover screws with four plastic caps.

Surface installation

- Replace black rubber seal by light-coloured round seal supplied.
- Insert unit into AP box.
- On Monitor L and LL, the »E« button is at the bottom right.
- On Monitor LLT, the text display is at the top of the module.
- Fasten using the four self-tapping screws supplied.
- Cover screws with four plastic caps.

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
Data gatherer

Preliminary installation

The Dräger Gas Management data gatherer is designed for installation on DIN mounting rails.

Recommended mounting location:
AP housing 360

AP housing complete: G 50 748
Housing dimensions:
WxLxH = 160x360x100 mm

- Fit housing to wall at place of use.

Prepare connecting cable

Operating voltage

either 24V AC or 24V DC

Wiring

- Only one lead should be allocated to each terminal of the data gatherer and upgrade units.

Cable type

- In principle, only flexible cables should be used. Very rigid cables put excessive mechanical stress on the terminals.
- To avoid wiring errors, only cables with colour-coded individual leads should be used.

Data and network cable

- E.g. J-Y (ST) Y number of leads x 2 x 0.8 mm

Number of leads according to volume of information to be processed, e.g. 16 x 2 x 0.8 mm

Power cable

- Remove sheath about 5 cm, strip leads.

When using the Dräger Gas Management system cable (G 50 725):
continue as described in chapter "Network cable".

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- Select step sequence for power cable and print out.
- Cut out step sequence with a blade to corresponding cable diameter.
- Slide step sequence over cable and press into starting position in housing.

Network cable

LON interface cable: multi-lead cable with stranded pairs, shielded.

E.g. GM network system cable (G 50 750); this cable comprises 2x2 pairs communications leads and 2 supply leads for the 24V supply to the emergency and operating alarm panels including shielding.

Alternatively:

– E.g. J-Y (ST) Y 2x2x0.8 mm


Other alternative cable types should be agreed with the design department for each project.

GM network system cable:

- Pair 1, shielded: La/Lb
 - Pair 2, shielded: Reserve
 - Pair 3, unshielded: 24V supply
- Remove sheath about 5 cm, separate shield from leads.
 - Remove shield except for shielded conductor, cover with shrinkdown plastic tubing, lead through shielded conductor.
 - Strip leads.
 - Select step sequence for network cable and print out.
 - Cut out step sequence with a blade to corresponding cable diameter.
 - Slide step sequence over cable and press into starting position in housing.

Data lines

- Remove sheath about 5 cm. Strip leads.
- Select step sequence for data lines and print out.

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- Cut out step sequence with a blade to corresponding cable diameter.
- Slide step sequence over cable and press into starting position in housing.

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
Final installation


Basic data gatherer

The basic data gatherer is an electronic module for installation on a mounting rail.

- Clamp basic data gatherer mounting rail module (G 50 716) on to mounting rail at left side of housing.

Data gatherer upgrade (DGU)

The data gatherer upgrade unit is an electronic module for mounting rail installation for expansion of the data gatherer by additional inputs. A basic data gatherer can be expanded by up to 3 upgrade units.

- Clamp data gatherer upgrade mounting rail module (G 50 717) to mounting rail to right of existing data gatherer.

- Connect existing data gatherer to upgrade unit by means of flat cable supplied.

Potential distributor terminals

The data gatherer has distributor terminals to distribute the ground potential.

Provision is made for connection of one distributor terminal to terminal 17 and one to terminal 18 of the data gatherer or upgrade units.

The distributor terminals should be fitted at the place of use.

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
Optional jumpering board / mounting rail

An optional jumpering board for mounting rail installation is available for jumpering of the 24V power supply and the network cable.

- Clamp jumpering board/top hat rail (G 50 731) on to mounting rail.

Connecting data lines

The data gatherers are designed to evaluate potential-free contacts.

- Connect lines to terminal of corresponding channel (1 to 16) and to associated potential distributor terminal as shown in planning documents.
- Potential distributor terminal at terminal 17 is allocated to channel 1 to 8. Potential distributor terminal at terminal 18 is allocated to channel 9 to 16.

Terminal designation for the data channels is identical on the basic data gatherer and upgrade modules.

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
The cable is connected directly to the basic data gatherer.

- Lead all subsequent leads to terminals 19 to 23 together through the ferrite core supplied.
- Strip leads.
- Connect network cable shield to terminal 21 (shield).
- Connect network cable to terminal 22 (La) and terminal 23 (Lb) (independent of polarity).
- 24V DC supply: Connect lead 0V to terminal 19 (GND) and lead +24V to terminal 20 (24V AC/DC). If polarity is reversed unit will not operate; in that case correct polarity.
- 24V AC supply: Connect supply leads to terminals 19 and 20 (24V AC/DC) (independent of polarity).
- Lay cable in housing so that there is no contact between electronics and shielded conductor.

Close housing

Place cover in position and fasten with screws supplied.

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System overview

Components of the Dräger Gas Management system network:

| Designation/Description | Part No. |
|--|----------|
| Monitor L | G 50 712 |
| Monitor LL | G 50 711 |
| Monitor LLT | G 50 710 |
| Basic data gatherer | G 50 716 |
| Data gatherer upgrade | G 50 717 |
| Accessories: | |
| HW/UP set 250 | G 50 741 |
| HW/UP set 310 | G 50 742 |
| HW/UP set 400 | G 50 743 |
| AP set 250 | G 50 745 |
| AP set 310 | G 50 746 |
| AP set 400 | G 50 747 |
| AP housing 360 | G 50 748 |
| Jumpering distributor / board | G 50 732 |
| Communications cable e.g. GM network system cable | G 50 750 |

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+49 451 8 82-20 80

<http://www.draeger.com>

Directive 93/42/EEC
concerning Medical Devices

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© Dräger Medical GmbH
Drägerwerk AG & Co. KGaA