Infinity Acute Care System

- Workstation Critical Care

- Workstation Neonatal Care

Evita V300

en	Supplement, page 2	hr	Dodatni dokument, stranica 39
da	Supplement, side 6	bg	Допълнение, стр. 43
no	Tilleggsdokument, side 10	sl	Priloga, stran 47
sv	Tilläggsinformation, sida 14	sr	Dodatak, stranica 51
fi	Lisälehti, sivu 18		
ru	Приложение, стр. 22		
pl	Suplement, strona 27		
lt	Priedas, 31 puslapis		
lv	Papildinājums, 35. lappuse		

WARNING

**To properly use this medical device,
read and comply with the instructions
for use and this supplement.**

Supplement to the instructions for use

WARNING

To properly use this medical device, read and comply with the instructions for use and this supplement.

Instructions for use	Part number	Edition
Infinity Acute Care System - Workstation Critical Care, Workstation Neonatal Care	9052077	6 – 2016-01 and higher
Evita V300 SW 2.n	9052995	3 – 2016-04 and higher

- Keep this supplement with the instructions for use. This supplement updates the information of the instructions for use in the following chapters.

For your safety and that of your patients

The following text replaces the sections:

- "Networking" (9052077)
- "Device combinations" (9052995)

Device combinations

This device can be operated in combination with other Dräger devices or with devices from third-party manufacturers. Follow the accompanying documents of the individual devices.

If a device combination is not approved by Dräger, the safety and the functional integrity of the individual devices may be compromised. The operating organization must ensure that the device combination complies with the applicable editions of the relevant standards for medical devices.

Device combinations approved by Dräger meet the requirements of the following standards:

- IEC 60601-1 (electrical safety, mechanical safety, software)
- IEC 60601-1-2 (EMC)
- IEC 60601-1-8 (alarm systems)

General safety information

The following text replaces the sections:

- "Note on EMC/ESD risk for the device function" (9052077)
- "Information on electromagnetic compatibility" (9052995)

Electromagnetic compatibility (EMC)

Medical electrical equipment is subject to special precautionary measures concerning electromagnetic compatibility. During installation and before initial operation, follow the information in section: "EMC declaration" (page 4).

This device can be affected by other electrical devices.

WARNING

Risk due to electrostatic discharge

Malfunctions that endanger the patient may occur if no protective measures against electrostatic discharge are employed in the following situations:

- **When touching the pins of connectors that carry the ESD warning symbol.**
- **When establishing connections with these connectors.**

To prevent malfunctions, observe the following measures and train the relevant personnel:

- **Observe the ESD protective measures. Such measures may include wearing antistatic clothing and shoes, touching a potential equalization pin before and while making the connection, or using electrically insulating and antistatic gloves.**
- **Observe the requirements for the electromagnetic environment. Observe the following section: "Electromagnetic environment" (page 4).**

WARNING

Risk due to electromagnetic disturbance

Wireless communication devices (e.g., cellular phones) and medical electrical equipment (e.g., defibrillators, electrosurgical devices) emit electromagnetic radiation. When such devices are operated too close to this device or its cables, the functional integrity of this device may be compromised by electromagnetic disturbances. As a result, the patient could be put at risk.

Maintain a distance of at least 0.3 m (1.0 ft) between this device and wireless communication devices, to ensure that the essential performance of this device is fulfilled.

Maintain an adequate distance between this device and other medical electrical equipment.

Technical data

The following text replaces the sections:

- "EMC Declaration" (9052077)
- "EMC Declaration" (9052995)

EMC declaration

General information

This device was tested for electromagnetic compatibility using accessories from the list of accessories. Other accessories may only be used if they do not compromise the electromagnetic compatibility. The use of non-compliant accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the device.

This device may be used in the direct vicinity of other devices only if Dräger has approved this device arrangement. If no approval has been given by Dräger, it must be ensured that this device functions correctly in the desired arrangement before use. The instructions for use for the other devices must be followed.

Electromagnetic environment

This device may only be used in environments specified in the following sections of the instructions for use:

- "Environment of use" (9052077)
- "Environment of use" (9052995)

Emissions	Compliance
Radiated emissions	Class A, group 1 (30 MHz to 1 GHz)
Conducted emissions	Class A, group 1 (150 kHz to 30 MHz)

NOTE

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Immunity against	Test level and required electromagnetic environment
Electrostatic discharge (ESD) (IEC 61000-4-2)	Contact discharge: ± 8 kV
	Air discharge: ± 15 kV
Fast transient electrical disturbances (bursts) (IEC 61000-4-4)	Power cable: ± 2 kV
	Longer signal input lines/output lines: ± 1 kV
Impulse voltages (surges) (IEC 61000-4-5)	Voltage, external conductor – external conductor: ± 1 kV
	Voltage, external conductor – protective ground conductor: ± 2 kV
Magnetic fields at mains frequency (IEC 61000-4-8)	50 Hz: 30 A/m
Voltage dips and short interruptions in the supply voltage (IEC 61000-4-11)	Voltage dips of 30 % to 100 %, 8.3 ms to 5 s, different phase angles
Radiated high-frequency disturbances (IEC 61000-4-3)	80 MHz to 2.7 GHz: 3 V/m
Conducted high-frequency disturbances (IEC 61000-4-6)	150 kHz to 80 MHz: 3 V, ISM bands: 6 V
Electromagnetic fields in the vicinity of wireless communication devices	Various frequencies from 385 MHz to 5785 MHz: 9 V/m to 28 V/m

Recommended separation distances from wireless communication devices

To ensure that the functional integrity of this device is maintained, there must be a separation distance of at least 1.0 m (3.3 ft) between this device and wireless communication devices.

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© Drägerwerk AG & Co. KGaA
Drägerwerk AG & Co. KGaA