Druckluft-Steuerschrank / Compressed air control panel

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**Instructions for Use
Series DWSBM / DWSRM / DWSRD**

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elektro stuhr gmbh

For your and your patients safety

Strictly follow the Instructions for Use

Any use of the pressure control panel requires full understanding and strict observation of these instructions for use. The apparatus is only to be used for the purposes specified here.

Maintenance

The cabinet is maintenance-free.
A regular operational test should be carried out together with the maintenance of the plant (maintenance interval of the compressor).
Maintenance of the cabinet only by authorized personnel.
Only original parts are to be used for maintenance.

CAUTION !

Not for use in areas of explosion hazard.
The cabinet is not approved or certified for use in areas of explosion hazard.

Safe connection with other electrical equipment

No electrical connections to equipment which is not listed in these instructions for use

Liability for proper function or damage

The liability for proper function of the cabinet is irrevocably transferred to the owner or operator if the cabinet is improperly serviced or repaired by personnel not authorized by Service or if the cabinet is used in a manner not conforming to the intended use.
The manufacturer is not liable for damage which will be caused by non-compliance with the recommendations given above.
The warranty and liability provisions of the commercial terms of sale and delivery of the manufacturer are not modified by the recommendations given above.

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Intended use

Compressed air control panel for fully-automatic operation and monitoring of compressed air systems with two or three compressors of the type DWSBM, DWSRM or DWSRD.

Operation is assured even if one compressor fails – or if two compressors fail in a system with three compressors.

CAUTION !


The equipment must not be used in the presence of inflammable gases – risk of fire!

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Operation

- 1 Switch on power supply = Set main switch (1)
to "ON"

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
If the system does not work:

- Refer to the schedule "Fault-Cause-Remedy".

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Function of the buttons:

- 1 = DEL: delete (not required)
- 2 = ALT: time/date switching
- 3 = push buttons: select points of the menu
- 4 = OK: switch, store
- 5 = ESC: switch back, cancel
- 6 = * : advice of status
- 7 = red LED: collective malfunction message
- 8 = green LED: operation

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

The Möller control system has the following functions for:

- Advice of status in automatic running
- Advice of disturbance
- Switching base load switching
- System pressure parameter adjustment

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The following information are to read:

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C1: off (compressor is switched off)
 C1: Automatic off (compressor is in automatic function and has not received the switch-on point)
 C1: Automatic on (compressor is in automatic function and is running)
 C1: malfunction (compressor is not operable)
 C3: inactive (only in double operation)

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

Pressure

Pressure: the actual pressure in bar
Switch off pressure: adjusted switch off pressure.

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K1: 0000 h (Betriebstunden)

You can read the following information:

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
Failure emergency service

The system switches to emergency service at falling below the minimal pressure. During the emergency service all pumps which are in automatic function stay operable. In order to except an oversized system load only two compressor will be operable during the emergency service. At achieving the minimal pressure the control switches back to the normal mode. All 3 compressors start to obtain the system pressure.

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

Failure maximum pressure exceeded

This status message will be indicated if the system pressure is 0,5 bar over the set switch-off pressure. The failure can be reset by the button "OK" at falling below the maximum pressure

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
Failure system pressure overview

The message will be released if the system pressure does not achieve the requested switch-off pressure at running compressors within the adjusted time. To reset the message push the button "OK".

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Fault – Cause - Remedy

The following table is intended as a guide for locating and remedying the causes of any faults which may arise.

Only qualified specialists may carry out work for which the equipment must be opened.

Fault	Cause	Remedy
Motor does not start (manual mode)	Power supply not connected to the motor Motor damaged	Main switch "ON" Check main voltage Check circuit-breaker Check motor protection relay Check mode switch
Motor does not start (automatic mode)	Pressure wrongly set on PLC	Set pressostat switching points correctly
Pressure does not build up	2 Phases wrongly connected	Restore correct phase sequence Call for Service
Motor fuses are tripped	Motor current exceeds rated amperage of fuse	Check rated amperage of fuse in accordance with circuit diagram. Check motor data in accordance with circuit. Check motor. Call for Service.
Compressors do not switch between basic load / peak load / reserve load	Failure of Control system	Check control panel Call for Service
Red indicator lights (control voltage)	Failure of Control voltage	Install replacement fuse Call for Service

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Overview Switch cabinet in general

See circuit-diagram, enclosed to the cabinet.

Maintenance

The following operational test should be made together with the maintenance of the plant (maintenance interval of compressor).

Check cutting-in / cutting-off points of the Möller control system.

Check setting of pressure switch for tank pressure alarm.

The corresponding instructions for use must be observed when servicing the compressor.

Disposal

At the end of the equipment's useful service life:

Dispose of the equipment in accordance with national waste disposal regulations,

or

return the equipment to an authorized disposal company for correct disposal.

Further information can be obtained from the local public authorities.

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© elektro stuhr gmbh 2019