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
A steep portion in the inspiratory pressure volume loop may, however, indicate both “healthy” tidal ventilation as well as “unhealthy” intra-tidal recruitment. More recent theoretical and clinical studies have found evidence of alveolar recruitment above the upper inflection point [88, 89]. It was shown that it was possible to shift the upper inflection point farther upwards, suggesting that the UIP does not necessarily signify the beginning of hyperinflation but may instead hold recruitment potential [88]. Both phenomena may occur at the same time in different lung regions, partly masking each other, making therapeutic conclusions difficult. Nonetheless, both experimental and clinical studies demonstrated good results with PEEP settings of LIP + 2 cm H₂O [60, 90].

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Referencias

- [1] Somerson SJ, Sicilia MR. Historical perspectives on the development and use of mechanical ventilation. *J Am Assoc Nurse Anesth* 1992; 60: 83-94.
- [2] Lewandowski K. Extracorporeal membrane oxygenation for severe acute respiratory failure. *Crit Care*. 2000;4(3):156-68. Epub 2000 Apr 12. Review.
- [3] Phua J, Badia JR, Adhikari NK, Friedrich JO, Fowler RA, Singh JM, Scales DC, Stather DR, Li A, Jones A, Gattas DJ, Hallett D, Tomlinson G, Stewart TE, Ferguson ND. Has mortality from acute respiratory distress syndrome decreased over time? A systematic review. *Am J Respir Crit Care Med*. 2009 Feb 1;179(3):220-7. Epub 2008 Nov 14. Review.
- [4] Slutsky AS, Tremblay LN: Multiple System Organ Failure. Is Mechanical Ventilation a Contributing Factor? *Am J Respir Crit Care Med* Vol 157. pp 1721–1725, 1998
- [5] Cook DJ, Walter SD, Cook RJ, et al. Incidence of and risk factors for ventilator-associated pneumonia in critically ill patients. *Ann Intern Med* 1998;129:433–440
- [6] Guidelines for the Management of Adults with Hospital-acquired, Ventilator-associated, and Healthcare-associated Pneumonia. *Am J Respir Crit Care Med* Vol. 171.pp 388-416, 2005
- [7] Dreyfuss D, Saumon G. Ventilator-induced lung injury: lessons from experimental studies. *Am J Respir Crit Care Med*. 1998 Jan;157(1): 294-323.

- [8] Ricard JD, Dreyfuss D, Saumon G. Ventilator-induced lung injury. *European Respiratory Journal*, 1 August 2003, vol. 22, no. Supplement 42, pp. 2S-9S
- [9] Slutsky AS. Lung injury caused by mechanical ventilation. *Chest*. 1999 Jul;116(1 Suppl):9S-15S.
- [10] Chu EK, Whitehead T, Slutsky AS. Effects of cyclic opening and closing at low- and high-volume ventilation on bronchoalveolar lavage cytokines. *Crit Care Med*. 2004 Jan;32(1):168-74.
- [11] Beale R, Grover ER, Smithies M, Bihari D. Acute respiratory distress syndrome (“ARDS”): no more than a severe acute lung injury? *BMJ*. 1993 Nov 20;307(6915):1335-9. Review.
- [12] Marini JJ, Gattinoni L. Ventilatory management of acute respiratory distress syndrome: a consensus of two. *Crit Care Med*. 2004 Jan;32(1):250-5.
- [13] Gattinoni L, Carlesso E, Cadringer P, Valenza F, Vagginelli F, Chiumello D. Physical and biological triggers of ventilator-induced lung injury and its prevention. *Eur Respir J Suppl*. 2003 Nov;47:15s-25s.
- [14] Murphy DB, Cregg N, Tremblay L, Engelberts D, Laffey JG, Slutsky AS, Romaschin A, Kavanagh BP. Adverse Ventilatory Strategy Causes Pulmonary-to-Systemic Translocation of Endotoxin. *Am J Respir Crit Care Med* Vol 162. pp 27-33, 2000
- [15] Terragni PP, Rosboch G, Tealdi A, Corno E, Menaldo E, Davini O, Gandini G, Herrmann P, Mascia L, Quintel M, Slutsky AS, Gattinoni L, Ranieri VM. Tidal hyperinflation during low tidal volume ventilation in acute respiratory distress syndrome. *Am J Respir Crit Care Med*. 2007 Jan 15;175(2):160-6. Epub 2006 Oct 12.

- [16] Gattinoni L, Pesenti A, Avalli L, Rossi F, Bombino M. Pressure-volume curve of total respiratory system in acute respiratory failure. Computed tomographic scan study. *Am Rev Respir Dis.* 1987 Sep;136(3):730-6
- [17] Bidani A, Tzouanakis AE, Cardenas VJ Jr, Zwischenberger JB. Permissive hypercapnia in acute respiratory failure. *JAMA*1994, 272:957-962.
- [18] Hickling KG, Walsh J, Henderson S, Jackson R. Low mortality rate in adult respiratory distress syndrome using low-volume, pressure-limited ventilation with permissive hypercapnia: a prospective study. *Crit Care Med* 1994, 22:1568-1578.
- [19] Putensen C, Mutz N, Putensen-Himmer G, Zinserling J. Spontaneous breathing during ventilatory support improves ventilation-perfusion distributions in patients with acute respiratory distress syndrome. *Am J Respir Crit Care Med* 1999;159:1241-1248.
- [20] Neumann P, Wrigge H, Zinserling J, Hinz J, Maripuu E, Andersson LG, Putensen C, Hedenstierna G. Spontaneous breathing affects the spatial ventilation and perfusion distribution during mechanical ventilatory support. *Crit Care Med.* 2005 May;33(5):1090-5.
- [21] Wrigge H, Zinserling J, Neumann P, Muders T, Magnusson A, Putensen C, Hedenstierna G. Spontaneous breathing with airway pressure release ventilation favors ventilation in dependent lung regions and counters cyclic alveolar collapse in oleic-acid-induced lung injury: a randomized controlled computed tomography trial. *Crit Care.* 2005;9(6):R780-9. Epub 2005 Nov 16.
- [22] Putensen C, Hering R, Muders T, Wrigge H. Assisted breathing is better in acute respiratory failure. *Curr Opin Crit Care.* 2005 Feb;11(1):63-8. Review.

- [23] Priori SG, Aliot E, Blomstrom-Lundqvist C, et al. Task Force on Sudden Cardiac Death of the European Society of Cardiology. *Eur Heart J* 2001; 16: 1374-1450.
- [24] Wysocki M, Cracco C, Teixeira A, Mercat A, Diehl JL, Lefort Y, Derenne JP, Similowski T. Reduced breathing variability as a predictor of unsuccessful patient separation from mechanical ventilation. *Crit Care Med.* 2006 Aug;34(8):2076-83.
- [25] Funk DJ, Graham MR, Girling LG, Thliveris JA, McManus BM, Walker EK, Rector ES, Hillier C, Scott JE, Mutch WA. A comparison of biologically variable ventilation to recruitment maneuvers in a porcine model of acute lung injury. *Respir Res.* 2004 Nov 24;5:22.
- [26] McMullen MC, Girling LG, Graham MR, Mutch WA. Biologically variable ventilation improves oxygenation and respiratory mechanics during one-lung ventilation. *Anesthesiology* 2006; 105:91-7
- [27] Arold SP, Mora R, Lutchen KR, Ingenito EP, Suki B. Variable tidal volume ventilation improves lung mechanics and gas exchange in a rodent model of acute lung injury. *Am J Respir Crit Care Med* 2002; 165:366-71
- [28] Mutch WA, Harms S, Ruth Graham M, Kowalski SE, Girling LG, Lefevre GR. Biologically variable or naturally noisy mechanical ventilation recruits atelectatic lung. *Am J Respir Crit Care Med.* 2000 Jul; 162(1):319-23.
- [29] Gama de Abreu M, Spieth PM, Pelosi P, Carvalho AR, Walter C, Schreiber-Ferstl A, Aikele P, Neykova B, Hübler M, Koch T. Noisy pressure support ventilation: A pilot study on a new assisted ventilation mode in experimental lung injury. *Crit Care Med* 2008; 36:818-27

- [30] Spieth PM, Carvalho AR, Pelosi P, Hoehn C, Meissner C, Kasper M, Hübler M, von Neindorff M, Dassow C, Barrenschee M, Uhlig S, Koch T, de Abreu MG. Variable tidal volumes improve lung protective ventilation strategies in experimental lung injury. *Am J Respir Crit Care Med*. 2009 Apr 15;179(8):684-93. Epub 2009 Jan 16.
- [31] Spieth PM, Carvalho AR, Güldner A, Pelosi P, Kirichuk O, Koch T, de Abreu MG. Effects of different levels of pressure support variability in experimental lung injury. *Anesthesiology*. 2009 Feb;110(2):342-50.
- [32] Suki B, Alencar AM, Sujeer MK, Lutchen KR, Collins JJ, Andrade JS, Ingenito EP, Zapperi S, Stanley HE. Life-support system benefits from noise. *Nature* 1998; 393:127-8
- [33] Mutch WA, Graham MR, Girling LG, Brewster JF. Fractal ventilation enhances respiratory sinus arrhythmia. *Respir Res*. 2005 May 9;6:41.
- [34] Chang HK. Mechanisms of gas transport during ventilation by high-frequency oscillation. *J Appl Physiol*. 1984 Mar;56(3):553-63. Review.
- [35] Ferguson ND, Slutsky AS. Point. High-frequency ventilation is the optimal physiological approach to ventilate ARDS patients. *J Appl Physiol*. 2008 Apr;104(4):1230-1. Epub 2007 Nov 29. Links
- [36] Putensen, Zech, Wrigge. Longterm effects of spontaneous breathing during ventilatory support in patients with acute lung injury. *Am J Respir Crit Care Med* 2001; 164:43-49
- [37] Habashi NM. Other approaches to open-lung ventilation: airway pressure release ventilation. *Crit Care Med*. 2005 Mar;33 (3 Suppl):S228-40. Review.

- [38] Putensen, Zech, Wrigge. Longterm effects of spontaneous breathing during ventilatory support in patients with acute lung injury. *Am J Respir Crit Care Med* 2001; 164:43-49
- [39] Froese AB, Kinsella JP. High-frequency oscillatory ventilation: lessons from the neonatal/pediatric experience. *Crit Care Med*. 2005 Mar;33 (3 Suppl):S115-21.
- [40] Keszler M, Durand DJ. Neonatal high-frequency ventilation. Past, present, and future. *Clin Perinatol*. 2001 Sep;28(3):579-607. [Links](#)
- [41] Lapinsky SE, Mehta S. Bench-to-bedside review: Recruitment and recruiting maneuvers. *Crit Care*. 2005 Feb;9(1):60-5. Epub 2004 Aug 18.
- [42] Rimensberger PC, Pristine G, Mullen BM, Cox PN, Slutsky AS. Lung recruitment during small tidal volume ventilation allows minimal positive end-expiratory pressure without augmenting lung injury. *Crit Care Med*. 1999 Sep;27(9):1940-5.
- [43] Rimensberger PC, Cox PN, Frndova H, Bryan AC. The open lung during small tidal volume ventilation: concepts of recruitment and “optimal” positive end-expiratory pressure. *Crit Care Med*. 1999 Sep;27(9):1946-52.
- [44] Carvalho AR, Jandre FC, Pino AV, Bozza FA, Salluh JI, Rodrigues R, Ascoli FO, Giannella-Neto A. Positive end-expiratory pressure at minimal respiratory elastance represents the best compromise between mechanical stress and lung aeration in oleic acid induced lung injury. *Crit Care* 2007, 11:R86

- [45] Nielsen J, Nilsson M, Fredén F, Hultman J, Alström U, Kjaergaard J, Hedenstierna G, Larsson A. Central hemodynamics during lung recruitment maneuvers at hypovolemia, normovolemia and hypervolemia. A study by echocardiography and continuous pulmonary artery flow measurements in lung-injured pigs. *Intensive Care Med.* 2006 Apr;32(4):585-94. Epub 2006 Mar 7.
- [46] Odenstedt H, Lindgren S, Olegård C, Erlandsson K, Lethvall S, Aneman A, Stenqvist O, Lundin S. Slow moderate pressure recruitment maneuver minimizes negative circulatory and lung mechanic side effects: evaluation of recruitment maneuvers using electric impedance tomography. *Intensive Care Med.* 2005 Dec;31(12):1706-14. Epub 2005 Sep 22.
- [47] Villagrà A, Ochagavía A, Vatua S, Murias G, Del Mar Fernández M, Lopez Aguilar J, Fernández R, Blanch L. Recruitment maneuvers during lung protective ventilation in acute respiratory distress syndrome. *Am J Respir Crit Care Med.* 2002 Jan 15;165(2):165-70.
- [48] Musch G, Harris RS, Vidal Melo MF, O'Neill KR, Layfield JD, Winkler T, Venegas JG. Mechanism by which a sustained inflation can worsen oxygenation in acute lung injury. *Anesthesiology.* 2004 Feb;100(2):323-30. Erratum in: *Anesthesiology.* 2004 May;100(5):1336.
- [49] Grasso S, Mascia L, Del Turco M, Malacarne P, Giunta F, Brochard L, Slutsky AS, Marco Ranieri V. Effects of recruiting maneuvers in patients with acute respiratory distress syndrome ventilated with protective ventilatory strategy. *Anesthesiology.* 2002 Apr;96(4):795-802.
- [50] Riva DR, Oliveira MB, Rzezinski AF, Rangel G, Capelozzi VL, Zin WA, Morales MM, Pelosi P, Rocco PR. Recruitment maneuver in pulmonary and extrapulmonary experimental acute lung injury. *Crit Care Med.* 2008 Jun;36(6):1900-8.

- [51] Gattinoni L, Pelosi P, Suter PM, Pedoto A, Vercesi P, Lissoni A. Acute respiratory distress syndrome caused by pulmonary and extrapulmonary disease: Different syndromes? *Am J Respir Crit Care Med* 1998; 158:3-11
- [52] Thille AW, Richard JC, Maggiore SM, Ranieri VM, Brochard L. Alveolar recruitment in pulmonary and extrapulmonary acute respiratory distress syndrome: comparison using pressure-volume curve or static compliance. *Anesthesiology*. 2007 Feb;106(2):212-7.
- [53] Brower RG, Morris A, MacIntyre N, Matthay MA, Hayden D, Thompson T, Clemmer T, Lanken PN, Schoenfeld D. ARDS Clinical Trials Network, National Heart, Lung, and Blood Institute, National Institutes of Health. Effects of recruitment maneuvers in patients with acute lung injury and acute respiratory distress syndrome ventilated with high positive end-expiratory pressure. *Crit Care Med*. 2003 Nov;31(11):2592-7. Erratum in: *Crit Care Med*. 2004 Mar;32(3):907.
- [54] Borges JB, Okamoto VN, Matos GF, Carames MP, Arantes PR, Barros F, Souza CE, Victorino JA, Kacmarek RM, Barbas CS, Carvalho CR, Amato MB. Reversibility of lung collapse and hypoxemia in early acute respiratory distress syndrome. *Am J Respir Crit Care Med*. 2006 Aug 1;174(3):268-78. Epub 2006 May 11.
- [55] Kacmarek RM. Strategies to optimize alveolar recruitment. *Curr Opin Crit Care*. 2001 Feb;7(1):15-20.
- [56] Bugedo G, Bruhn A, Hernandez G, Rojas G, Varela C, Tapia JC, Castillo L. Lung computed tomography during a lung recruitment maneuver in patients with acute lung injury. *Intensive Care Med*. 2003 Feb;29(2):218-25. Epub 2003 Jan 18.

- [57] Grasso S, Mascia L, Del Turco M, Malacarne P, Giunta F, Brochard L, Slutsky AS, Marco Ranieri V. Effects of recruiting maneuvers in patients with acute respiratory distress syndrome ventilated with protective ventilatory strategy. *Anesthesiology*. 2002 Apr;96(4):795-802.
- [58] Odenstedt H, Lindgren S, Olegård C, Erlandsson K, Lethvall S, Aneman A, Stenqvist O, Lundin S. Slow moderate pressure recruitment maneuver minimizes negative circulatory and lung mechanic side effects: evaluation of recruitment maneuvers using electric impedance tomography. *Intensive Care Med*. 2005 Dec;31(12):1706-14. Epub 2005 Sep 22.
- [59] Patroniti N, Foti G, Cortinovis B, Maggioni E, Bigatello LM, Cereda M, Pesenti A. Sigh improves gas exchange and lung volume in patients with acute respiratory distress syndrome undergoing pressure support ventilation. *Anesthesiology*. 2002 Apr;96(4):788-94.
- [60] Amato, Barbas, Medeiros et al. Effect of a protective ventilation strategy on mortality in the ARDS. *New England Journal of Medicine* 1998; 338:347-354.
- [61] Brower RG, Matthay MA, Morris A, Schoenfeld D, Taylor Thompson B, Wheeler A. The Acute Respiratory Distress Syndrome Network: Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome. *N Engl J Med*. 2000 May 4;342(18):1301-8.
- [62] Stewart TE, Meade MO, Cook DJ, Granton JT, Hodder RV, Lapinsky SE, Mazer CD, McLean RF, Rogovein TS, Schouten BD, Todd TR, Slutsky AS. Evaluation of a ventilation strategy to prevent barotrauma in patients at high risk for acute respiratory distress syndrome. Pressure- and Volume-Limited Ventilation Strategy Group. *N Engl J Med*. 1998 Feb 5;338(6):355-61.

- [63] Brochard L, Roudot-Thoraval F, Roupie E, Delclaux C, Chastre J, Fernandez-Mondéjar E, Clémenti E, Mancebo J, Factor P, Matamis D, Ranieri M, Blanch L, Rodi G, Mentec H, Dreyfuss D, Ferrer M, Brun-Buisson C, Tobin M, Lemaire F. Tidal volume reduction for prevention of ventilator-induced lung injury in acute respiratory distress syndrome. The Multicenter Trail Group on Tidal Volume reduction in ARDS. *Am J Respir Crit Care Med*. 1998 Dec;158(6):1831-8.
- [64] Brower RG, Shanholtz CB, Fessler HE, Shade DM, White P Jr, Wiener CM, Teeter JG, Dodd-o JM, Almog Y, Piantadosi S. Prospective, randomized, controlled clinical trial comparing traditional versus reduced tidal volume ventilation in acute respiratory distress syndrome patients. *Crit Care Med*. 1999 Aug;27(8):1492-8.
- [65] Brower RG, Lanken PN, MacIntyre N, et al. Higher versus lower positive end-expiratory pressures in patients with the acute respiratory distress syndrome. *N Engl J Med*. 2004;351(4):327-336.
- [66] Meade MO, Cook DJ, Guyatt GH, Slutsky AS, Arabi YM, Cooper DJ, Davies AR, Hand LE, Zhou Q, Thabane L, Austin P, Lapinsky S, Baxter A, Russell J, Skrobik Y, Ronco JJ, Stewart TE. Lung Open Ventilation Study Investigators. Ventilation strategy using low tidal volumes, recruitment maneuvers, and high positive end-expiratory pressure for acute lung injury and acute respiratory distress syndrome: a randomized controlled trial. *JAMA*. 2008 Feb 13;299(6):637-45.
- [67] Mercat A, Richard JC, Vielle B, Jaber S, Osman D, Diehl JL, Lefrant JY, Prat G, Richecoeur J, Nieszkowska A, Gervais C, Baudot J, Bouadma L, Brochard L. Expiratory Pressure (Express) Study Group. Positive end-expiratory pressure setting in adults with acute lung injury and acute respiratory distress syndrome: a randomized controlled trial. *JAMA*. 2008 Feb 13;299(6):646-55.

- [68] Vieillard-Baron A, Charron C, Jardin F. Lung “recruitment” or lung overinflation maneuvers? *Intensive Care Med.* 2006 Jan;32(1):177-8. Epub 2005 Nov 11.
- [69] Lim CM, Soon Lee S, Seoung Lee J, Koh Y, Sun Shim T, Do Lee S, Sung Kim W, Kim DS, Dong Kim W. Morphometric effects of the recruitment maneuver on saline-lavaged canine lungs. A computed tomographic analysis. *Anesthesiology.* 2003 Jul;99(1):71-80.
- [70] Caramenz MP, Kacmarek RM, Helmy M, Miyoshi E, Malhotra A, Amato MB, Harris RS. A comparison of methods to identify open-lung PEEP. *Intensive Care Med.* 2009 Apr;35(4):740-7. Epub 2009 Jan 31.
- [71] Carvalho AR, Jandre FC, Pino AV, Bozza FA, Salluh JI, Rodrigues R, Ascoli FO, Giannella-Neto A. Positive end-expiratory pressure at minimal respiratory elastance represents the best compromise between mechanical stress and lung aeration in oleic acid induced lung injury. *Crit Care* 2007, 11:R86.
- [72] Schreiter D, Reske A, Stichert B, Seiwerts M, Bohm SH, Kloepfel R, Josten C. Alveolar recruitment in combination with sufficient positive end-expiratory pressure increases oxygenation and lung aeration in patients with severe chest trauma. *Crit Care Med.* 2004 Apr;32(4):968-75.
- [73] Badet M, Bayle F, Richard JC, Guérin C. Comparison of optimal positive end-expiratory pressure and recruitment maneuvers during lung-protective mechanical ventilation in patients with acute lung injury/acute respiratory distress syndrome. *Respir Care.* 2009 Jul;54(7):847-54.

- [74] Hickling KG. Best compliance during a decremental, but not incremental, positive end-expiratory pressure trial is related to open-lung positive end-expiratory pressure: a mathematical model of acute respiratory distress syndrome lungs. *Am J Respir Crit Care Med*. 2001 Jan;163(1):69-78
- [75] Toth I, Leiner T, Mikor A, Szakmany T, Bogar L, Molnar Z. Hemodynamic and respiratory changes during lung recruitment and descending optimal positive end-expiratory pressure titration in patients with acute respiratory distress syndrome. *Crit Care Med*. 2007 Mar;35(3):787-93.
- [76] Albaiceta GM, Piacentini E, Villagr a A, Lopez-Aguilar J, Taboada F, Blanch L. Application of continuous positive airway pressure to trace static pressure-volume curves of the respiratory system. *Crit Care Med*. 2003 Oct;31(10):2514-9.
- [77] Suter PM, Fairley B, Isenberg MD. Optimum end-expiratory airway pressure in patients with acute pulmonary failure. *N Engl J Med* 1975;292(6):284-289.
- [78] Brochard L. What is a pressure-volume curve? *Crit Care*. 2006;10(4):156. Review.
- [79] Albaiceta GM, Luyando LH, Parra D, Menendez R, Calvo J, Pedreira PR, Taboada F. Inspiratory vs. expiratory pressure-volume curves to set end-expiratory pressure in acute lung injury. *Intensive Care Med*. 2005 Oct;31(10):1370-8. Epub 2005 Aug 10.
- [80] Mankikian B, Lemaire F, Benito S, Brun-Buisson C, Harf A, Maillot JP, Becker J. A new device for measurement of pulmonary pressure-volume curves in patients on mechanical ventilation. *Crit Care Med*. 1983 Nov; 11(11): 897-901.

- [81] Blanc Q, Sab JM, Philit F, Langevin B, Thouret JM, Noel P, Robert D, Guerin C. Inspiratory pressure-volume curves obtained using automated low constant flow inflation and automated occlusion methods in ARDS patients with a new device. *Intensive Care Med.* 2002 Jul;28(7):990-4. Epub 2002 Jun 12.
- [82] Gama AM, Meyer EC, Gaudencio AM, Grunauer MA, Amato MB, de Carvalho CR, Barbas CS. Different low constant flows can equally determine the lower inflection point in acute respiratory distress syndrome patients. *Artif Organs.* 2001 Nov; 25(11): 882-9.
- [83] Servillo G, Svantesson C, Beydon L, Roupie E, Brochard L, Lemaire F, Jonson B. Pressure-volume curves in acute respiratory failure: automated low flow inflation versus occlusion. *Am J Respir Crit Care Med.* 1997 May;155(5):1629-36.
- [84] Lu Q, Vieira SRR, Richecoeur J, Puybasset L, Kalfon P, Coriat P, Rouby JJ. A simple automated method for measuring pressure-volume curves during mechanical ventilation. *Am J Respir Crit Care Med* 1999;159:275-282.
- [85] Piacentini E, Wysocki M, Blanch L. A new automated method versus continuous positive airway pressure method for measuring pressure-volume curves in patients with acute lung injury. *Intensive Care Med.* 2009 Mar;35(3):565-70. Epub 2008 Oct 14.
- [86] Artigas A, Bernard GB, Carlet J, et al. The American-European Consensus Conference on ARDS, part 2. Ventilatory, pharmacologic, supportive therapy, study design strategies, and issues related to recovery and remodeling. *Am J Respir Crit Care Med* 1998; 157: 1332± 1347.

- [87] Matamis D, Lemaire F, Harf A, et al. Total respiratory pressure-volume curves in the adult respiratory distress syndrome. *Chest* 1984; 86: 58± 66.
- [88] Jonson B, Richard JC, Straus C, Mancebo J, Lemaire F, Brochard L. Pressure-volume curves and compliance in acute lung injury: evidence of recruitment above the lower inflection point. *Am J Respir Crit Care Med*. 1999 Apr; 159(4 Pt 1): 1172-8.
- [89] Hickling KG. The pressure-volume curve is greatly modified by recruitment. A mathematical model of ARDS lungs. *Am J Respir Crit Care Med*. 1998 Jul;158(1):194-202.
- [90] Takeuchi M, Goddon S, Dolhnikoff M, Shimaoka M, Hess D, Amato MB, Kacmarek RM. Set positive end-expiratory pressure during protective ventilation affects lung injury. *Anesthesiology*. 2002 Sep;97(3):682-92.
- [91] Kallet RH. Pressure-volume curves in the management of acute respiratory distress syndrome. *Respir Care Clin N Am*. 2003 Sep; 9(3): 321-41.
- [92] Creamer KM, McCloud LL, Fisher LE, Ehrhart IC. Closing pressure rather than opening pressure determines optimal positive end-expiratory pressure and avoids overdistention. *Chest*. 1999 Jul;116 (1 Suppl):26S-27S.
- [93] Rimensberger PC, Cox PN, Frndova H, Bryan AC. The open lung during small tidal volume ventilation: concepts of recruitment and "optimal" positive end-expiratory pressure. *Crit Care Med*. 1999 Sep;27(9):1946-52.

- [94] Hickling KG. Best compliance during a decremental, but not incremental, positive end-expiratory pressure trial is related to opening positive end-expiratory pressure: a mathematical model of acute respiratory distress syndrome lungs. *Am J Respir Crit Care Med.* 2001 Jan;163(1):69-78
- [95] Nishida T, Suchodolski K, Schettino GP, Sedek K, Takeuchi M, Kacmarek RM. Peak volume history and peak pressure-volume curve pressures independently affect the shape of the pressure-volume curve of the respiratory system. *Crit Care Med.* 2004 Jun;32(6):1358-64.
- [96] Takeuchi M, Sedek KA, Schettino GP, Suchodolski K, Kacmarek RM. Peak pressure during volume history and pressure-volume curve measurement affects analysis. *Am J Respir Crit Care Med.* 2001 Oct 1;164(7):1225-30.
- [97] Hickling KG. The pressure-volume curve is greatly modified by recruitment. A mathematical model of ARDS lungs *Am J Respir Crit Care Med.* 1998 Jul;158(1):194-202
- [98] Mehta S, Stewart TE, MacDonald R, Hallett D, Banayan D, Lapinsky S, Slutsky A. Temporal change, reproducibility, and interobserver variability in pressure-volume curves in adults with acute lung injury and acute respiratory distress syndrome. *Crit Care Med.* 2003 Aug; 31(8): 2118-25.
- [99] Harris RS, Hess DR, Venegas JG. An objective analysis of the pressure-volume curve in the acute respiratory distress syndrome. *Am J Respir Crit Care Med* 2000, 161:432-9.
- [100] Venegas JG, Harris RS, Simon BA. A comprehensive equation for the pulmonary pressure-volume curve. *J Appl Physiol.* 1998 Jan; 84(1): 389-95.

- [101] Heller H, Brandt S, Schuster KD. Development of an algorithm for improving the description of the pulmonary pressure-volume curve. *J Appl Physiol* 2002, 92:1770.
- [102] Henzler D, Orfao S, Rossaint R, Kuhlen R. Modification of a sigmoidal equation for the pulmonary pressure-volume curve for asymmetric data. *J Appl Physiol* 2003, 95:2183-4. author reply 2184
- [103] Albaiceta GM, Garcia E, Taboada F. Comparative study of four sigmoid models of pressure-volume curve in acute lung injury. *Biomed Eng Online*. 2007 Feb 14;6:7.
- [104] Nielsen J, Ostergaard M, Kjaergaard J, et al. Lung recruitment maneuver depresses central haemodynamics in patients following cardiac surgery. *Intensive Care Med* 2005; 31:1189-1194
- [105] Nielsen J, Nilsson M, Fredén F, Hultman J, Alström U, Kjaergaard J, Hedenstierna G, Larsson A. Central hemodynamics during lung recruitment maneuvers at hypovolemia, normovolemia and hypervolemia. A study by echocardiography and continuous pulmonary artery flow measurements in lung-injured pigs. *Intensive Care Med*. 2006 Apr;32(4):585-94. Epub 2006 Mar 7.
- [106] Dall'ava-Santucci J, Armaganidis A, Brunet F, Dhainaut JF, Chelucci GL, Monsallier JF, Lockhart A. Causes of error of respiratory pressure-volume curves in paralyzed subjects. *J Appl Physiol*. 1998 Jan;64(1):42-9
- [107] Henzler D, Mahnken A, Dembinski R, Waskowiak B, Rossaint R, Kuhlen R. Repeated generation of the pulmonary pressure-volume curve may lead to derecruitment in experimental lung injury. *Intensive Care Med*. 2005 Feb;31(2):302-10. Epub 2004 Dec 9.

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Dräger Medical Hispania S.A.
Dräger Medical México, S.A. de C.V.
Dräger Medical GmbH se convertirá
en Drägerwerk AG & Co. KGaA.