Supplement

Oxylog 3000 *plus*

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WARNING
To properly use this medical device,
read and comply with the instructions
for use and this supplement.

Emergency and Transport Ventilator

Supplement to the instructions for use

WARNING

To properly use this medical device, read and comply with the instructions for use and this supplement.

Instructions for use	Part number	Edition
Oxylog 3000 <i>plus</i> en	5705305	12 – 2020-01
Oxylog 3000 <i>plus</i> zh	5705324	17 – 2020-11
Oxylog 3000 <i>plus</i> es	5705309	11 – 2020-01
Oxylog 3000 <i>plus</i> ptBR	9052937	09 – 2020-01

- Keep this supplement with the instructions for use of the medical device.

This supplement completes the information in the instructions for use in the following chapters:

For Your Safety and that of Your Patients

Installing accessories

WARNING

Risk of incorrect operation and of incorrect use
Follow the instructions for use of all accessories carefully.

Abbreviations

Abbreviation	Explanation
IP34	Level 3 protection against foreign bodies and contact, Level 4 protection against ingress of liquids

Assembly

Connecting the disposable breathing circuit for pediatric patients

WARNING

Risk of patient injury
If the tidal volume is less than 250 mL, use a pediatric breathing circuit.
If the tidal volume is more than 250 mL, use an adult breathing circuit.

Preparation

Zero calibration before ventilation

WARNING

Risk of patient injury
The CO₂ measurement can be inaccurate and lead to a risk of patient injury.

- Before the CO₂ filter check can be performed, perform a CO₂ zero calibration.
- Do not breathe on the CO₂ sensor during the CO₂ zero calibration and the CO₂ filter check.
- Only perform the CO₂ zero calibration and the CO₂ filter check at ambient temperatures between 10 and 40 °C (50 and 104 °F).
- Do not perform the CO₂ zero calibration when exposed to direct sunlight.

Monitoring

Performing zero calibration during ventilation

WARNING

Risk of patient injury
The CO₂ measurement can be inaccurate and lead to a risk of patient injury.

- Before the CO₂ filter check can be performed, perform a CO₂ zero calibration.
- Do not breathe on the CO₂ sensor during the CO₂ zero calibration and the CO₂ filter check.
- Only perform the CO₂ zero calibration and the CO₂ filter check at ambient temperatures between 10 and 40 °C (50 and 104 °F).
- Do not perform the CO₂ zero calibration when exposed to direct sunlight.

Configuration

Performing zero calibration

WARNING

Risk of patient injury
The CO₂ measurement can be inaccurate and lead to a risk of patient injury.

- Before the CO₂ filter check can be performed, perform a CO₂ zero calibration.
- Do not breathe on the CO₂ sensor during the CO₂ zero calibration and the CO₂ filter check.
- Only perform the CO₂ zero calibration and the CO₂ filter check at ambient temperatures between 10 and 40 °C (50 and 104 °F).
- Do not perform the CO₂ zero calibration when exposed to direct sunlight.

Problem solving

Messages in the alarm message field

Alarm	Cause	Remedy	Rank
<i>! CO₂ inaccurate or failed</i>	The CO ₂ measurement is inaccurate.	Use external CO ₂ monitoring. Contact DrägerService.	42
	The CO ₂ sensor is defective.	Replace the faulty CO ₂ sensor.	
	The zero calibration or filter check of the CO ₂ sensor was activated and no valid measured values were generated.	Wait until the CO ₂ zero calibration/filter check is completed.	
<i>!! etCO₂ low</i>	The lower alarm limit for the end-tidal CO ₂ concentration has been exceeded.	Check patient condition.	27
		Check alarm limits.	
		Adjust the alarm limit, if necessary.	
	The CO ₂ sensor has fallen off the cuvette.	Reattach the CO ₂ sensor to the cuvette.	
<i>! No int. battery charging</i>	Internal battery cannot be charged due to a faulty battery or too hot or cold environment.	Press the Alarm Reset key to confirm the alarm. Remove the internal battery and reinsert or replace it.	44

Technical data

Device specifications

Protection class against ingress of liquids	IP34, the device is splash-proof from all directions and protected against foreign objects >2.5 mm (0.1 in)
---	---

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© Drägerwerk AG & Co. KGaA
Drägerwerk AG & Co. KGaA