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In the chapter "For Your Safety and that of Your Patients" the following additional information must be observed:

Emergency method

The operator of the health facility must develop an emergency procedure to respond to catastrophic failure of the connected pipeline system when the medical gas supply to all medical devices is interrupted simultaneously.

Room requirements

Supply systems with gas cylinders and cylinders containing non-cryogenic liquids must not be located in the same room as medical air compressors, oxygen concentrators or vacuum supply systems.

In the chapter "Operation" the following additional information must be observed:

Pressure display

WARNING

The connectors of the individual gas cylinders of the manifold are equipped with a non-return valve preventing complete emptying of the manifold in case of a leakage at one connector.

The use of the non-return valves has the effect that the pressure gauge of the RS 20/RS 80 only shows the pressure of the gas cylinder with the highest pressure in the manifold. Single empty gas cylinders could remain undiscovered.

Particularly the individual gas cylinders of standby manifolds must therefore be checked for correct filling at regular intervals.

To make sure that the complete reserve volume can be provided at any time, the standby manifold can be used regularly approx. every 3 months for the normal supply and afterwards be replaced by full gas cylinders.

Checking the standby manifold

CAUTION

After attaching the gas cylinder and on a regular basis in the normal operation, the gas cylinder manifold must be checked for leakages to prevent unnoticed emptying particularly of the standby manifold.

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© Drägerwerk AG & Co. KGaA
Drägerwerk AG & Co. KGaA