
RS 80/2 Reducing Station Redukční stanice

**For Medical Gas Pipeline Systems
Instructions for Use**

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For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

Maintenance

The apparatus must be inspected and serviced regularly by trained service personnel at six monthly intervals. Repair and general overhaul of the apparatus may only be carried out by trained service personnel. We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Only authentic Dräger spare parts may be used for maintenance. Observe chapter "Maintenance Intervals".

Accessories

Do not use accessory parts other than those in the order list.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by DrägerService or if the apparatus is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

Drägerwerk AG & Co. KGaA

Regulations

The national laws and regulations governing pressure reducers and compressed-gas systems are to be observed.

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Intended Use

RS 80/2 – Reducing Station for Medical Gas Pipeline Systems, for reducing and monitoring gas from a loop system.


Suitable for:

Oxygen, nitrous oxide, carbon dioxide, nitrogen as well as for inert, non-corrosive gases.

May not be used for:


Corrosive, aggressive or toxic gases, acetylene, propane, butane and others flammable gases.

RS 80/2
for O₂ and
not flammable gases (NB)

RS 80/2
pro O₂ a
nehořlavé plyny (NB)

RS 80/2
for AIR 5 bar*
and AIR 8 bar

RS 80/2
pro AIR 5 bar*
a AIR 8 bar

* 1 bar = 1 kPa x 100

* 1 bar = 1 kPa x 100

Testing and Commissioning

Pressure reducing stations may not be operated until they have been tested and commissioned by trained and qualified personnel.

National regulations must be observed.

In Europe, e. g. EN 737-3.

When installation or maintenance procedures have been carried out, a comprehensive testing and commissioning programme must be undertaken by trained personnel.

These tests determine:

1. whether the safety requirements for the protection of patients and staff have been fulfilled –

and

2. whether the performance characteristics of the medical gas pipeline system are met.

Written records of the tests must be kept.

Following this testing programme, the pipeline system which is ready for operation and the documentation are handed over to the user, and users are then given instructions.

The handing-over procedure must be formally recorded.

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Installation

- Only by trained service personnel and according to Installation Instructions "Reducing Station RS 80/2"

Installation of additional accessories

- Observe the installation instructions for the accessories.
- Pressure gauges and pressure switch gauges can be replaced under pressure. A restriction valve in the connection prevents large-scale leaking.

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
Operation

Do not allow equipment for oxygen or other oxidizing gases to come into contact with oil and grease, or other lubricants which are not certified for use with oxygen by the appropriate authority.

In the Federal Republic of Germany:
BAM, Berlin (Federal Institute for Material Research and Testing).

- 1,2** Pressure reducers in use.
- 3** Valves open.
- 4** Valves open depending on which outlets were selected.
- 5** NIST connection closed.

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
Emergency Supply Unit

- Inform the operator.
- 1 Connect emergency supply 5 bar or 8 bar of the relevant gas type to the NIST connection.
- 2 Close valves.

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
Emergency Access

In an emergency the device can also be opened without a key.

- Hit button with fist – remove cover.

Replace button into cover:

- Open lock using key.
- Insert button into cover and lock device using key.

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
Shut-down

The cylinder manifold may only be shut down if no patients are going to be put at risk by turning off the compressed gas supply.

Follow internal regulations and any legal requirements in Europe, such as EN 737-3.

Attach information labels.

- 1 Close valves.
- 2 If connected:
Remove the emergency supply from the NIST connector.
- 3 By manually operating the safety valve, allow the gas to escape into the open air.

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Fault-Cause-Remedy

Fault	Cause	Remedy
Breakdown in operating pressur	Breakdown of gas supply	Establish emergency supply. Reduce gas consumption. Have repaired by experts.
Operating pressure too low (P2 < 4 bar) bzw. (P2 < 6.4 bar)	Inadequate gas supply	Establish emergency supply. Reduce gas consumption. Pressure OK after these measures? If not, have repaired by experts.
Operating pressure too high (P2 > 6 bar) bzw. (P2 > 9.6 bar)	Increased gas consumption	Have repaired by experts.
Increased gas consumption	Safety valve blows off	Observe pressure gauge. Have repaired by experts.
Increased gas consumption	Minor leakage	Have repaired by experts.
Increased gas consumption Inadequate gas supply Gas accumulates in the room	Major Leckage	No fire, do not smoke! Establish emergency supply. Reduce gas consumption. If possible, shut off respective range of gas supply. Have repaired by experts.
Pressure gauge does not show required pressure (5 bar / 8 bar)	Pressure setting too low or too high	Re-adjust pressure using regulating screw of pressure reducer, then relock regulating screw.
	Pressure gauge defective	Have repaired by experts.
Delivery pressure increases at zero delivery	Safety valve blows off	Have repaired by experts.

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Maintenance Intervals

Regularly, at least every six months

- Inspection of the entire system by trained service personnel.
We recommend that a service contract is obtained with DrägerService.

Regularly, at least once a year

- Check pressure reducers for set values and subsequent pressure rise.
- Check internal alarm thresholds / alarm messages.
- Check function and condition of safety valve.
- Make leak test of the entire system.

Every six years

- Complete overhaul / replacement of the following:
 - pressure reducers (as obligatory for safety reasons).
 - safety valves and blow-off valves.
 - sintered metal filter inserts.
 - all Elastomer parts, such as O-rings, pressure reducer diaphragm, closing bolts.
- Use only original Dräger replacement parts for all maintenance procedures.

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Service



Replace safety valve

The safety valve can be replaced without interrupting operation.

- 1 Under no circumstances should the adaptor be loosened, as this would lead to gas escaping. For safety reasons the adaptor is secured with Loctite.
- 2 Unscrew the cap nut (a small amount of gas may escape)

While safety valve is being replaced the mains fusing is interrupted.


- 3 Observe operating pressure gauge. The pressure must not exceed 6 or 9 bar.
- Reinsert the new or repaired safety valve as quickly as possible and check functioning by turning the knurled screw.

Replace pressure reducer

Replace the pressure reducers individually. The second pressure reducer must remain operational.

- 4 Replace pressure reducer:
 - 5 Close valves.
- 6 Replace pressure reducer:
 - 7 Close valves.
- Attach to base plate and unscrew cap nuts on left and right-hand sides of pressure reducer. Pull out pressure reducer and replace.

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
Co je co

RS 80/2 O₂ a N₂O


- 1 Pipeline system outlet
- 2 Pressure gauge O₂ and/or non-flammable gases
- 3 Connections for valve block (for retrofitting)
- 4 Safety valve
- 5 Vent pipe
- 6 Non-return valve
- 7 Emergency supply unit NIST-connection O₂ or N₂O
- 8 Pressure reducer, right
- 9 Screw connection (input)
- 10 Pressure reducer, left
- 11 Ball valve
- 12 Pressure switch gauge (exceeded and falling below required pressure)
- 13 Valve set
- 14 Terminals for pressure switches

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RS 80/2 AIR 5 and AIR 8

RS 80/2 AIR 5 a AIR 8

- 1 Pipeline system outlet
- 2 Pressure gauge non-flammable gases
- 3 Connections for valve block (for retrofitting)
- 4 Safety valve
- 5 Non-return valve
- 6 Emergency supply unit NIST-connection AIR
- 7 Pressure reducer, right
- 8 Screw connection (input)
- 9 Pressure reducer, left
- 10 Ball valve
- 11 Pressure switch gauge (exceeded and falling below required pressure)
- 12 Valve set
- 13 Terminals for pressure switches

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
Popis funkcí

Flow chart

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A,B Pressure reducers with test connector

- C Safety valve
- D Pressure switch gauge (for closing and opening)
- 1 Screw connection (input)
- 2 Screw connection for emergency supply
- 3–8 Ball valve
- 9–11 Connections for valve blocks
- 12,13 Operating pressure gauges
- 14 Vent pipe for safety valve (not applicable if using compressed air)
- 15 Pressure gauge for supply pressure
- P_H Delivery pressure
- P_V Supply pressure

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Pressure reducers A,B

These reduce the current supply pressure to 5 bar or 8 bar delivery pressure.

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Safety valve C

The function of the safety valve C is to protect the section of the unit following the valve from excess pressure. The gas released is lead off via pipe **14**. The pipe line must lead to the open air from the reducing station for gases other than compressed air, in order to ensure that gas can be blown off safely.

Pressure switch gauge D

In order to ensure safe operation, the pressure switch gauge monitors too high or too low pressures behind the pressure reducers. Alarms are passed on as a mal-function message to an emergency alarm panel or to the building management system (BMS).

The pressure switch gauge **D** opens in the event of pressure is falling below or is exceeding the required level (closed-circuit current principle).

Valves 3, 4, 5 and 6

are used to lock the pressure reducers **A** and **B** in the event of repairs. If **A** is being repaired or removed, then **3** and **5** must be shut, and for work on **B** valves **4** and **6** must be shut.

Connection 9-11

are for retrofitting of valve blocks for further gas distribution.

Pressure gauges 12, 13

display the current outlet pressure (delivery pressure) after the pressure reducers (normally 5 bar or 8 bar).

Pressure gauge 15

displays the current pressure (supply pressure) before each pressure reducer.

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Directive 93/42/EEC
concerning Medical Devices

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© Drägerwerk AG & Co. KGaA