Drägerwerk AG & Co. KGaA, 23558 Lübeck

To whom it may concern

Our reference

COVID-19 Supplement Handling of Dräger
Anesthesia Workstation

Phone

+49 451 882-0

April 2, 2020

**Supplement to “SARS-CoV-2 and handling of Dräger Anesthesia Workstation”, updated
March 24th, 2020**

(latest version can be obtained on:

https://www.draeger.com/en_corp/Corporate/Coronavirus-COVID-19)

**System setup recommendations for the use of Dräger anesthesia devices in the scope of
their intended use for neonatal and pediatric patients**

Please note: Additional and separate information is available for a potential off-label use for long-term ventilation. Refer to the URL stated above.

Background:

Specific patient groups cannot be sufficiently ventilated with mechanical filter (BSF or HMEF) at the Y-piece of the breathing circuit. Therefore, this results in the need to keep the dead space between the Y-piece and the airway connector to the patient at a low level given the small tidal volumes of those patients.

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System setup recommendation for neonatal patients (tidal volume approx. < 30 ml):

1. For neonatal patients, it is frequently necessary to sample breathing gas for the monitoring of CO₂, O₂, N₂O and anesthetic agents as close to the airways as possible. Otherwise the quality of the measured gas parameters might not be sufficient.
2. For neonatal patients, use an HME **at the Y-piece** for sufficient humidification of the airways.
e.g. Dräger HumidStar 2 (MP01745) for tidal volumes above 10ml
3. Use mechanical BSFs (adult) **at each the expiratory and the inspiratory port of the anesthesia device**.
e.g. Dräger SafeStar 55 (MP01790)


Figure 1: System setup recommendation for neonatal patients (tidal volume approx. < 30 ml)

- Check that all connections are securely fitted and tight.
- Perform a system test of the anesthesia device after the breathing circuit (breathing hoses, BSFs, HME, etc.) has been completely installed and prior to connecting a patient.
- Regularly check the inside of the BSFs and HME for liquids or visible soiling (secretion). If liquids or visible soiling are found, replace the BSFs and/or HME immediately.
Set the alarm limits for MV low and Paw high to suitable values.

Caution: If this is not done, there is a risk of pressure build-up and insufficient patient ventilation and of drawing air from the patient's lung

- After the completion of a case and before reprocessing the anesthesia workplace remove all components of the breathing circuit:
gas sample line, water trap, connector, Y-piece, HME, breathing hoses and BSFs.
- Dispose of all components safely.

System setup recommendation for pediatric patients (tidal volume approx. 30 - 300 ml):

1. For pediatric patients it is usually possible to sample breathing gas for monitoring of CO₂, O₂, N₂O and anesthetic agents at the HME/ HMEF (electrostatic) at the Y-piece.
For types and part numbers please see Dräger Accessory Catalogue.

Caution: The gas sampling port must be on the device side of the HME / HMEF. Otherwise breathing gas could get drawn from the patient's lung if the HME/ HMEF gets clogged.

2. Use mechanical BSFs (adult) e.g. Dräger SafeStar 55 (MP01790) **at each the expiratory and the inspiratory port of the anesthesia device.**


Figure 2: System setup recommendation for pediatric patients (tidal volume approx. 30 - 300 ml)

- Check that all connections are securely fitted and tight.
- Perform a system test of the anesthesia device after the breathing circuit (breathing hose, BSFs, HME / HMEF, etc.) has been completely installed and prior to connecting a patient.
- Regularly check the inside of the BSFs and HME / HMEF for liquids or visible soiling (secretion). The lower the freshgas flow and the greater the CO₂ production of the patient, the more humidity may accumulate in the breathing circuit.*
- If liquids or visible soiling are found, replace the BSF and/or HME / HMEF immediately.
Set the alarm limits for MV low and Paw high to suitable values.

* The accumulation of humidity within the breathing circuit can be reduced by increasing the FG flow.

Furthermore, the user can replace the soda lime before it is completely used up e.g. after the colour change Soda Lime 800+ indicates that 2/3 of the soda lime has been used up.

Caution: If this is not done, there is a risk of pressure build-up and insufficient patient ventilation.

- After completion of a case and before reprocessing the anesthesia workplace exchange all components of the breathing circuit:
Gas sample line, water trap, connector, Y-piece, HME/ HMEF, breathing hoses and BSFs.
- dispose all of these components safely

We kindly ask you to take this information into consideration and follow the respective Instructions for Use.

If you have any further questions or need any support, please contact Dräger's local service support or sales representative.

Kind Regards,

Meike Lessau
Global Product Manager
Business Unit Hospital
Consumables & Accessories
Medical Division

Matthias Hurlebaus
Product Manager
Business Area Anesthesiology
Business Unit Therapy
Medical Division

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