


D-11424-2016

## CO<sub>2</sub> measurement

During expiration, CO<sub>2</sub> is one of the most important parameter for analysing ventilation efficiency<sup>1</sup>. This is why precise and reliable monitoring of CO<sub>2</sub> concentrations of patients is so crucial.

“Within the expiration phase, CO<sub>2</sub> is one of the most important parameters, providing valuable information on the efficacy of ventilation, gas exchange and metabolism.”<sup>2</sup>

### Stabilize the patient and protect the lung as second step of the Respiration Pathway



#### MECHANICAL VENTILATION

As non-invasive as possible, as invasive as necessary. Along the Respiration Pathway a variance and diversity of treatment tools clearly improve the clinical decision-making.

**Continuous monitoring of CO<sub>2</sub>** concentrations can serve as an early warning system with regard to changes in the acuity level of critically ill patients.<sup>2</sup> The mainstream infrared absorption measurement (real time method) of CO<sub>2</sub> directly at the Y-piece provides reliable data which are displayed in real time on the ventilator screen.

## The CO<sub>2</sub> production VCO<sub>2</sub> is a measure of the physical stress on a patient.

#### WHY CO<sub>2</sub> MEASUREMENT?

The CO<sub>2</sub> option offers more than just CO<sub>2</sub> measurement.

- Display capnometry and capnography simultaneously, together with other patient parameters such as flow, volume/pressure curves and data<sup>6,7</sup>
- Enhance patient care with individually adjustable CO<sub>2</sub> alarm settings and optimized ventilation management
- Obtain accurate measurements even under challenging conditions, such as during active humidification
- Avoid gas leakage and resulting misrepresentation of ventilation settings
- Rapidly verify patient CO<sub>2</sub> values with direct sensor control via the reference filter
- Perform zero calibration without disconnecting the patient from the ventilator
- Combine the CO<sub>2</sub> option with optional mask ventilation (NIV/NIVplus) for real time CO<sub>2</sub> monitoring even during NIV
- Take full advantage of the automated weaning protocol SmartCare®/PS

#### THE PHYSIOLOGICAL CAPNOGRAM

A – B: Emptying of the upper dead space of the airways  
The CO<sub>2</sub> concentration in this section of the curve equals zero, as this is the first phase of expiration during which air from the upper airways, which has not been involved in the process of gas exchange, is analysed.

B – C: Gas from the lower dead space and alveoli  
The CO<sub>2</sub> concentration increases continuously, as the air being analysed comes partly from the upper airways and partly from the alveoli which are rich in CO<sub>2</sub>.

C – D: Alveolar gas


This phase is described as the “alveolar plateau”.  
The curve rises very slowly. The air being analysed comes mainly from the alveolar area.

D: End-tidal CO<sub>2</sub> partial pressure

Represents the highest possible concentration of exhaled CO<sub>2</sub> and is reached at the end of expiration. This point is described as end-tidal CO<sub>2</sub> (etCO<sub>2</sub>) and represents the final portion of air which was involved in the exchange of gases in the alveolar area. It thus represents under certain conditions a reliable index of CO<sub>2</sub> partial pressure in the arterial blood.

D – E: Inspiration

The CO<sub>2</sub> concentration falls rapidly, as fresh gas not containing CO<sub>2</sub> forces its way into the airways at the beginning of inspiration.




<sup>1</sup> St. John RE: End-tidal carbon dioxide monitoring Crit Care Nurs Vol 23, No. 4, August 2003; 83-88

<sup>2</sup> St. John RE: Exhaled gas analysis: technical and clinical aspects of capnography and oxygen consumption. Crit Care Nurs Clin N Am. 1989; 20:363-374

# Monitoring the dead space V<sub>ds</sub> reflects the current patient situation and indicates respiratory insufficiency.


- In intensive care, capnography is applied as a non-invasive way for evaluating a patient's ventilator status<sup>3</sup>.
- It can be used to assess changes in ventilation, pulmonary perfusion and metabolism to support optimization of ventilation settings<sup>2</sup>.
- Observing the arterial and etCO<sub>2</sub> difference respective gradient over a period of time can provide important information, related to either improved or worsening patient clinical status<sup>1</sup> and thus support increased patient safety.
- Capnography monitoring has been increasingly used in operating rooms, intensive care units and emergency departments to indicate incorrect intubation and to monitor cardiopulmonary resuscitation effectiveness<sup>4,5</sup>.
- CO<sub>2</sub> measurement can provide also ongoing measurement of dead space values and CO<sub>2</sub> production.



D-15892-2019



CO<sub>2</sub> curve and value (Evita V800)



Display of alarm limits (Evita V800): etCO<sub>2</sub> high, etCO<sub>2</sub> low



Positioning of CO<sub>2</sub> cuvette and CO<sub>2</sub> sensor

D-15892-2018

## YOUR BENEFITS WITH CO<sub>2</sub> MAINSTREAM MEASUREMENT:

- quick and easy non-invasive CO<sub>2</sub> measurement
- display of CO<sub>2</sub> values as end tidal CO<sub>2</sub> concentration and continuous CO<sub>2</sub> curves
- values such as VCO<sub>2</sub>, V<sub>ds</sub>, V<sub>ds</sub>/V<sub>te</sub>, VT<sub>CO<sub>2</sub></sub> and Slope phase III are also available\*

Evita V800/V600	Babylog VN800/VN600	Evita V500/V300	Babylog VN500	Savina 300 Select/Classic	Oxylog 3000 plus	Oxylog VE300
etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric	etCO <sub>2</sub> numeric
etCO <sub>2</sub> curve	etCO <sub>2</sub> curve	etCO <sub>2</sub> curve	etCO <sub>2</sub> curve	etCO <sub>2</sub> curve	etCO <sub>2</sub> curve	etCO <sub>2</sub> curve
VCO <sub>2</sub>	VCO <sub>2</sub>	VCO <sub>2</sub>	VCO <sub>2</sub>			
VT <sub>CO<sub>2</sub></sub>	VT <sub>CO<sub>2</sub></sub>	VT <sub>CO<sub>2</sub></sub>	VT <sub>CO<sub>2</sub></sub>			
V <sub>ds</sub> /V <sub>Te</sub>	V <sub>Te</sub>	V <sub>Te</sub>	V <sub>Te</sub>			
Slope phase III	Slope phase III	Slope phase III	Slope phase III			

\*depending on the device used

1 St. John RE.: End-tidal carbon dioxide monitoring Crit Care Nurs Vol 23, No. 4, August 2003; 83-88

2 St. John RE.: Exhaled gas analysis: technical and clinical aspects of capnography and oxygen consumption. Crit Care Nurs Clin N Am. 1989; 20:363-374

3 Bongard F, Sue D.: Pulse oximetry and capnography in intensive and transitional care units. West J. Med. 1992 Jan; 156(1); 57-64

4 AARC Guideline: Capnography/Capnometry during Mechanical Ventilation-2003 revision and update: Respiratory Care, May 2003 Vol. 48 No. 5

5 Behende et al.: Validity of a disposable and end-tidal CO<sub>2</sub> detection in verifying endotracheal tube placement in infants and children. Ann Erg Med 1992 31:142-5

”A volumetric capnogram contains extensive physiological information about metabolic production, circulatory transport and CO<sub>2</sub> elimination within the lungs. VCap\* is also the best clinical tool to measure dead spaces allowing a detailed analysis of the functional components of each tidal volume, thereby providing clinically useful hints about the lung’s efficiency of gas exchange.

Recent advances in VCap and our improved understanding of its clinical implications may help in overcoming the known limitations and reluctances to include expired CO<sub>2</sub> kinetics and dead space analysis in routine bedside monitoring. It is about time to start using this powerful monitoring tool to support decision making in the intensive care environment”.

Suarez-Sipmann F, et al; Volumetric capnography: the time has come., 2014 Jun;20(3):333-9

\*Vcap: Volumetric capnography

## ACCESSORIES

- Mainstream CO<sub>2</sub> sensor that measures continuously
- Reusable and disposable CO<sub>2</sub> cuvettes for adult and pediatric patients
- The patented design of the disposable CO<sub>2</sub> cuvette delivers the same performance as the reusable type, but without the additional cost of time consuming sterilization procedures and protocols



Mainstream CO<sub>2</sub> sensor for Evita Infinity V500, Evita V300, Savina 300

adult  
disposable



pediatric



reusable

Mainstream CO<sub>2</sub> measurement

End-tidal CO <sub>2</sub> concentration	etCO <sub>2</sub>
Range	0 to 100 mmHg or 0 to 13.2 Vol% (at 1013 mbar (or cmH <sub>2</sub> O)) or 0 to 13.3 kPa
Accuracy	According to ISO 80601-2-55 ±(3.3 mmHg +8 % of the measured value) ±(0.44 kPa +8 % of the measured value) ±(0.43 Vol% +8 % of the measured value)
Measurement conditions	Respiratory rate: 6 to 100 /min Inspiratory time: >250 ms Expiratory time: >250 ms
Drift of measurement accuracy	According to ISO 80601-2-55 <0.2 kPa (at 5.00 kPa) over 6 h
The measured values of the CO <sub>2</sub> measurement are barometrically pressure-compensated.	
T10...90	<35 ms
Response time, total	<200 ms
Time until specified accuracy is reached	<120 s (at 23 ±2 °C)
With reference to the displayed measured values the following dead spaces must be taken into account:	
CO <sub>2</sub> cuvette, adults (6870279, MP01062)	4.3 mL
CO <sub>2</sub> cuvette, pediatric patients (6870280, MP01063)	1.9 mL
Neonatal flow sensor ISO 15 (8411130)	0.9 mL
Neonatal flow sensor Y-piece (8410185)	1.7 mL

## CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

## Manufacturer:

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23542 Lübeck, Germany

## REGION DACH

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

## REGION EUROPE

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

## REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

## REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.  
25 International Business Park  
#04-20/21 German Centre  
Singapore 609916  
Tel +65 6308 9400  
Fax +65 6308 9401  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

## REGION NORTH AMERICA

Draeger, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 800 4DRAGER  
(+1 800 437 2437)  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

## REGION CENTRAL AND SOUTH AMERICA

Dräger Panama S. de R.L.  
59 East Street, Nuevo Paitilla,  
House 30, San Francisco Town  
Panama City, Panama  
Tel +507 377 9100  
Fax +507 377 9130  
[servicioalcliente@draeger.com](mailto:servicioalcliente@draeger.com)

Locate your Regional Sales Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

