


Dräger N95 Respirator Donning Instructions for Community Use



1. Normal workplace use requires each user to be tested to assure proper seal on the face (fit test). In the absence of a fit test, a user seal check is required to achieve an acceptable face-seal.



2. If you have long hair, it should be tied in a position that does not interfere with the head straps.



3. Position the mask under your chin and over your nose. Slip the first strap over your head to the back of your neck.



4. Pull the upper strap into position on the crown at the back of your head, ensure the strap sits above your ears.



5. Straps should not be twisted or crossed over.



6. Adjust the nose clip to minimize leakage and make sure the nose clip has a smooth fit around the nose and there are no gaps.



7. Cover the mask with both hands and exhale. The mask should inflate slightly.



8. If any air escapes at the side of the mask, readjust the mask until it fits properly and ensure that the nose clip is positioned correctly.



[Learn more about Dräger N95 Respirators](#)

*If you want to protect yourself and others, use the respirator without an exhalation valve.