Infinity Acute Care System

de	Gebrauchsanweisung, seite 2	tr	Kullanma Kılavuzu, sayfa 17
en	Instructions for Use, page 3	hu	Használati útmutató, 18. oldal
fr	Notice d'utilisation, page 4	sk	Návod na použitie, strana 19
es	Instrucciones de uso, página 5	cs	Návod k použití, strana 20
it	Istruzioni per l'uso, pagina 6	zh	使用说明书, 第 21 页
pt_ib	Instruções de utilização, página 7	ko	사용 지침, 22 참조
pt_br	Instruções de uso, página 8	ja	取扱説明書, 23 ページ
nl	Gebruiksaanwijzing, pagina 9	ro	Instrucțiuni de utilizare, pagina 24
da	Brugsvejledning, side 10	el	Οδηγίες Χρήσης, σελίδα 25
no	Bruksanvisning, side 11		
sv	Bruksanvisning, sida 12		
fi	Käyttöohje, sivu 13		
ru	Руководство по эксплуатации, стр. 14		
pl	Instrukcja obsługi, strona 15		
lt	Naudojimo instrukcija, 16 psl.		

WARNING

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TC

20 ms

Supplement to the Instructions for Use for Babylog VN500 SW 2.n with part number 9038982

This supplement must be kept with the Instructions for Use of the medical device.

This supplement updates the information in the Instructions for Use of the medical device in the following chapters:

For your safety and that of your patients, section "Product-specific safety information"

CAUTION

To ensure the accuracy of the pressure measurements during HFO, it is necessary to perform a breathing circuit check to determine the resistance and compliance values. Otherwise the mean airway pressure could deviate from the set values.

Technical data, section "Displayed measured values, Airway pressure measurement"

Accuracy does only apply for the measurement range specified.

Airway pressure measurement

Positive end-expiratory pressure

PEEP

Peak Inspiratory Pressure

PIP

Mean airway pressure

Pmean

Minimal airway pressure

Pmin

Lower pressure level in APRV

Plow

End-inspiratory pressure for mandatory breaths

EIP

Upper pressure level in APRV

Phigh

Range

–60 to 120 mbar (or hPa or cmH₂O)

Accuracy

±6 % of measured value, or ±0.5 mbar (or hPa or cmH₂O), whichever is greater

TO...90 (for Pmean)

33 s for intubated patients and 20 s with NIV

Airway pressure measurement in HFO

Pressure amplitude (peak-to-peak) in HFO

ΔPhf

Range

0 to 90 mbar (or hPa or cmH₂O)

Accuracy

±20 % of measured value, or ±5.5 mbar (or hPa or cmH₂O), whichever is greater, if the breathing circuit was calibrated

Mean airway pressure under high frequency ventilation

Pmean

Range

0 to 50 mbar (or hPa or cmH₂O)

Accuracy

±15 % of measured value, or ±1.9 mbar (or hPa or cmH₂O), whichever is greater, if **Ampl hf** <50 mbar (or hPa or cmH₂O), if the breathing circuit was calibrated

Accuracy

±20 % of measured value, or ±2.9 mbar (or hPa or cmH₂O), whichever is greater, if **Ampl hf** ≥50 mbar (or hPa or cmH₂O), if the breathing circuit was calibrated

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© Dräger Medical GmbH
Dräger Medical GmbH
Drägerwerk AG & Co. KGaA