

Dräger Babylog® VN500 Cuidados abiertos neonatales

Para las futuras generaciones. El Babylog® VN500 combina los años de experiencia de Dräger con la última tecnología. El resultado es una solución de ventilación completa e integrada para los pacientes más pequeños. Supere nuevas barreras hoy para estar preparado para el mañana.



Ventajas

Configurable user interface and monitoring tools

- Individual monitoring views that can be determined by the user
 - Standardised, intuitive and user-friendly graphical user interface
 - Online help including context-sensitive help functions
 - Extended monitoring functions and smart data visualizations
-

Decision-making tools that decrease cognitive workload

- Smart Pulmonary View provides a graphical display of the compliance and resistance, including spontaneous breathing
 - Trending, measured parameters, waveforms and loops
-

Workstation functions

- Configuration to suit your needs
 - Downloadable screenshots for training, research and knowledge transfer
 - Fast, standard configuration of all Babylog VN500 devices via USB
 - Ability to connect the C500 Cockpit display to an overhead projector
 - Several log export options to support learning and research
-

Advanced respiratory care

- HFOV including 'sigh' breaths for lung recruitment and Volume Guarantee
- Original Dräger leak adaptation and leak compensation technology to maintain sensitive triggering and volume targets
- PC-MMV can be used for weaning and promotes spontaneous breathing and automatically adapts support in line with patient needs
- Integrated non-invasive ventilation and high-flow oxygen therapy

Accesorios



MT-706-2001

Gafas nasales


Gafas nasales XS, referencia 8418415 (10 unidades)

Gafas nasales S, referencia 8418605 (10 unidades)

Gafas nasales M, referencia 8418416 (10 unidades)

Gafas nasales L, referencia 8418531 (10 unidades)

Gafas nasales XL, referencia 8418417 (10 unidades)



D-7384-2009

Accesorios para cuidados neonatales

Los accesorios de Dräger para la ventilación y la termoterapia neonatales ventilan, reducen el estrés y ayudan a promover el desarrollo del recién nacido de manera cuidadosa y efectiva con una amplia gama de accesorios diseñados específicamente para su uso con los pacientes más pequeños.



MT-6197-2006

BabyFlow desechable

Referencia 8418583 (20 unidades)

Especificaciones técnicas

Patient type	Neonatal and pediatric patients
Ventilation settings	
Ventilation modes	Pressure-controlled ventilation: <ul style="list-style-type: none"> - PC-CMV - PC-AC - PC-SIMV - PC-PSV - PC-MMV - PC-HFO - PC-APRV Support of spontaneous breathing: <ul style="list-style-type: none"> - SPN-CPAP/PS - SPN-CPAP/VS - SPN-PPS - SPN-CPAP
Enhancements	<ul style="list-style-type: none"> - Apnea ventilation - Flow trigger - Sigh - Volume Guarantee (VG) - Smart Pulmonary View - Automatic Tube Compensation® (ATC) - AutoRelease (in PC-APRV) - HFO-Sigh - Volume Guarantee for HFO (HFO-VG) - Leak Compensation - O₂-therapy
Special procedures	<ul style="list-style-type: none"> - Suction maneuver - Manual inspiration/hold - Medication nebulization
Therapy types	<ul style="list-style-type: none"> - Invasive ventilation (Tube) - Non-invasive ventilation (NIV) - O₂-therapy
Respiratory rate (RR)	0.5 to 150/min
Inspiratory time (Ti)	0.1 to 3 s
Maximum inspiratory time for flow cycled breaths (T _{imax})	Neonates 0.1 to 1.5 s Pediatric patients 0.1 to 4 s
Tidal volume (VT)	Neonates 2 to 100 mL Pediatric patients 20 to 300 mL
Inspiratory flow (Flow)	Neonates 2 to 30 L/min Pediatric patients 2 to 30 L/min
Tidal volume during Apnea Ventilation (VT _{apn})	Neonates 2 to 100 mL Pediatric patients 20 to 300 mL
Respiratory rate during Apnea Ventilation (RR _{apn})	2 to 150/min
Inspiratory pressure (P _{INSP})	1 to 80 mbar (or hPa or cmH ₂ O)
Inspiratory pressure limit (P _{MAX})	2 to 100 mbar (or hPa or cmH ₂ O)
Positive end-expiratory pressure (PEEP)	0 to 35 mbar (or hPa or cmH ₂ O)
Rise time for pressure support (Slope)	Neonates 0 to 1.5 s Pediatric patients 0 to 2 s
O ₂ concentration (FiO ₂)	21 to 100 Vol%
Trigger sensitivity (Flow trigger)	0.2 to 5 L/min

Especificaciones técnicas

Automatic Tube Compensation® (ATC)	<p>Inner diameter of the tube (Tube Ø)</p> <ul style="list-style-type: none"> - Endotracheal tube (ET) <ul style="list-style-type: none"> Pediatric patients 2 to 8 mm (0.08 to 0.31 in) Neonates 2 to 5 mm (0.08 to 0.2 in) - Tracheostomy tube (Trach.) <ul style="list-style-type: none"> Pediatric patients 2.5 to 8 mm (0.1 to 0.31 in) - Degree of compensation 0 to 100%
High Frequency Oscillation (PC-HFO)	<ul style="list-style-type: none"> - Mean airway pressure (MAPhf) 5 to 50 mbar (or hPa or cmH₂O) - Frequency of oscillation (fhf) 5 to 20 Hz - I to E (I:Ehf) 1:1 to 1:3 - Pressure amplitude (Ampl hf) 5 to 90 mbar (or hPa or cmH₂O) - Maximum pressure amplitude (Ampl hf max) in HFO (VG) 5 to 90 mbar (or hPa or cmH₂O) - Tidal volume (VThf) 0.2 to 40 mL - Sigh pressure (Psigh) 6 to 80 mbar (or hPa or cmH₂O) - Respiratory rate of sigh (RRsigh) 0 to 30/min - Sigh pressure rise time (Slope sigh) <ul style="list-style-type: none"> Pediatric patients 0 to 2 s, Neonates 0 to 1.5 s - Sigh inspiratory time (Tisigh) 0.1 to 3 s
Leakage compensation	<ul style="list-style-type: none"> - On / Off - On: full compensation active - Off: trigger compensation active
O ₂ -therapy	<p>Continuous Flow 2 to 50 L/min</p> <p>O₂ concentration FiO₂ 21 to 100 Vol%</p>
Displayed measured values	
Airway pressure measurement	<p>Positive end-expiratory pressure (PEEP)</p> <p>Peak inspiratory pressure (PIP)</p> <p>Mean airway pressure (Pmean)</p> <p>Minimum airway pressure (Pmin)</p> <p>Lower pressure level in APRV (Plow)</p> <p>End-inspiratory pressure for mandatory breaths (EIP)</p> <p>Upper pressure level in APRV (Phigh)</p> <p>Range -60 to 120 mbar (or hPa or cmH₂O)</p>
Flow measurement (proximal)	
Minute volume measurement	<p>Expiratory minute volume (MVe)</p> <p>Inspiratory minute volume (MVi)</p> <p>Mandatory expiratory minute volume (MVemand)</p> <p>Spontaneous expiratory minute volume (MVspan)</p> <p>Spontaneous expiratory minute volume (MV)</p> <p>Range 0 to 30 L/min BTPS</p>
Tidal volume measurement	<p>Tidal volume (VT)</p> <p>Inspiratory tidal volume (not leakagecompensated) of mandatory breaths (VTimand)</p> <p>Expiratory tidal volume (not leakagecompensated) of mandatory breaths (VTemand)</p> <p>Inspiratory tidal volume (not leakagecompensated) of spontaneous breaths (VTispon)</p> <p>Range 0 to 1,000 mL BTPS</p>

Especificaciones técnicas

Respiratory rate measurement	Respiratory rate (RR) Mandatory respiratory rate (RRmand) Spontaneous respiratory rate (RRspon) Range 0 to 300/min
O ₂ measurement (inspiratory side)	Inspiratory O ₂ concentration (in dry air) (FiO ₂) Range 18 to 100 Vol%
CO ₂ measurement in mainstream (pediatric patients only)	End-expiratory CO ₂ concentration (etCO ₂) Range 0 to 100 mmHg
Displayed calculated values	
Compliance (C)	Range 0 to 650 mL/mbar
Resistance (R)	Range 0 to 1,000 mbar/(L/s)
Leakage minute volume (MVleak)	Range 0 to 30 L/min BTPS
Rapid Shallow Breathing (RSB)	Range Neonates 0 to 300 (/min/mL) Range Pediatric patients 0 to 9999 (/min/L)
Curve displays	Airway pressure Paw (t) -30 to 100 mbar Flow (t) -40 to 40 L/min Volume V (t) 2 to 300 mL CO ₂ (t) 0 to 100 mmHg
Alarms/Monitoring	
Expiratory minute volume (MVe)	High/Low
Airway pressure (Paw)	High/Low
Insp. O ₂ concentration (FiO ₂)	High/Low (automatic)
End-expiratory CO ₂ concentration (etCO ₂)	High/Low
Respiratory rate (RR)	High
Volume monitoring (VT)	Low (automatic)
Apnea alarm time (T _{apn})	5 to 60 seconds, Off
Disconnect alarm delay time (T _{disconnect})	0 to 60 seconds
Performance characteristics	
Control principle	time-cycled, pressure-controlled, volume-constant
Inspiratory flow (BTPS)	max. 60 L/min
Base flow, neonates	6 L/min
Base flow, pediatric patients	3 L/min
Safety valve	Opens if medical compressed air supply fails (supply gas flow is not sufficient to provide the inspiratory flow required), enables spontaneous breathing with ambient air.
Endotracheal suction	
Disconnection detection	automatic
Reconnection detection	automatic
Initial oxygen enrichment	max. 3 minutes
Active suction phase	max. 2 minutes
Final oxygen enrichment	max. 2 minutes
Oxygen enrichment for suction maneuver	Factor for neonates 1 to 2 Factor for pediatric patients 1 to 2
Maneuver settings	
Sigh pressure (Δ _{int} PEEP)	0 to 20 mbar (or hPa or cmH ₂ O)
Time interval between sighs (Interval sigh):	20 s to 180 min
Number of cycles for a sigh (Cycles sigh):	1 to 20 exhalations
Medication nebulization	for 5,10,15,30 minutes
Control principle	time-cycled, pressure-controlled, volume-constant

Especificaciones técnicas

Inspiratory flow (BTPS)	Pediatric patients 2 to 30 L/min Neonates 2 to 30 L/min
Operating data	
Mains power supply	
Mains power connection	100 V to 240 V, 50/60 Hz
Current consumption	
at 230 V	max. 1.4 A
at 230 V	max. 3.0 A
Inrush current	approx. 8 to 24 A peak approx. 6 to 17 A quasi RMS
Power consumption	
maximum	300 W
during ventilation, without charging the battery	approx. 100 W ventilation unit with Medical Cockpit approx. 180 W with GS500
Digital machine output	Digital output and input via an RS232 C interface Dräger MEDIBUS and MEDIBUS.X
Gas supply	
O ₂ gauge pressure	2.7 to 6.0 bar (or 270 to 600 kPa or 39 to 87 psi)
Air gauge pressure	2.7 to 6.0 bar (or 270 to 600 kPa or 39 to 87 psi)
Physical Specifications	
Dimensions (W x H x D)	
Babylog® VN500 and Infinity® C500	420 mm x 685 mm x 410 mm (16.5 in x 27.0 in x 16.1 in)
Babylog® VN500 and Infinity® C500 on trolley	577 mm x 1,420 mm x 687 mm (22.7 in x 55.1 in x 26.7 in)
Weight	
Babylog® VN500 and Infinity® C500	approx. 25 kg (55.1 lbs)
Babylog® VN500 and Infinity® C500 on trolley	approx. 59 kg (130 lbs)
GS500	approx. 10.5 kg (23 lbs)
PS500	approx. 27 kg (59.5 lbs)
Adapter for 38 mm pole	approx. 2.35 kg (5.18 lbs)
Infinity® C500	
Diagonal screen size	17" TFT color touch screen
Input / Output ports	<ul style="list-style-type: none"> - RS232 (9-pin) connectors - USB ports for data collection - 1 DVI for digital video output - RJ 45 Ethernet connectors (service purposes only)
<p>ATC®, trademarked by Dräger. AutoFlow®, trademarked by Dräger</p> <p>BTPS – Body Temperature Pressure Saturated. Measured values relating to the conditions of the patient lung (98.6 °F), steamsaturated gas, ambient pressure.</p> <p>1 mbar = 100 Pa</p>	
Some functionalities are available as an option.	

Observaciones

No todos los productos, características o servicios están disponibles para la venta en todos los países. Las marcas comerciales citadas están registradas en ciertos países únicamente y no necesariamente en el país en el que se publique este material. Visite www.draeger.com/trademarks para conocer el estado actual.

SEDE PRINCIPAL

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Alemania
www.draeger.com

Fabricante:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23542 Lübeck, Alemania

Localice a su representante
de ventas regional en:
www.draeger.com/contacto



ARGENTINA

Dräger Argentina S.A.
Colectora Panamericana Este
1717B, 1607BLF San Isidro,
Buenos Aires
Tel +54 11 48 36 8300 / Fax -8321

BRASIL

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06460-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

CHILE

Dräger Chile Ltda.
Av. Presidente Eduardo
Frei Montalva 6001-68
Complejo Empresarial El Cortijo,
Conchalí, Santiago
Tel +56 2 2482 1000 / Fax -1001

COLOMBIA

Draeger Colombia S.A.
Oficina: Calle 100 #13-21
Oficina 503 Bogotá D.C.
Tel +60 17 94 50 50

ESPAÑA

Dräger Hispania, S.A.
C/ Xaudaró 5, 28034 Madrid
Tel +34 90 011 64 24
Fax +34 91 358 36 19
atencionalcliente@draeger.com

MÉXICO

Dräger Medical México,
S.A. de C.V., German Centre
Av. Santa Fe, 170 5-4-14
Col. Lomas de Santa Fe
01210 México D.F.
Tel +52 55 52 61 43 37
Fax +52 55 52 61 41 32

PERÚ

Draeger Perú SAC
Av. San Borja Sur 573-575
Lima 41
Tel +511 626 95-95 / Fax -73

PORTUGAL

Dräger Portugal, Lda.
Rua Nossa Senhora da
Conceição, n.º 3, R/c
2790-111 Carnaxide
Tel +351 21 155 45 86
Fax +351 21 155 45 87
clientesportugal@draeger.com