


Dräger DrugCheck® 3000 On-the-spot Drug Screening Device

Use the Dräger DrugCheck® 3000 to find out within minutes if a person has recently consumed certain drugs. This compact oral fluid-based drug test is for presumptive or preliminary screening. The device does not require electricity and can be used practically anywhere.

In the United States, the DrugCheck 3000 is not a medical device, is not cleared by the United States Food and Drug Administration as a medical device, and is only made available for law enforcement use and purposes.



Benefits

Obtain samples easily and safely

The Dräger DrugCheck® 3000 is a single use screener, consisting of two components: an absorbent collector (swab) for obtaining an oral fluid sample and a test cassette for the analysis. The test cassette contains a liquid buffer solution and a window with two test strips that display the control and test lines. The test can be performed in a few easy steps: 1) take a saliva sample using the swab, 2) insert the swab into the funnel opening of the test cassette, 3) shake the test kit, 4) begin the test analysis, 5) read the results shown in the display window. A color indicator on the swab disappears as soon as it has absorbed enough oral fluid for a test.

As soon as the control lines appear in the window, you can read the results. If the test result is negative, a line will appear alongside the respective substance class (drug). This means that none of the target substance was detected in the sample. If a line fails to appear next to one of the substance classes, then the result for this substance is positive.

Fast and sensitive analysis

Check individuals for up to six substance classes with the Dräger DrugCheck 3000: cocaine, opiates, amphetamines, methamphetamines/designer drugs (e.g. Ecstasy, MDMA), benzodiazepines (e.g. in prescription medicines), and cannabis (THC). Of all these substance classes, cannabis is the drug consumed most frequently – and also the most difficult of all compounds to identify (THC = delta-9- tetrahydrocannabinol). This is why Dräger optimized the DrugCheck 3000 to detect THC – and now offers two measurement options. The fast option (25 ng/ml THC threshold) takes less than 3 minutes and the sensitive option (15 ng/ml THC threshold) takes less than 5 minutes to show a positive or negative result on the test window.

Convenient, hygienic on-the-spot drug detection

The pocket-sized design makes the DrugCheck 3000 easy to transport. It can be made ready for use quickly and easily. The DrugCheck® 3000 test kit is tamper-proof and hygienic to use.

Maintenance Free

Because the DrugCheck 3000 has no electrical parts, it's safe to use in hazardous areas. The portable test kit can be easily disposed of after use and there is no calibration or servicing required - which makes using the device a hassle-free experience.

Backed by decades of experience

Dräger has decades of experience in alcohol measuring equipment and drug detection methods. For the DrugCheck® 3000, Dräger employed the testing principle of the proven Dräger DrugTest® 5000 system, which reliably detects even minute traces of THC.

Benefits

Dräger DrugCheck® App

Dräger DrugCheck® App: drug detection just got even more convenient and efficient. Scan, test and record the data: mobile data logging on your smartphone is that simple. The new Dräger DrugCheck® App allows you to scan, analyze, log and archive your test results in just a few steps. Video instructions guide you safely through the process of handling the DrugCheck 3000, analysing the test result.

Infographic

| Dräger DrugCheck® 3000 2 Measurement Options for Detecting THC | |
|--|--|
| FAST  | SENSITIVE  |
| 25ng/ml | 15ng/ml |
| Results in < 3 minutes | Results in < 5 minutes |

When using the test kit, please reference the Quick IFU document for the correct procedure of detecting the desired THC threshold.

Services



[Dräger DrugCheck® App in the App Store / Google Play Store](#)

Scan the QR code to go directly to the download page for the Dräger DrugCheck® App in the App Store / Google Play Store

Related Products



Dräger DrugTest® 5000

No pipetting, no drips, no timing: With the Dräger DrugTest® 5000, drug test is carried out simply and quickly. The collected saliva sample can be analyzed immediately—for accurate results on the spot.



Dräger SSK 5000

Ideal for surface drug testing applications, the Dräger SSK 5000 sampling system is quick and easy to use. When used with the Dräger DrugTest® 5000, the sample can be transferred to a test cassette and screened for possible drug contamination. The sample can also be stored and transported for further testing.

Related Products



D-4542E-2015

Dräger Alcotest® 5820

The Dräger Alcotest® 5820 eliminates guesswork by providing a reliable breath alcohol concentration reading every time. The device ensures precise results by using the same measurement technology used by police in over 30 million breath alcohol tests a year.



D-30476-2017

Dräger Alcotest® 5000

The Dräger Alcotest® 5000 is a professional breath-alcohol tester which detects the presence of alcohol. The high-speed breathalyzer lets you perform numerous tests in minimal time. Its special funnel reduces the back-flow of expired air to a minimum, preventing the risk of infection to subsequent test subjects.

Technical Data

Test for detecting drugs in oral fluid

| | |
|--------------------------------|--|
| Dimensions (W x H x D) | 32 x 111 x 57 mm 1.26 x 4.37 x 2.24 inches |
| Weight | <30 g (0.07 lb) |
| Operating range | Operation: 5 °C to 40 °C (41 °F to 104 °F) at 5% to 95% relative humidity |
| Storage/Transportation | 4 °C to 25 °C (39 °F to 77 °F) |
| Time of measurement | Fast measurement: analysis <3 minutes Sensitive measurement: analysis <5 minutes |
| Selecting the measurement mode | Depending on the selected preincubation time and the desired THC cutoff, fast or more sensitive reading. |
| Licenses/Standards | Outside of the European Union for non-medical use (Non-IVD) or for forensic use (only USA). |

Ordering Information

| Dräger DrugCheck® 3000 | Description | Ordering no. |
|---|------------------------|--------------|
| Test for detecting drugs in oral fluid | | |
| Dräger DrugCheck® 3000 STK 6-Non-IVD | Packaging unit 20 kits | 83 27 960 |
| Dräger DrugCheck® 3000 detects the following substances: | | |
| <ul style="list-style-type: none"> - Cocaine (COC) - Opiates (OPI) - Cannabis (THC) - Amphetamines (AMP) - Methamphetamines (MET) - Benzodiazepines (BZO) | | |

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS
Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Allemagne
www.draeger.com

USA
Dräger, Inc.
7256 S. Sam Houston Parkway W.,
Suite 100
Houston, TX 77085
1 800 4DRAGER
(1 800 437 2437)

Locate your Regional
Sales Representative at:
www.draeger.com/contact

