


Quick start guide – HI-Flow Star for adults


1a Connection to a blender

1. Connect the O₂ hose (A) to the blender (B) and the water chamber (C).
2. Connect the heated inspiratory hose to the second opening of the water chamber.
3. Attach the heating cable adaptors and temperature sensors to the humidifier and hose system.




1b Connecting to a ventilator

1. Connect the inspiratory hose (A) to the inspiratory port of the ventilator (B) and to the water chamber (C).
2. Connect the heated inspiratory hose to the second opening of the water chamber.
3. Connect the heating cable adaptors and temperature sensors to the humidifier and hose system.



2

1. Use the sizing gauge to establish the size of the prong required. The sizing gauge shows the prong sizes available.
2. The prong should only close the nostrils to a maximum of 50%.
3. Connect the nasal cannula to the inspiratory hose.
4. Start the device as per the instructions for use and set the parameters. Check that there is a respiratory flow coming out of the prong.



3

1. The nasal cannula hose may run to the left or right of the patient.
2. Turn the nasal cannula so that the hose from the nasal cannula runs in the direction of the device.
3. Carefully place the prong in the patient's nose and put the headgear on.
4. The size of the headgear can be adjusted using the hook and loop fasteners.

