# Supplement (MS32633)

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Tilleggsdokument for Infinity Acute Care System VG4.n – Norsk – side 27

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## WARNING

For a full understanding of the performance characteristics of this device, the user should carefully read this supplement and the related instructions for use before use of the device.

## Infinity Acute Care System VG4.n

This page has been left blank intentionally.

## Supplement for Infinity Acute Care System VG4.n

---

This document supplements the instructions for use of Infinity Acute Care System (IACS) Monitoring Applications Software VG4.n, and Infinity M540 patient monitor Software VG4.n.

### Compatible Medibus.X devices

---

The Cockpit device connectivity interface allows data from various standalone devices to display parameters, settings, and waveforms on the Cockpit.

The following table lists which devices and corresponding software versions are supported with Medibus.X version 1.0.3.

Supported device	Supported software version
<b>Ventilators</b>	
Dräger Evita V500	2.31
Dräger Babylog VN500	2.31
Savina 300	4.02
Carina	3.21
Dräger Evita V300	2.31
Oxylog 3000+	1.04
<b>Anesthesia machines</b>	
Primus family, Apollo	4.5
Infinity Perseus A500	1.11
Zeus IE	1.04
Fabius family	3.35

## ST monitoring specifications

---

In the ST segment analysis table, the updated ST input accuracy is as follows:

### ST segment analysis

ST input accuracy	<p><math>\pm 0.5</math> mm (.05 mV) or 15% of the measured value is valid for STI, STII, STIII, aVR, aVL, aVF, V, V+, V1 to V6, dV1, dV3, dV4 and dV6.</p> <p><math>\pm 3.2</math> mm (0.32 mV) for STVM and STCVM.</p> <p><b>CAUTION</b></p> <p><b>ST accuracy may be impacted if ECG double-counting occurs or with highly irregular heart rates.</b></p> <p><b>NOTE</b></p> <p>In an unlikely scenario where the heart rate is low (e.g., 30 bpm, 45 bpm), and ST-segment deviation is high (<math>\geq \pm 13</math> mm), the claimed ST accuracy might not be met on one or two leads.</p>
-------------------	---

## Security recommendations

---

Dräger makes the following security recommendations:

- Physical security of the patient monitors is recommended and is the responsibility of the operating organization.
- Physical security of the telecommunications closet is recommended and is the responsibility of the operating organization.
- Dräger recommends that operating organizations restrict physical access to unused Ethernet ports on the IACS.
- Dräger recommends that operating organizations restrict physical access to unused USB and serial ports on the IACS.
- Dräger relies on the medical device isolation mechanism of the VLANs and the proper configuration, implementation, and use of the operating organization's security measures to prevent the introduction of malware onto the Infinity network.

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These Instructions for Use only apply to  
**Infinity Acute Care System VG4.n**  
with the Serial No.:

If no Serial No. has been filled in by Dräger,  
these Instructions for Use are provided for  
general information only and are not intended for  
use with any specific machine or unit.


This document is provided for customer  
information only, and will not be updated or  
exchanged without customer request.

Directive 93/42/EEC  
concerning Medical Devices

 Manufacturer

 **Draeger Medical Systems, Inc.**  
3135 Quarry Road  
Telford, PA 18969-1042  
U.S.A.  
 (215) 721-5400  
(800) 4DRAGER  
(800 437-2437)  
FAX (215) 723-5935  
 <http://www.draeger.com>

 In Europe, Middle East, Africa, Latin  
America, Asia Pacific distributed by

 **Drägerwerk AG & Co. KGaA**  
Moislinger Allee 53 – 55  
23542 Lübeck  
Germany  
 +49 451 8 82-0  
FAX +49 451 8 82-20 80  
 <http://www.draeger.com>

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Edition: 3 – 2018-04

(Edition: 1 – 2015-12)

Dräger reserves the right to make modifications  
to the equipment without prior notice.

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© Drägerwerk AG & Co. KGaA