

## ComfortStar® Disposable anesthesia face masks

More patient comfort and easy handling.



Patient comfort and ease of use were the top priorities when developing the Dräger disposable ComfortStar anesthesia face masks.

- Patient comfort is supported by a teardrop or a round-shaped mask, which provides an anatomical fit with effective seal.
- Non-slip ridges and a flexible mask body provide healthcare professionals an easy and comfortable grip on the mask.
- All ComfortStar masks are compatible with all Dräger anesthesia devices.



MT-4197-2008

**ComfortStar® masks** for adult application (sizes 4-6) are equipped with ridges on the mask shell to provide a non-slip grip.



MT-1176-2005

**Two-sided pump ball** for air management to inflate/deflate mask cushion.



MT-3938-2008

**Color coded hook rings** for an easier size identification.

## ComfortStar® anesthesia face masks are for disposable use

The adjustable mask cushion combined with the round or teardrop shape provides anatomical fit and thus more patient comfort and an effective seal.

The inflatable mask cushion and the flexibility of the mask give the physician a comfortable grip. In addition, the ComfortStar masks sizes 4-6 are equipped with slip-stop ridges to provide a safe and non-slip grip.

The transparent mask body allows for continuous observation of the patient's condition. All masks are pre-equipped with 'Traffic Light' colored hook rings. They are color coded for easier size identification. Smaller sizes (0-1) have reddish color, mid-sizes (2-3) yellowish color and large sizes (4-6) greenish color. Small and large headgears are also available for easy and secure fixation.


A potential side effect for patients who have undergone surgical procedure under general anesthetic is post-operative nausea and vomiting (PONV). This often delays patient recovery. Research shows that added aromas can help minimizing these unpleasant side effects<sup>1)</sup>.

ComfortStar scented anesthesia face masks have the same excellent benefits of the ComfortStar standard masks, but in addition have added aromas.

They are available with strawberry or bubble gum scents in sizes 1-4, and mint scent in sizes 3-6. The appropriate color scheme provides optical stimulation and is intended to reduce pre-operative anxiety in patients.


<sup>1)</sup> Anderson, L.A., Gross, J.B.:  
Aromatherapy with peppermint, isopropyl alcohol, or placebo is equally effective in relieving postoperative nausea, 2004, in: Journal of PeriAnesthesia Nursing, Vol. 39, No 1, 29-55

## ComfortStar® anesthesia face masks and accessories




MT-4205-2008

ComfortStar® standard masks, sizes 0-6



MT-4170-2008

ComfortStar® strawberry scented masks, sizes 1-4



MT-4126-2008

ComfortStar® bubble gum scented masks, sizes 1-4




MT-4149-2008

ComfortStar® mint scented masks, sizes 3-6




MT-2859-2006

Headgear, Silicone, reusable, small and large



MT-9656-2006

Headgear, Soft Fabric, disposable, small and large



MT-9654-2006

Headgear, Soft Rubber, disposable, small and large



MT-11178-2006

Two-sided pump ball for air management

## ORDER INFORMATION

Description	Package of	Order No.	Description	Package of	Order No.
<b>A</b> ComfortStar anesthesia face mask, disposable, standard, with hook ring			<b>D</b> ComfortStar anesthesia face mask, disposable, mint scented, with hook ring		
<b>A0</b> size 0, neonate	20 pcs.	MP01510	<b>D3</b> size 3, toddler	20 pcs.	MP01543
<b>A1</b> size 1, new born	20 pcs.	MP01511	<b>D4</b> size 4, small adult	20 pcs.	MP01544
<b>A2</b> size 2, infant	20 pcs.	MP01512	<b>D5</b> size 5, medium adult	20 pcs.	MP01545
<b>A3</b> size 3, toddler	20 pcs.	MP01513	<b>D6</b> size 6, large adult	20 pcs.	MP01546
<b>A4</b> size 4, small adult	20 pcs.	MP01514			
<b>A5</b> size 5, medium adult	20 pcs.	MP01515	<b>E</b> Headgear, Silicone, reusable		
<b>A6</b> size 6, large adult	20 pcs.	MP01516	<b>E1</b> small	1 pc.	MP01557
			<b>E2</b> large	1 pc.	MP01558
<b>B</b> ComfortStar anesthesia face mask, disposable, strawberry scented, with hook ring			<b>F</b> Headgear, Soft Fabric, disposable		
<b>B1</b> size 1, new born	20 pcs.	MP01521	<b>F1</b> small	1 pc.	MP01559
<b>B2</b> size 2, infant	20 pcs.	MP01522	<b>F2</b> large	1 pc.	MP01560
<b>B3</b> size 3, toddler	20 pcs.	MP01523			
<b>B4</b> size 4, small adult	20 pcs.	MP01524	<b>G</b> Headgear, Soft Rubber, disposable		
			<b>G1</b> small	1 pc.	MP01561
<b>C</b> ComfortStar anesthesia face mask, disposable, bubble gum scented, with hook ring			<b>G2</b> large	1 pc.	MP01562
<b>C1</b> size 1, new born	20 pcs.	MP01531			
<b>C2</b> size 2, infant	20 pcs.	MP01532	<b>H</b> Two-sided pump ball for air management	1 pc.	MP01590
<b>C3</b> size 3, toddler	20 pcs.	MP01533			
<b>C4</b> size 4, small adult	20 pcs.	MP01534			

**HEADQUARTERS**  
 Draegerwerk AG & Co. KGaA  
 Moislinger Allee 53–55  
 23558 Lübeck, Germany

www.draeger.com

**REGION EUROPE CENTRAL  
 AND EUROPE NORTH**  
 Draeger Medical GmbH  
 Moislinger Allee 53–55  
 23558 Lübeck, Germany  
 Tel +49 451 882 0  
 Fax +49 451 882 2080  
 info@draeger.com

**REGION EUROPE SOUTH**  
 Draeger Médical S.A.S.  
 Parc de Haute  
 Technologie d'Antony 2  
 25, rue Georges Besse  
 92182 Antony Cedex, France  
 Tel +33 1 46 11 56 00  
 Fax +33 1 40 96 97 20  
 dlmfr-contact@draeger.com

**REGION MIDDLE EAST, AFRICA,  
 CENTRAL AND SOUTH AMERICA**  
 Draeger Medical GmbH  
 Branch Office Dubai  
 Dubai Healthcare City, P.O. Box 505108  
 Dubai, United Arab Emirates  
 Tel + 971 436 24 762  
 Fax + 971 436 24 761  
 contactuae@draeger.com

**REGION ASIA / PACIFIC**  
 Draeger Medical  
 South East Asia Pte Ltd  
 25 International Business Park  
 #04-27/29 German Centre  
 Singapore 609916, Singapore  
 Tel +65 6572 4388  
 Fax +65 6572 4399  
 asia.pacific@draeger.com

**REGION NORTH AMERICA**  
 Draeger Medical, Inc.  
 3135 Quarry Road  
 Telford, PA 18969-1042, USA  
 Tel +1 215 721 5400  
 Toll-free +1 800 437 2437  
 Fax +1 215 723 5935  
 info.usa@draeger.com

Manufacturer:  
 Draeger Medical GmbH  
 23542 Lübeck, Germany  
 The quality management system at  
 Draeger Medical GmbH is certified  
 according to ISO 13485, ISO 9001  
 and Annex II.3 of Directive 93/42/EEC  
 (Medical devices).