

Drägersorb® 800+

Soda lime is essential for CO₂ absorption in inhalation anaesthesia machines with rebreathing systems.



MT-3887-2008

Dräger is the only manufacturer of anaesthesia machines that develops and produces its own soda lime. Therefore, Drägersorb® 800+ addresses the specific needs of today's anaesthesia technology.

- High CO₂ absorption capacity due to uniform utilization, pellet form prevents channel formation, longer lifetime reduces costs
- Little dust generation creates a safer environment for your personnel and protects your equipment
- Additional safety for patients, few undesirable reaction products if accidentally dried out
- Customer-friendly packaging ensures easy handling

DRÄGERSORB® 800+ – ECONOMIC EXCELLENCE

Thanks to its unique pellet form, Drägersorb® 800+ does not undergo channel formation. Our soda lime absorbs CO₂ evenly and is effectively utilized. Its pellet form and high porosity create an optimal surface. Drägersorb® 800+ is therefore able to absorb more CO₂ and lasts longer than conventional soda lime.

DRÄGERSORB® 800+ – CREATING A BETTER ENVIRONMENT

The shape and composition of the Drägersorb® 800+ pellets enhance its strength and resistance to abrasion. This helps protect the lungs of both the patient and anaesthesia staff. Dust emissions are harmful to high-tech anaesthesia machines which is another reason why Drägersorb® 800+ was developed – to protect your valuable investments and reduce ongoing costs.



MT-1988-2008

Disposable absorber

MT-562-99



5 L - Container

MT-2004-2008



Disposable Absorber

D-566-2009

Drägersorb® CLIC
Adapter

DRÄGERSORB® 800+ – ADDITIONAL SAFETY

Drägersorb® 800+ significantly reduces the generation of decomposition products if the soda lime accidentally dries out prior to use.¹ Drägersorb® 800+ may be used with all halogenated inhalation anaesthetic agents. Of course, only fresh soda lime guarantees optimal use for patients. A change in color from white to violet blue clearly indicates when the soda lime needs to be changed.

DRÄGERSORB® 800+ – USER-FRIENDLY PACKAGING

The Drägersorb® CLIC disposable absorber for use with all modern Dräger anaesthesia systems delivers a quick, clean and uncomplicated disposable solution.

¹ Results confirmed by the independent institution Spectral Service Köln

Advantages of this system include:

- even greater convenience for anaesthesia staff, no contact with soda lime
- absolutely no dust contact
- more complete utilization reduces costs; change absorber at any time, even during surgery, thanks to the non-return valve.

The Drägersorb® CLIC disposable absorber for Dräger anaesthesia systems is attached to the usual absorber holder by means of an adapter.

Drägersorb® 800+ is also available in convenient 5 L containers.

ORDER INFORMATION

Description	Order No.
5 L - Container	MX00001
Drägersorb® 1.2 L CLIC 800+ – Disposable Absorber (6 per package)	MX00004
Drägersorb® CLIC Adapter	MX50090

Additional CLIC accessories on request.

HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlimfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA

Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel + 971 436 24 762
Fax + 971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:

Dräger Medical GmbH
23542 Lübeck, Germany
The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).