

## MonoLead ECG Lead-Wire Set

More time with your patients, less time with wires.




D-18305-2010

Frustrating tangles. Lost time. Patient discomfort. Managing the “spaghetti” that results from traditional ECG lead wires is a tedious, time-consuming distraction that takes focus away from your patient. As you attach the cables, you need to untangle and route the wires. Until now.

Dräger Medical, a leader in cable management, introduces a better way: the MonoLead™ ECG lead-wire set.

MonoLead is a patent-pending, one-wire ECG lead-set solution that directly replaces multi-wire systems. The single, tangle-proof lead set accommodates standard electrodes and electrode placement, while eliminating the tangled clutter that frequently occurs with traditional multi-lead sets. Ideal for both adult and pediatric patients, MonoLead is compatible with all Dräger Medical Infinity® patient monitors.



D-7684-2010

MonoLead 6 ECG Cable



## More comfortable for your patients, easier for you

Because MonoLead is a single cable rather than a complicated clutter of wires, it's less intimidating for your patient and their loved ones. Additionally it is more physically comfortable. That's because MonoLead is designed with Z-Snap™ zero-insertion force electrode snaps, which connect easily and hold securely, with no pressure or discomfort for your patient.

MonoLead is fast and easy to apply, and uses the standard lead placement that's familiar to nurses everywhere. It's flexible, durable and easy to clean. And the bright green color makes the cable easy to identify – even under a sheet.

So spend less time with wires and more time with your patients. Discover MonoLead, the one-wire ECG cable solution that's better for you, better for your patients.

### Consistent Volume Guarantee

- Tangle proof – replaces 3-, 5- and 6 electrode lead sets with a single wire
- Fast and easy to apply – has clearly labeled, color-coded snaps
- Same placement and results – supports the standard (Mason-Likar) electrode placement on the torso that is familiar to clinicians worldwide
- More comfortable – incorporates a patented Z-Snap™ zero-insertion force electrode snaps
- Suitable for both pediatric and adult patients
- Easy to clean – wipes clean with standard disinfectants

**ORDERING INFORMATION**

|   | 3368391; 5950196<br>MultiMed 5, 2,5m<br>Multimed 5, 1,5m | 5191221<br>MultiMed 6 | MS20093<br>MultiMed plus | MS20094<br>MultiMed plus OR | 5589663<br>MultiMed 12 |
|---|--|-----------------------|--------------------------|-----------------------------|------------------------|
| <b>Pods:</b>  |  |                       |                          |                             |                        |
| <b>MonoLeads T-connector</b>                          |  |                       |                          |                             |                        |
| MS14555 MonoLead 3 ECG CBL<br>T-connector IEC1        | X (a)  | X (a)                 | X                        | X                           |                        |
| MS14559 MonoLead 5 (limb) ECG CBL<br>T-connector IEC1 | X (a)  | X (a)                 | X                        | X                           |                        |
| MS14683 MonoLead 6 ECG CBL<br>T-connector IEC1        |  | X (a)                 | X                        | X                           |                        |
| <b>Adapters:</b>                                      |  |                       |                          |                             |                        |
| MS14679 Adapter MonoLead<br>to MultiMed 5             | X  |                       |                          |                             |                        |
| MS14680 Adapter MonoLead<br>to MultiMed 6             |  | X                     |                          |                             |                        |
| <b>HF-protection for OR</b>                           |  |                       |                          |                             |                        |
| 5947226 ESU block<br>for Multimed 5                   | X  |                       |                          |                             |                        |
| 7486140 ESU block<br>for Multimed 6                   |  | X                     |                          |                             |                        |
| <b>MonoLead Dualpin</b>                               |  |                       |                          |                             |                        |
| MS16160 MonoLead 3 ECG CBL<br>D/PIN IEC1              | X  | X                     |                          |                             |                        |
| MS16161 MonoLead 5 (limb) ECG CBL<br>D/PIN IEC1       | X  | X                     |                          |                             | X                      |
| MS16232 MonoLead 5 (chest) ECG CBL<br>D/PIN IEC1      |  |                       |                          |                             | X                      |

X (a) Adapter required

MT-5523-2004



MT-5530-2004



MT-5535-2004



Patented Z-Snap™ Zero-Insertion Force Electrode Snaps

**HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

Effective 2010-09:  
Dräger Medical AG & Co. KG  
changes to  
**Dräger Medical GmbH**

**REGION EUROPE CENTRAL  
AND EUROPE NORTH**

Dräger Medical AG & Co. KG  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION EUROPE SOUTH**

Dräger Médical S.A.S.  
Parc de Haute  
Technologie d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[d1mfr-contact@draeger.com](mailto:d1mfr-contact@draeger.com)

**REGION MIDDLE EAST, AFRICA,  
CENTRAL AND SOUTH AMERICA**

Dräger Medical AG & Co. KG  
Branch Office Dubai  
Dubai Healthcare City, P.O. Box 505108  
Dubai, United Arab Emirates  
Tel + 971 436 24 762  
Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION ASIA / PACIFIC**

Draeger Medical  
South East Asia Pte Ltd  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**Manufacturer:**

Dräger Medical AG & Co. KG  
23542 Lübeck, Germany  
The quality management system at  
Dräger Medical AG & Co. KG is  
certified according to ISO 13485,  
ISO 9001 and Annex II.3 of Directive  
93/42/EEC (Medical devices).