


Date / Datum 2025-07-30	Confirmation of Compatibility (CCR) EU 2017/745 Annex I, 14.1	
Document ID CCR108-032-2507-009-0		Page 1 of 2

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck
Germany

hereby confirms under its sole responsibility, that the medical device

VentStar Helix

may be used in accordance with its intended purpose together with

Product name	part number	Savina 300	IACS Workstation Critical Care (Evita V300/V500)	Babylog 8000	IACS WS Neonatal Care (VN500)	Evita V600	Evita V800	Babylog VN600	Babylog VN800	F&P MR850 Humidifier
VentStar Helix Dual Heated	MP02606	X	X							X
VentStar Helix Heated	MP02607	X	X							X
VentStar Helix Dual Heated N plus	MP02650	X	X	X	X	X	X	X	X	X
VentStar Helix Heated N plus	MP02608	X	X	X	X	X	X	X	X	X

- without losing or compromising the ability to perform as intended
- without the need for modification or adaption of any part of the combined devices
- without conflict/interference or adverse reaction

Combination of the listed products is done according to the current instructions for use including mounting instructions of the respective devices in therein indicated configurations only.

Any deviation of the indicated configurations not authorized by Drägerwerk AG & Co. KGaA will invalidate this confirmation.

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
Postal address:
23542 Lübeck, Germany
Tel +49 451 882-0
Fax +49 451 882-2080
info@draeger.com
www.draeger.com

Bank details:
Commerzbank AG, Lübeck
IBAN: DE95 2304 0022 0014 6795 00
Swift-Code: COBADE FF 230

Sparkasse zu Lübeck
IBAN: DE15 2305 0101 0001 0711 17
Swift-Code: NOLADE21SPL

Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7903 HL
General partner:
Drägerwerk Verwaltungs AG
Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7395 HL
UID-Nr. DE135082211

Chairwoman of the Supervisory Board for Drägerwerk AG & Co. KGaA: Maria Dietz

Chairman of the Supervisory Board for Drägerwerk Verwaltungs AG: Stefan Lauer

Executive Board:
Stefan Dräger (chairman)
Stefanie Hirsch
Rainer Klug
Gert-Hartwig Lescow
Dr. Reiner Piske
Anton Schrofner

Date / Datum 2025-07-30	Confirmation of Compatibility (CCR) EU 2017/745 Annex I, 14.1	Dräger
Document ID CCR108-032-2507-009-0		Page 2 of 2

Date / Datum
2025-07-30

Date / Datum
2025-07-30



Head of Regulatory Affairs
Business Unit Hospital Consumables &
Accessories
Medical Division

Head of Product Area Systems
Business Unit Hospital Consumables &
Accessories
Medical Division

Holger Wagner

Theresa Strampfer

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
Postal address:
23542 Lübeck, Germany
Tel +49 451 882-0
Fax +49 451 882-2080
info@draeger.com
www.draeger.com

Bank details:
Commerzbank AG, Lübeck
IBAN: DE95 2304 0022 0014 6795 00
Swift-Code: COBA DE FF 230

Sparkasse zu Lübeck
IBAN: DE15 2305 0101 0001 0711 17
Swift-Code: NOLADE21SPL

Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7903 HL
General partner:
Drägerwerk Verwaltungs AG
Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7395 HL
UID-Nr. DE135082211

Chairwoman of the Supervisory Board for
Drägerwerk AG & Co. KGaA: Maria Dietz

Chairman of the Supervisory Board for
Drägerwerk Verwaltungs AG: Stefan Lauer

Executive Board:
Stefan Dräger (chairman)
Stefanie Hirsch
Rainer Klug
Gert-Hartwig Lescow
Dr. Reiner Piske
Anton Schrofner