

1 Supplement to the instructions for use

1.1 Innovian Anesthesia software compatibility

Microsoft updates

Dräger recommends the regular installation of applicable software updates for Microsoft products on Innovian Anesthesia servers and workstations. Updates may be installed via Windows Update, Windows Server Update Services (WSUS), or System Center Configuration Manager (SCCM).

Innovian Anesthesia 8.0/8.0 SP1

Innovian Anesthesia Software	Compatibility
Server software	Microsoft Windows Server 2016 Microsoft Windows Server 2019
Database software	Microsoft SQL Server 2016 Microsoft SQL Server 2017 Microsoft SQL Server 2019
SQL Server Reporting Services reports	Microsoft SQL Server 2016 Microsoft SQL Server 2017 Reporting Services Microsoft SQL Server 2019 Reporting Services
HL7 Interface software	Rhapsody Integration Engine 5.5 (Rhapsody support ended 2023) Rhapsody Integration Engine 6.0-6.5
Workstation (Innovian Application Suite)	Microsoft Windows 10 Pro x64 Microsoft Windows 10 Enterprise x64 Microsoft Windows 11 Enterprise x64 (non-clinical)

Innovian IQ (Innovian Anesthesia 8.0 SP1 compatible)

Innovian IQ Software	Compatibility
Server software	Microsoft Windows Server 2016 Microsoft Windows Server 2019
Database software	Microsoft SQL Server 2016 Microsoft SQL Server 2019
SQL Server Reporting Services reports	Microsoft SQL Server 2016 Microsoft SQL Server 2019 Reporting Services

Innovian Anesthesia 7.0 / 7.0 SP1

Innovian Anesthesia Software	Compatibility
Server software	Microsoft Windows Server 2012 R2 (Microsoft support ended 2023) Microsoft Windows Server 2016 Microsoft Windows Server 2019
Database software	Microsoft SQL Server 2012 (Microsoft support ended 2022) Microsoft SQL Server 2014 (Microsoft support ended 2024) Microsoft SQL Server 2016 Microsoft SQL Server 2017 Microsoft SQL Server 2019
SQL Server Reporting Services reports	Microsoft SQL Server 2012 Microsoft SQL Server 2014 Microsoft SQL Server 2016 Microsoft SQL Server 2017 Reporting Services Microsoft SQL Server 2019 Reporting Services
HL7 Interface software	Rhapsody Integration Engine 5.5 (Rhapsody support ended 2023) Rhapsody Integration Engine 6.0-6.5
Workstation (Innovian Application Suite)	Microsoft Windows 10 Pro x64 Microsoft Windows 10 Enterprise x64

1.2

Trademarks

Trademarks owned by third-party manufacturers

Trademark	Trademark owner
HL7	Health Level Seven International
Microsoft	Microsoft Corp.
Microsoft Windows	
SQL Server	
Rhapsody	InterOperability Bidco, Inc., dba Rhapsody

This page has been left blank intentionally.



Manufacturer



Draeger Medical Systems, Inc.

3135 Quarry Road
Telford, PA 18969
U.S.A.



(215) 721-5400
(800) 4DRAGER

FAX

(215) 723-5935



<http://www.draeger.com>

European representative



Drägerwerk AG & Co. KGaA



Moislinger Alle 53 – 55
D-23542 Lübeck
Germany



+49 451 8 82-0

FAX

+49 451 8 82-2080



<http://www.draeger.com>



Distributed in the US by



Draeger, Inc.

3135 Quarry Road
Telford, PA 18969-1042
U.S.A.



(215) 721-5400
(800) 4DRAGER

FAX

(215) 723-5935



<http://www.draeger.com>



Distributed in Canada by



Draeger Medical Canada, Inc.

2425 Skymark Avenue, Unit 1
Mississauga, Ontario, L4W4Y6
Canada



(905) 212-6600
(866 343-2273)

FAX

(905) 763-1890



Canada.support@draeger.com

9511339 – en
© Draeger Medical Systems, Inc.
Edition: 2 – 2024-09
(Edition: 1 – 2020-05)
Dräger reserves the right to make changes to
the product without prior notice.

