

## Supplement to the instructions for use

---

This supplement applies to the medical devices listed in the table.

Medical device	Part number	Edition
Fabius GS	9038685	Edition 1 and higher
Fabius GS <i>premium</i>	9038923	Edition 1 and higher
Fabius Tiro	9038686	Edition 1 and higher
Fabius plus	9039054	Edition 1 and higher

- Keep this supplement with the instructions for use of the medical device.

### **WARNING**

**To properly use this medical device, read and comply with the instructions for use and this supplement.**

The supplement updates the information of the instructions for use in the following chapter:

## Technical Data


---

### Machine Data

#### Power supply

Rating non-configurable, with optional Dräger multiple socket outlet (see instructions for use for Multiple socket outlet for medical devices, 9038774)

100 to 240 VAC, 50/60 Hz, 2.3 Amax.


 Manufacturer:


**Dräger Medical GmbH**

 Moislinger Allee 53 – 55

D-23542 Lübeck

Germany

 +49 451 8 82-0

FAX  +49 451 8 82-20 80

 <http://www.draeger.com>

**9053746** – en

© Dräger Medical GmbH

Edition: 2 – 2015-01

(Edition: 1 – 2012-04)

Dräger reserves the right to make modifications to the medical device without prior notice.



As of 2015-08:  
Dräger Medical GmbH  
changes to  
**Drägerwerk AG & Co. KGaA**