



# LiteStar®

## Consumables and Accessories

---

Providing high patient comfort during anesthesia and improving handling for caregivers: The teardrop shape and recessed housing of our anesthesia mask LiteStar support an anatomical fit with effective sealing and provide easier and comfortable grip for your caregiver routines.

# Benefits

## LiteStar

The Dräger mask LiteStar is intended for single patient use in anesthesia. The new design with the recessed grip supports a perfect control for the physician. In addition, the transparent mask body allows for continuous observation of the patient's condition.

## Color-Coding

Every mask size is visualized by a different label color. This labeling supports a quick and easy identification during picking which saves time during the daily routine.

## Additional mask accessories

Easy application is supported through additional mask accessories such as a hook ring and headgears for fixation.

Two different headgears in small and large designs offer patients and physicians a soft and elastic support and holding position.

# Details



D-37691-2021

Recessed housing side-face




D-37702-2021

Recessed housing front-face




D-44293-2021

Hook ring for LiteStar mask



MT-9654-2006

Headgear, soft rubber



MT-9658-2006

Headgear, soft fabric

# Ordering Information

<b>LiteStar anesthesia face mask, disposable, without hook ring</b>	<b>Packaging size</b>	<b>Part number</b>
size 1, small infant	50 pcs.	MP11001
size 2, large infant	50 pcs.	MP11002
size 3, child	50 pcs.	MP11003
size 4, small adult	50 pcs.	MP11004
size 5, medium adult	50 pcs.	MP11005
size 6, large adult	50 pcs.	MP11006
size 7, extra large adult	50 pcs.	MP11007
<b>Hook ring for LiteStar masks</b>		
with 22 mm ID connector (size 3-7)	30 pcs.	MP11008
<b>Headgear, Soft Fabric, disposable</b>		
small	1 pcs.	MP01559
large	1 pcs.	MP01560
<b>Headgear, Soft Rubber, disposable</b>		
small	1 pcs.	MP01561
large	1 pcs.	MP01562

Not all products, features, or services are for sale in all countries.  
Mentioned Trademarks are only registered in certain countries and not necessarily in the country  
in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

**Corporate Headquarters**  
**Drägerwerk AG & Co. KGaA**  
Moislinger Allee 53-55  
23558 Lübeck  
🌐 [www.draeger.com](http://www.draeger.com)

**Region Europe**  
**Drägerwerk AG & Co. KGaA**  
Moislinger Allee 53-55  
23558 Lübeck, Germany  
☎ +49 451 882 0  
☎ +49 451 882 2080  
✉ [info@draeger.com](mailto:info@draeger.com)

**Region Asia Pacific**  
**Draeger Singapore Pte. Ltd.**  
61 Science Park Road  
The Galen #04-01  
Singapore 117525  
☎ +65 6872 9288  
☎ +65 6259 0398  
✉ [asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**Region Middle East, Africa**  
**Dräger Safety AG & Co. KGaA**  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
☎ +971 4 4294 600  
☎ +971 4 4294 699  
✉ [contactuae@draeger.com](mailto:contactuae@draeger.com)

**Manufacturer**  
**Drägerwerk AG & Co. KGaA**  
Moislinger Allee 53-55  
23558 Lübeck

**Region Central and South America**  
**Dräger Indústria e Comércio Ltda.**  
Al. Pucurui - 51 - Tamboré  
06460-100 - Barueri - São Paulo  
☎ +55 (11) 4689-4900  
✉ [relacionamento@draeger.com](mailto:relacionamento@draeger.com)



Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)