



Catheter Mounts Breathing Circuits

Catheter mount tube extensions and airway connectors make it easier for you to safely move adult and pediatric patients. They are very easy to use and help keep your hands free, giving you more room to provide care. ErgoStar® catheter mounts and catheter connectors are just one part of Dräger's complete line of breathing system equipment.

Highlights

ErgoStar® catheter mounts

Dräger smooth bore catheter mounts are flexible and lightweight, and provide low resistance to gas flow, while the external coil reduces kinking and blockage. Clinicians can enjoy the ultimate in flexibility and ease-of-use with the double swivel connector that minimizes traction and torque on the tube. The dual cap design allows you to maintain PEEP during suction or bronchoscopy and the device-sided 22F/15M connector serves as a 2-in-1 solution with less dead space. Each flexible catheter mount and airway connector is also equipped with a red safety cap so that the system is completely closed and protected until immediately prior to being used on the patient. There is also a Luer-Lock sample port at a practical 45° angle, allowing for comfortable gas sampling as close to the patient as possible. Enjoy improved workflow thanks to the single-use design and the standard connectors, allowing you to securely and safely connect to other components. All Dräger catheter mounts are also PVC- and latex-free, eliminating the risk of softeners such as DEHP.

Connective Products



Dräger Babyroo® TN300

- Configurable open care warmer for the labour and delivery unit and NICU
- Provides stable temperature with two combined heat sources
- Includes optional Respiratory Support Module with AutoBreath
- Meets resuscitation guidelines and supports intuitive and ergonomic handling



Dräger Evita V800

- Seamless transition between O₂ therapy, NIV, and invasive ventilation modes.
- Individualised lung protective ventilation for all patient groups with advanced diagnostic tools.
- Support of a synchronised way for quick and efficient weaning.
- Streamlined user interface allows intuitive, efficient operation and guidance.
- Natively implemented SDC-standard allows secure alarm distribution



Dräger Evita® V600

- Early Mobilization and Patient Transport
- Lung Protective Ventilation
- Automatic weaning with SmartCare®/PS
- Future-proof open connectivity



Dräger Savina® 300 Select

- Turbine-driven ventilation
- Intuitive and simple operation
- Sophisticated ventilation modes

Not all products, features, or services are for sale in all countries. Trademarks mentioned herein are the property of its respective owner. Trademarks may be owned by Drägerwerk AG & Co. KGaA (Dräger) or its affiliates in certain countries and not necessarily in the country in which this material is released. Visit www.draeger.com/trademarks for the current status of Dräger's trademarks.

Dräger South Africa (Pty) Ltd.

Johannesburg: 2 Ruacana Street, Waterfall Commercial District, (cnr Beatty Street), Buccleuch, 2090
Tel: +27(0)11 059 4200

Dräger South Africa (Pty) Ltd.

Cape Town: Unit 2, 7 on London, 7 London Circle, Brackengate Business Park Brackenfell, 7561
Tel: +27(0) 21 819 3000

Dräger South Africa (Pty) Ltd.

Durban: 17 Thynk Industrial Park Cnr Brickworks and Queen Nandi Drive, Briardene, Durban North, 7561
Tel: +27(0) 31 812 7900

Dräger Kenya Ltd

Nairobi: 8th Floor Office 2, Sanlam Tower, Westlands
Tel: +25(0) 47 462 40336

Corporate Headquarters

Drägerwerk AG & Co. KGaA

Moislinger Allee 53-55
23558 Lübeck
Germany

