

## SediStar® - Anesthesia face masks Consumables and Accessories

We help you by providing high patient comfort during anesthesia and improved handling for caregivers. The teardrop shape of our SediStar supports an anatomical fit with an effective seal, making your procedures easier and more comfortable. Beyond that the masks are compatible with all our anesthesia devices, they are free of harmful chemicals to help protect you and your patients.



## Benefits

---

### Supporting infection prevention to keep patients safe

During anesthetic procedures, every breath your patient takes is critical. We share a common goal: to protect your patients from infection during anesthesia. To help you achieve this goal and break the chain of infection, our SediStar is a disposable, single patient use anesthesia mask. The teardrop shape makes it easy to put our SediStar on without touching the sealing lip for quick and secure application, protecting both staff and patients from healthcare-associated infections.

---

### Supports patient safety by providing a clear view of the patient's face

Full visibility of your patient's face during anesthesia-induced procedures is critical for safe and continuous observation. SediStar masks have a transparent mask body that allows you to observe your patient's condition.

---

### Combining comfort and ease of use

Safe anesthesia delivery means having the right equipment that not only keeps your patient safe, but also makes the clinician's job easier. Our SediStar anesthesia masks are ergonomically designed to give you optimal handling, effective seal and comfortable fit for different facial anatomies. Additional mask accessories, such as different headgears, can help to make handling easy, fast and comfortable for you and your patients.

---

### Easy to find the best fit

The correct fit of an anesthesia mask is essential for patient comfort and effective anesthesia procedures. Our SediStar masks offer an optimised fit for different facial anatomies through their ergonomic design and sealing lip. To make it easy to identify the right mask, each mask size is identified by a different colour label. This labelling supports quick and easy identification during picking, saving time in your daily routine and additionally, help your patient's comfort during anaesthesia.

---

### Helping to shape a more sustainable future

As a manufacturer of medical technology for over 130 years, we are committed to provide excellence in the quality of our products. And we are also committed and believe in taking responsibility in creating sustainable products. We designed our SediStar anaesthesia masks to combine patient safety and sustainability. That's why our SediStar masks are free of natural latex and do not contain PVC or phthalates, e.g. DEHP (which is possibly carcinogenic to humans<sup>1</sup>). Additionally they do not contain BPA, which is potentially toxic.<sup>2</sup> This is our commitment, this is the responsibility we have for a healthier environment.

<sup>1</sup> World Health Organization. *Guidelines for drinking-water quality*. World Health Organization. 2017. 424-425. ISBN 978-92-4-154995-0.

<sup>2</sup> Risk Management for Bisphenol A (BPA) | US EPA | Bisphenol A | EFSA (europa.eu)

## Details



D-4741-2019

Housing side-face




D-4737-2019

Housing front-face



D-25951-2020

Hook ring for SediStar mask



MT-9654-2006

Headgear, soft rubber



MT-9658-2006

Headgear, soft fabric

## Ordering Information

	Packaging size	Part number
SediStar anesthesia face mask, disposable, without hook ring - Available in following sizes:		
size 1, small infant	30 pcs.	MP01701
size 2, large infant	30 pcs.	MP01702
size 3, pediatric	30 pcs.	MP01703
size 4, small adult	30 pcs.	MP01704
size 5, medium adult	30 pcs.	MP01705
size 6, large adult	30 pcs.	MP01706
size 7, extra large adult	30 pcs.	MP01707
Hook ring for SediStar masks with 22 mm ID connector (size 3-7)	30 pcs.	MP01708
Headgear, Soft Fabric, disposable		
small	1 pcs.	MP01559
large	1 pcs.	MP01560
Headgear, Soft Rubber, disposable		
small	1 pcs.	MP01561
large	1 pcs.	MP01562

Not all products, features, or services are for sale in all countries.  
Mentioned Trademarks are only registered in certain countries and not necessarily in the country  
in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

### CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
[www.draeger.com](http://www.draeger.com)

### Manufacturer:

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23542 Lübeck, Germany

### REGION EUROPE

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.  
61 Science Park Road  
The Galen #04-01  
Singapore 117525  
Tel: +65 6872 9288  
Fax: +65 6259 0398  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.  
Al. Pucurui - 51 - Tamboré  
06460-100 - Barueri - São Paulo  
Tel. +55 (11) 4689-4900  
[relacionamento@draeger.com](mailto:relacionamento@draeger.com)

Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

