



Silicone Masks Consumables and Accessories

Making hospital stays as comfortable as possible for patients is our common goal. To effectively enhance patient comfort our silicone masks are made of medical grade silicone that has superior ability to conform to the patient's anatomy.

Benefits

Silicone Masks

The medical grade silicone is resistant to high or low temperature and durable. Our silicone masks can be autoclaved and thus, they can be reused on multiple patients.

There are different mask designs available. The smaller sizes (sizes 00-2) have a rather round shape in order to provide a comfortable fit and an effective seal for infants with more round shaped face contours. The bigger sizes (sizes 3-5) have an anatomically formed teardrop shape to ensure a convenient fit and seal to more adult patients' facial characteristics. The silicone masks have a flexible, soft, and transparent mask body.

The mask body design allows for a comfortable rest of the thumb. The form and material create an easy grip for the physician. The transparent mask body allows for continuous visual observation of the patient's condition.


Easy application is supported through additional mask accessories such as hook rings and headgears for fixation. Two different headgears in small and large designs offer patients and physicians a soft and elastic support and holding position.

Details



D-7361-2011

Hook ring for silicone mask



MT-9654-2006

Headgear, soft rubber



MT-9658-2006

Headgear, soft fabric

Ordering Information

Anesthesia silicone mask, reusable, without hook ring	Packaging size	Part number
size 00, premature baby	1 pcs.	MP02900
size 0, small neonate	1 pcs.	MP01910
size 1, large neonate	1 pcs.	MP01911
size 2, small child	1 pcs.	MP01912
size 3, large child	1 pcs.	MP01913
size 4, small adult	1 pcs.	MP01914
size 5, large adult	1 pcs.	MP01915

Hook ring		
reusable, small, for mask sizes 00-1 with 15 mm OD connector	10 pcs.	MP01916
reusable, large, for mask sizes 2-5 with 22 mm ID connector	10 pcs	MP01917

Headgear, Soft Fabric, disposable		
small	1 pcs.	MP01559
large	1 pcs.	MP01560

Headgear, Soft Rubber, disposable		
small	1 pcs.	MP01561
large	1 pcs.	MP01562

Not all products, features, or services are for sale in all countries. Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

Corporate Headquarters

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
■ www.draeger.com

Manufacturer

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

Region Europe

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
☎ +49 451 882 0
☎ +49 451 882 2080
✉ info@draeger.com

Region Middle East, Africa

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
☎ +971 4 4294 600
☎ +971 4 4294 699
✉ contactuae@draeger.com

Region Asia Pacific

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
☎ +65 6872 9288
☎ +65 6259 0398
✉ asia.pacific@draeger.com

Region Central and South America

Dräger Indústria e Comércio Ltda.
Al. Pucuruí 51 - Tamboré
06460-100 - Barueri - São Paulo
☎ +55 (11) 4689-4900
✉ relacionamento@draeger.com



Locate your Regional Sales
Representative at:
www.draeger.com/contact