






















### Accessories & Consumables - **FABIUS Family**

Family	Part No	Main Description	Units per pack
	<b>MX00004</b>	Drägersorb® 800+ CLIC, disposable absorber, 6 pcs	6
	<b>MX50100</b>	Drägersorb® CLIC Free, disposable absorber, 6 pcs	6
	<b>MX00012</b>	Draegersorb CLIC Adapter	1
	<b>MX00013</b>	Draegersorb CLIC Intermediate Plate (Fabius Adapter plate not needed for Cosy 2.6, only earlier Cosy's that have the big refillable canister.	1
	<b>MX00001</b>	Draegersorb 800+ Soda Lime (5 Litre)	1
	<b>MX50050</b>	Draegersorb Free (5-Litre)	1
	<b>MKO1900</b>	SpiroLife sensor -Autoclavable Reusable	1
	<b>8403735</b>	Spirolog flow sensor, 5 pcs	5
	<b>6872130</b>	Watertrap for SCIO (Set of 12)	12











### Accessories & Consumables - **FABIUS Family**

Family	Part No	Main Description	Units per pack
	8290286	Sample Line Set Single Use	10
	6850645	Oxygen sensor (O2 capsule)	1
<a href="#">To see full range of Breathing Filters - please see attached</a>			
	MP01770	Filter CareStar 30 (PK/50) <i>Bacterial/viral filter, compressible volume: 30 ml; recommended tidal volume: 100 - 1500 ml</i>	50
	MP01790	Filter SafeStar 55 'Mechanical' (PK/50) <i>HEPA breathing filter, compressible volume: 55 ml; recommended tidal volume: 200 - 1500 ml</i>	50
	MP01805	Filter / HME TwinStar 55 (PK/50) <i>Bacterial/viral filter and HME, compressible volume: 55 ml; recommended tidal volume: 200 - 1500 ml</i>	50
	MP00321	Anesthesia Set Basic, Latex Free, 2 corrugated EVA-hoses <i>1.5 m (4.9 ft), Y-piece with Luer Lock, breathing bag with connector, 2L latex free bag.</i>	25
	MP00325	Vent Set Flex, Latex-Free <i>Standard Flex breathing circuit: 2 extendable PP-hoses up to 1.8 m (5.9 ft), Y-piece with Luer Lock.</i>	25
	MP00320	Vent Set Basic, Latex-Free <i>Standard breathing circuit: 2 corrugated EVA-hoses 1.5 m (4.9 ft), Y-piece with Luer Lock.</i>	25
	MP00317	Vent Set COAX, Latex-Free. <i>Coaxial breathing circuit (hose-in- hose) 1.5 m (4.9 ft), space saving and effective. Inspiration gas is warmed up by the expiration gas .</i>	20
		<b>ComfortStar® Standard</b> sizes: 0-6 Size 0: Neonate (# MP01510) Size 1: New Born (# MP01511) Size 2: Infant (# MP01512) Size 3: Toddler(# MP01513) Size 4: Small Adult (# MP01514) Size 5: Medium Adult (# MP01515) Size 6: Large Adult (# MP01516)	20


## Accessories & Consumables - **FABIUS Family**

Family	Part No	Main Description	Units per pack
		<b>ComfortStar® Mint</b> Scented sizes: 3-6 For adult application. Size 3: Toddler (# MP01543) Size 4: Small Adult (# MP01544) Size 5: Medium Adult (# MP01545) Size 6: Large Adult (# MP01546)	20
		<b>ComfortStar® Strawberry</b> Scented sizes: 1-4 For paediatric application. Size 1: New Born (# MP01521) Size 2: Infant (# MP01522) Size 3: Toddler (# MP01523) Size 4: Small Adult (# MP01524)	20
		<b>ComfortStar® Bubblegum</b> Scented sizes: 1-4 For paediatric application. Size 1: New Born (# MP01531) Size 2: Infant (# MP01532) Size 3: Toddler (# MP01533) Size 4: Small Adult (# MP01534)	20
		<b>LiteStar® Mask</b> sizes: 1-6 Size 1: New Born (# MP01501) Size 2: Infant (# MP01502) Size 3: Toddler (# MP01503) Size 4: Small Adult (# MP01504) Size 5: Medium Adult (# MP01505) Size 6: Large Adult (# MP01506)	30
	<b>MP01840</b>	ErgoStar® CM 40 Catheter mount, double-swivel, double cap connector, 10 cm.	50
	<b>MP01845</b>	ErgoStar® CM 45 Catheter mount, double-swivel, double cap connector, 15 cm.	50
	<b>MP01850</b>	ErgoStar® CM 50 Extendible Catheter mount, double-swivel, double cap connector, 10-16cm.	50
	<b>MP01855</b>	ErgoStar® CM 55 Catheter mount, 15cm	50
	<b>MP01860</b>	ErgoStar® CM 60, Catheter mount, double-swivel, double cap connector, 10 cm. (16ml low deadspace for pediatric application)	50

### Accessories & Consumables - **FABIUS Family**

Family	Part No	Main Description	Units per pack
	<b>MP01890</b>	ErgoStar® AC 90, double-swivel, double cap connector	100
	<b>MP01895</b>	ErgoStar® AC 95, 90°-connector with luer port for gas sampling	100
	<b>2166054</b> <b>2166127</b> <b>2166062</b> <b>2166070</b>	Silicone breathing bag 0.5 l Silicone breathing bag 1.5 l Silicone breathing bag 2.3 l Silicone breathing bag 3.0 l	1
	<b>M30146</b>	Reusable breathing hose set for standard anaesthesia ventilation.	1
	<b>8404985</b>	Hose water trap.	1
	<b>2165856</b>	Breathing hose, paediatric, 40 cm, Reusable	1
	<b>2165821</b>	Breathing hose, paediatric, 60 cm, Reusable	1
	<b>2165651</b>	Breathing hose, paediatric, 110 cm, Reusable	1
	<b>2165961</b>	Breathing hose, paediatric, 150 cm, Reusable	1
	<b>M25649</b>	Elbow connector	1

**Accessories & Consumables - FABIUS Family**

Family	Part No	Main Description	Units per pack
	8403075	Y-piece, child, 90°	1
	8401645	Cap, 7 mm (0.3 in) (To seal the opening of the Y-piece)	1