

ComfortStar® - Anesthesia face masks Consumables and Accessories

During anaesthesia-induced surgical procedures, every breath your patient takes is critical. With its unique teardrop or round-shape design, our ComfortStar® anaesthesia face masks provide a superior anatomical fit and seal performance. Comfortable to wear and easy to use, it is also compatible with all Dräger anaesthesia devices. ComfortStar® supports every patient age to breathe easier.



Benefits

Comfort and flexibility - all in one single use: ComfortStar anaesthesia face masks

The safe delivery of anaesthetics involves having the right equipment that not only keeps your patient secure, but also facilitates the clinician's work. To give you optimal flexibility of use, our ComfortStar masks are designed with an inflatable cushion that gives you a comfortable, steady grip. When choosing ComfortStar mask sizes 4-6, you have the additional feature of built-in ridges in the body of the mask, which provide you with a safe, non-slip grip.

Moreover, the complete visibility of your patient's face during anaesthesia-induced procedures is critical for safe and continuous observation. ComfortStar masks have transparent mask body that allows you to observe every detail of your patient's condition. Similarly, for quick and easy size identification, all masks are pre-equipped with 'Traffic Light' coloured hook rings. Smaller sizes (0-1) are red in colour, medium sizes (2-3) are yellow, and large sizes (4-6) are green. You can also include small and large headgears for a simple and secure fix.

Scented options for more comfort: ComfortStar anaesthesia face masks

A potential side effect for patients who have undergone surgical procedure under general anaesthesia is postoperative nausea and vomiting (PONV). This often delays patient recovery. Research shows that added aromas can help minimise these unpleasant side effects ¹⁾.

To minimise postoperative discomfort, ComfortStar scented anaesthesia face masks have the same excellent benefits and features of the ComfortStar standard masks. They are available in either strawberry or bubble gum scents when choosing sizes 1-4, and mint scent in sizes 3-6. For added optical stimulation, we have also designed our ComfortStar masks with an appropriate colour scheme which is intended to reduce pre-operative anxiety in patients.

¹⁾ Anderson, L.A., Gross, J.B.: Aromatherapy with peppermint, isopropyl alcohol, or placebo is equally effective in relieving postoperative nausea, 2004, in: Journal of PeriAnesthesia Nursing, Vol. 39, No 1, 29-55

Ordering Information

Description	Package of	Order No.
ComfortStar anesthesia face mask, disposable, standard, with hook ring		
size 0, neonate	20 pcs.	MP01510
size 1, new born	20 pcs.	MP01511
size 2, infant	20 pcs.	MP01512
size 3, toddler	20 pcs.	MP01513
size 4, small adult	20 pcs.	MP01514
size 5, medium adult	20 pcs.	MP01515
size 6, large adult	20 pcs.	MP01516
ComfortStar anesthesia face mask, disposable, strawberry scented, with hook ring		
size 1, new born	20 pcs.	MP01521
size 2, infant	20 pcs.	MP01522
size 3, toddler	20 pcs.	MP01523
size 4, small adult	20 pcs.	MP01524
ComfortStar anesthesia face mask, disposable, bubble gum scented, with hook ring		
size 1, new born	20 pcs.	MP01531
size 2, infant	20 pcs.	MP01532
size 3, toddler	20 pcs.	MP01533
size 4, small adult	20 pcs.	MP01534
ComfortStar anesthesia face mask, disposable, mint scented, with hook ring		
size 3, toddler	20 pcs.	MP01543
size 4, small adult	20 pcs.	MP01544
size 5, medium adult	20 pcs.	MP01545
size 6, large adult	20 pcs.	MP01546
Headgear, Soft Fabric, disposable		
small	1 pc.	MP01559
large	1 pc.	MP01560
Headgear, Soft Rubber, disposable		
small	1 pc.	MP01561
large	1 pc.	MP01562
Two-sided pump ball for air management	1 pc.	MP01590

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel: +65 6872 9288
Fax: +65 6259 0398
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06460-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

CANADA

Draeger Medical Canada Inc.
2425 Skymark Avenue, Unit 1
Mississauga, Ontario, L4W 4Y6
Tel +1 905 212 6600
Toll-free +1 866 343 2273
Fax +1 905 212 6601
Canada.support@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

