

## Dräger SpO<sub>2</sub> Sensors Consumables and Accessories

Reusable and Disposable SpO<sub>2</sub> Sensors.



## Benefits

---

### Dräger SpO<sub>2</sub> Sensors

- Dräger SpO<sub>2</sub> Sensors are compatible with Dräger Infinity® Patient Monitors.
- Dräger SpO<sub>2</sub> Sensors are validated for all Dräger Infinity Patient Monitors.
- Dräger SpO<sub>2</sub> Sensors are compatible with the Nellcor® Oxisensor II.

## Technical Data

LED specification red 660 nm typical 1.8 mW typical

Infrared 905 nm typical 2.0 mW typical

SpO<sub>2</sub> Accuracy adults or pediatrics: 70-100% ± 2 digits

SpO<sub>2</sub> Accuracy in neonates: 70-100% ± 3 digits

The accuracy of the pulse rate depends on the monitor.

## Ordering Information

### Dräger Sensors

	Order number
SpO <sub>2</sub> Sensor, Adults, Reusable	MS 13 235
SpO <sub>2</sub> Sensor, Adults, Vinyl, Disposable	MS 16 449
SpO <sub>2</sub> Sensor, Adults, Foam, Disposable	MS 16 445
SpO <sub>2</sub> Sensor, Pediatrics, Vinyl, Disposable	MS 16 448
SpO <sub>2</sub> Sensor, Pediatrics, Foam, Disposable	MS 16 444
SpO <sub>2</sub> Sensor, Infants, Fabric, Disposable	MS 16 447
SpO <sub>2</sub> Sensor, Neonates, Fabric, Disposable	MS 16 446

### Adapter cable

To connect Dräger/ Nellcor SpO <sub>2</sub> Sensors to MultiMed 5/6, MultiMed 12 Pod and NeoMed Pod	
Extension cable SpO <sub>2</sub> cable, 1 m	33 68 433
Extension cable SpO <sub>2</sub> cable, 2 m	33 75 834
To connect Dräger/ Nellcor SpO <sub>2</sub> Sensors to MultiMed plus, MultiMed plus OR Pods and M300	
SpO <sub>2</sub> cable, 1.2 m	MS 18 683
SpO <sub>2</sub> cable, 3.0 m	MS 17 330

## Notes

### CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
[www.draeger.com](http://www.draeger.com)

### Manufacturer:

Dräger Medical GmbH  
Moislinger Allee 53–55  
23558 Lübeck, Germany

### As of August 2015

Dräger Medical GmbH changes  
to Drägerwerk AG & Co. KGaA

### REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION EUROPE SOUTH

Dräger Médical S.A.S.  
Parc de Haute Technologie  
d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[d1mfr-contact@draeger.com](mailto:d1mfr-contact@draeger.com)

### REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial  
S. de R.L.  
Complejo Business Park,  
V tower, 10th floor  
Panama City  
Tel +507 377 9100  
Fax +507 377 9130  
[contactcsa@draeger.com](mailto:contactcsa@draeger.com)

### REGION ASIA / PACIFIC

Draeger Medical  
South East Asia Pte Ltd.  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### REGION NORTH AMERICA

Draeger Medical, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 215 721 5400  
Toll-free+1 800 437 2437  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

Locate your Regional  
Sales Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

