

MiniMe® 2 NIV pediatric nasal mask Consumables and Accessories

Dräger offers the SleepNet MiniMe® 2 mask available in two sizes for a greater range of comfort for children between ages 2 and 12. Enhanced Custom Fit Technology™ enables the mask to be easily molded to a child's face, maximizing comfort. The MiniMe® 2 also features Advanced Cushion Technology (A.C.T.), an ergonomic cushion technology that delivers an effective seal yet minimizes pressure on small faces.

AIR°gel

Maximum comfort, minimum weight. Patented Air°gel makes the MiniMe 2 one of the most gentle fitting masks on the market.



Touchless Spacebar™

Unique design eliminates any contact with the patient's forehead.

Size Indicator

Both Vented and Non-Vented come in two sizes.

Custom Fit Technology

The shapeable shell mask can be easily molded to fit individual needs.
Flex, Fit, Forget.

Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility

Vented and Non-Vented Versions Available

For CPAP and Bi-level therapy in children



Benefits

Custom Fit Technology™

A unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit and comfort.

Advanced Cushion Technology™

The ergonomic AIR gel® cushion design provides the utmost in gentle cushioning to protect the child's delicate skin from pressure points and irritations which can lead to infection while maintaining an effective seal.

Cleaning

The MiniMe® mask should be cleaned daily in warm water with a mild detergent. Rinse well and allow to air dry away from direct sunlight. Single-patient use only.

- Available in two sizes (Small and Large) with 4 point simple and breathable headgear.
 - Products do not contain natural rubber latex, PVC, or DEHP.
- With the addition of the MiniMe® 2 mask, Dräger offers a complete line of single-patient use NIV masks from neonate to adult supporting the continuum of care.

Technical Data

MINIMI® 2 MASKS WITH HEADGEAR SPECIFICATIONS

Product Number	Style	Headgear Size	Leak Port	Mask Size	Mask Weight	Dead Space	Mask Dimensions
MP01553	Vented	Small	√	Small	0.069 lbs	22 ml	2.32" D x 2.00" W x 1.71" H
MP01554	Vented	Large	√	Large	0.089 lbs	40 ml	2.59" D x 2.69" W x 2.14" H
MP01555	Non-Vented	Small	X	Small	0.069 lbs	22 ml	2.32" D x 2.00" W x 1.71" H
MP01556	Non-Vented	Large	X	Large	0.089 lbs	40 ml	2.59" D x 2.69" W x 2.14" H

Info about Sleepnet

- 100% designed and manufactured in the USA
- Made with environmentally conscious production practices



Does not contain natural rubber latex

Scan code to view MiniMe 2 Quick Fit Guide



Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country
in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Sleepnet Corporation
5 Merrill Industrial Drive,
Hampton, New Hampshire 03842
USA

USA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969-1042
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

CANADA

Draeger Medical Canada Inc.
2425 Skymark Ave
Mississauga, Ontario, L4W 4Y6
Tel +1 905 212-6600
Toll-free +1 866 343 2273
Fax +1 905 212-6601
Canada.support@draeger.com