

MultiMed[®] plus and MultiMed[®] plus OR Consumables and Accessories

The Infinity[®] MultiMed[®] plus and plus OR cables combine ECG, pulse oximetry and temperature monitoring into a single lightweight package – reducing cable clutter at the point of care.



Benefits

MultiMed plus and MultiMed plus OR

Both cables support Infinity patient monitors and are compatible with bedside and patient-worn monitoring accessories. MultiMed plus and MultiMed plus OR cables are durable, easy to use, and help facilitate more efficient patient transport. For the Intensive Care environment, the MultiMed plus provides 3, 5 and 6-lead capability and can be used with our convenient MonoLeads® or traditional ECG-lead set.

For the OR, the MultiMed plus OR provides 3, 5 and 6-lead capability, plus an integrated electro surgical unit (ESU) filter that helps to eliminate High Frequency surgery artifacts for improved reliability and reduced weight.

- Lets you standardize on one ECG cable system for patient monitors, Infinity TruST® Telemetry* and Infinity M300.
- Supports 3, 5 and 6-lead ECG configurations

* MonoLead is not compatible with InfinityTruST Telemetry.

Standardize on one ECG-cable connection

With the new MultiMed plus and MultiMed plus OR cable, you can use the same ECG leads for OR, intensive care, and telemetry. The single-plug connection lets you quickly and easily swap out ECG-leads between MultiMed plus/MultiMed plus OR, TruST Telemetry and the Infinity M300, saving valuable time during patient transfer. In addition, the durable single-pin design helps ensure long life and reliable, positive connection.

Integrated ESU filter within MultiMed plus OR

Eliminate the need for separate ESU filters when connecting the ECG cables to remove HF surgery artifacts. Now, Dräger offers a single integrated solution that reduces weight and increases efficiency while simplifying processes.

Technical Data

Size (H x W x D)	20.5 x 82 x 120.9 mm
Weight	498.9 g (17.6 oz)
Connections:	Infinity Gamma, Gamma XL, Gamma X XL, Delta, Delta XL, Kappa, Kappa XLT, Vista or Vista XL
2.5 m connection cable (MultiMed to Monitor)	

Environmental Requirements

Humidity (operating):	20% to 90%, non-condensing
Storage:	10% to 95% with packaging
Altitude range (operating):	525 to 795 mmHg
Storage:	375 to 795 mmHg
Temperature range (operating):	10 °C to 40 °C
Storage:	-20 °C to 50 °C
Water resistance IPX1, drip proof	

Standards

IEC 60601-1, Isolated patient connection, Type Cardiac Floating, CF

ECG Monitoring

Patient category	adult and pediatric
Available leads	I, II, III, aVR, aVF, aVL, V1 to V6

Lead labeling on cable IEC1 and IEC 2 color coding. MultiMed® plus and plus OR are intermediate cables. For other relevant specifications see applicable monitor datasheet. For a complete list of Dräger approved accessories, please contact your local Dräger representative.

Ordering Information

MultiMed® plus Cable	MS20093
MultiMed® plus OR Cable	MS20094

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA

Locate your Regional
Sales Representative at:
www.draeger.com/contact



REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).