

## Oxylog® Trolley

With the easy-to use functional Oxylog® trolley system intensive care patients can quickly and smoothly be transported throughout the hospital without interrupting therapy.



D-34704-2011

### THE COMPREHENSIVE TRANSPORT SOLUTION FOR INTENSIVE CARE PATIENTS

Equipped with all the accessories needed to transport intensive care patients, the space-saving trolley can easily be attached to the patient bed where it is safely installed and keeps all the important equipment within reach:

- a transport ventilator, such as Oxylog® 3000 plus or Oxylog® VE300
- a transport monitor, such as Infinity® M540
- infusion bags
- a suction system
- a connection to the bed
- a pressure reducer
- a cylinder holder for two O<sub>2</sub> cylinders
- a power socket strip

It is no longer necessary to attach equipment parts to the bed for transport: The whole transport equipment can be installed on the trolley and is always ready for use.

### TRANSPORT WITH INTELLIGENT ACCESSORIES

Combined with the Dräger systems for therapy and monitoring and other accessories the Oxylog trolley is the comprehensive solution for interhospital patient transport. The transport ventilator Dräger Oxylog 3000 plus is able to ensure continuous ventilation at a high quality during the whole transport. If combined with, for example, Infinity® monitoring by Dräger, this allows continuous patient monitoring.

Dräger offers the Oxylog trolley in five different versions. You can add a device holder or a carrying system holder to your trolley.

**Version A:** Trolley with device holder for the Oxylog 3000 plus

**Version B:** Trolley with carrying system holder for the Oxylog 3000 plus

**Version C:** Trolley with device holder for the Oxylog VE300

**Version D:** Trolley with carrying system holder for the Oxylog VE300

**Version E:** Without holder for ventilator

D-76118-2013



Oxylog® Trolley with device holder

D-76125-2013



Oxylog® Trolley with carrying system holder

## ORDERING INFORMATION

|   | Order no.      |
|---|----------------|
| <b>Oxylog® trolley</b>                              | <b>5790150</b> |
| Connection to the bed                               | 8419800        |
| Gas cylinder holder                                 | G93110         |
| Universal holder with standard rail                 | G93140         |
| Swivel arm for patient monitor                      | G25709         |
| Swivel arm for accessories                          | G93111         |
| Cylinder pack for bronchial suction                 | 2M85056        |
| Rail support with claw                              | M26842         |
| Oxylite pressure reducer                            | MP00100        |
| Connection tube mount                               | 2M85446        |
| Power pack holder for standard rail                 | MC00155        |
| Transport holder and clamp for patient monitor M540 | MS28144        |
| Swivel arm for transport holder                     | MS20427        |
| Power socket strip (configuration)                  | 2609000        |
| Power socket strip attachment kit                   | 5790936        |



Not all products, features, or services are for sale in all countries.

Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

### CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

### Manufacturer:

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23542 Lübeck, Germany

### REGION DACH

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION EUROPE

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
Tel +971 4 4294 600  
Fax +971 4 4294 699  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### REGION ASIA PACIFIC

Dräger Singapore Pte. Ltd.  
61 Science Park Road  
The Galen #04-01  
Singapore 117525  
Tel +65 6872 9288  
Fax +65 6259 0398  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### REGION CENTRAL AND SOUTH AMERICA

Dräger Panamá Comercial  
S. de R.L.  
59 East Street, Nuevo Paitilla,  
House 31, San Francisco Town  
Panama City, Panama  
Tel +507 377-9100 / Fax -9130  
[servicioalcliente@draeger.com](mailto:servicioalcliente@draeger.com)

Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

