


Quickstop hinged arms Consumables and Accessories

Keeping everything close at hand. To leave your hands free to care for the patient, all the joints of the Quickstop hinged arm can be fixed into place and released with just one quick hand movement.



MT-0804-2008

Benefits

Quickstop hinged arms

The hinged arm is totally flexible, allowing it to be positioned precisely where it is needed. All cables, lines and hose systems can be arranged in such a way as to more or less completely prevent accidental disconnection, making the Quickstop hinged arm an essential accessory in anaesthesia, intensive care, surgery or indeed just about any other area of the hospital. The arm is provided with a universal rail bracket which can be used to attach it to all kinds of different objects (operating tables, beds, wall-mounted rails and devices). Thanks to the high quality of its materials, you will appreciate the benefits of this durable and maintenance-free Dräger product for years to come. Quickstop hinged arms - keeping everything close at hand.

- Quick
- Precise
- Flexible
- Universal
- Robust

Ordering Information

Quickstop hinged arm 1	Order no.
Length 1067 mm, contains:	2M 85 705
3 x hose clamps (2 pcs.)	
1 x hose clamps, coaxial (2 pcs.)	
Load capacity: max. 2 kg	
Quickstop hinged arm 2	
Length 1078 mm, contains:	2M 85 706
2 x hose clamps (2 pcs.)	
1 x hose clamps, coaxial (double)	
1 x hose clamp with ball	
Load capacity: max. 2 kg	
Optional accessories	
Holder for Y-piece (for ball mounting)	84 09 746
Quickstop hinged arm Pod	
Length 460 mm	2M 86 464
For the fixation of Dräger monitoring pods.	
Load capacity: max. 1,5 kg	
Optional accessories	
Quickstop set for the fixation of MultiMed and NeoMed pods.	MP 00 653
Quickstop hinged arm, short	
Length 403 mm, contains: 2M 86 463	2M 86 463
1 x universal tube holder plate for several cable- and tube dimensions	
Load capacity: max. 4 kg	
Spare Parts	
Hose clamps, coaxial (double)	MP00493
Hose clamp (double)	8406138
Hose clamp with ball	8409841
Holder for Y-piece	8409746
Universal tube holder plate (for ball mounting)	2M86465
Quickstop set for the fixation of MultiMed and NeoMed pods.	MP00653

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

Locate your Regional
Sales Representative at:
www.draeger.com/contact



The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).