



# Savina® & More

---

**Dräger**

Technology for Life



# Improving acute care by considering every detail

---

Dear customer,

State-of-the art medical devices, innovative digital solutions, comprehensive services as well as quality consumables and accessories – these are the cornerstones that make it easier for you to provide dedicated healthcare to your patients.

Providing consumables and accessories that enable high-quality therapy, ensure staff satisfaction and at the same time optimise the cost of care is what we do. Despite their apparent simplicity, these products play an integral role in administering successful and effective medical treatment by connecting medical devices with patients. Whether it's supporting infection prevention and control with a broad portfolio of disposable and reusable accessories, or improving workplace efficiency by streamlining tasks with intelligently designed and easy-to-use accessories – our consumables and accessories can help you overcome your challenges. You can focus on your patients, knowing that you can rely on the quality and performance of the products you use in your daily routines.

Moreover, we always keep patient comfort and safety in mind to make patient and staff experience as good as possible. And we go one step further. That's why we invest heavily in the design and development of sustainable accessories that help us shape a more sustainable future.


Improved clinical, operational, and economic outcomes through the entire patient pathway are our common priority – our goal is to support you all the way with our consumables and accessories as they are **always an essential part of your acute care.**



# Protecting the respiratory system

As non-invasive as possible, as invasive as necessary: Discover the variability and diversity of our treatment tools for use along the respiration pathway. They enable administration of protective mechanical ventilation therapy in your ICU.

pathway-2020-en



## PREVENT

Avoiding intubation as long as possible helps to decrease the possibility of patients developing ventilator-associated-pneumonia (VAP).

## STABILISE

Stabilising your patient is crucial to protecting their lungs from potential conditions such as ventilator-associated lung injury (VALI), which could result in acute lung injury (ALI) or acute respiratory distress syndrome (ARDS).

## WEAN

Weaning your patient efficiently reduces incidences of ALI or ARDS.

## RECOVER

Ensuring successful weaning and avoiding reintubation improves your patient's recovery.



- A Savina®
- B Savina® 300
- C Savina® 300 Classic
- D Savina® 300 Select
- E Savina® 300 NIV



# On this way you'll find all your accessories for the Savina® family

<b>1. APPLICATIONS FOR RETROFITTING</b>	<b>08</b>	<b>6. VENTILATION ACCESSORIES</b>	<b>23</b>
<b>2. HOSES FOR GAS SUPPLY</b>	<b>10</b>	6.1 Disposable breathing circuits	23
<b>3. ACCESSORIES FOR BASIC UNIT</b>	<b>11</b>	6.2 Reusable breathing circuits/sets and related accessories	29
3.1 Hinged arms	12	6.3 Breathing system filters and HMEs	37
3.2 Sensors	13	6.4 Catheter mounts and airway connectors	41
3.3 Expiratory valves and filters	14	6.5 Non-invasive patient interfaces and accessories	43
3.4 CO <sub>2</sub> Monitoring	15	6.6 High-Flow Therapy Nasal Cannula and accessories	49
3.5 Trolley and accessories for the Trolley	16	6.7 Pack2Go®	50
<b>4. BREATHING GAS HUMIDIFIERS</b>	<b>17</b>	6.8 Set2go	50
4.1 Humidifier Aquapor H300	17	6.9 Accessories for manual Ventilation	51
4.2 Fisher & Paykel Humidifier and related Accessories	18		
<b>5. MEDICATION NEBULIZER</b>	<b>21</b>		
5.1 Pneumatic medication nebulizer	21		
5.2 Aerogen Nebulizer	22		

## APPLICATIONS FOR RETROFITTING

### 1. APPLICATIONS FOR RETROFITTING

		SAVINA® 300	SAVINA® 300 CLASSIC	SAVINA® 300 SELECT	SAVINA® 300 NIV	SAVINA®
<b>Software update SW 05.02.01</b>	<b>8422604</b>	■	■	■		
<b>Software update SW 3.11</b>	<b>8421497</b>				■	
<b>PC-APRV*</b>	<b>8420887</b>	■	■	■		
<b>PC-BIPAP**/PC-SIMV+</b>	<b>8414131</b>	■	■	■	■	
<b>PC-AC</b>	<b>8420888</b>	■	■	■		
<b>VC-MMV*</b>	<b>8420889</b>	■	■	■		
<b>AutoFlow®</b>	<b>8422700</b>	■	■	■		
<b>AutoFlow® with software SW 3.11</b>	<b>8422700</b>				■	
<b>NIV with software SW 3.11</b>	<b>8415595</b>				■	
<b>Non-Invasive Ventilation, NIV</b>	<b>8422628</b>	■	■	■		
<b>MonitoringPlus (loops, trends)***</b>	<b>8420879</b>	■	■	■		
<b>PaediatricPlus*</b>	<b>8421313</b>	■	■	■		
<b>Special maneuvers</b>	<b>8421315</b>	■	■	■		
<b>Automatic tube compensation, ATC®*</b>	<b>8422626</b>	■	■	■		
<b>O<sub>2</sub> therapy</b>	<b>8422627</b>	■	■	■		
<b>CO<sub>2</sub> monitoring kit including sensor and software option***</b>	<b>8420766</b>	■	■	■		
<b>Nurse call kit including connector</b>	<b>8420760</b>	■	■	■		
<b>Sound reduction kit</b>	<b>8421038</b>	■	■	■		

\*Retrofit applications are not available for the product variant Savina® 300 Classic and Savina® 300 NIV.

\*\*Trademark used under license.

\*\*\*Retrofit application are not available for the product variant Savina® 300 NIV.

## APPLICATIONS FOR RETROFITTING

<b>A</b>	<b>External batteries set for Savina® 300 trolley, consisting of battery box and 2 batteries</b>	<b>8417191</b>
	<b>External batteries set for compact trolley, consisting of battery box and 2 batteries</b>	<b>8422736</b>
<b>B</b>	<b>Batteries 12 V/17 Ah (2 batteries included without box)</b>	<b>8422634</b>
	<b>DC battery cable kit for Savina®</b>	<b>8421430</b>
<b>C</b>	<b>Transport supply unit (TSU) including transport coupling for Savina® 300 trolley (not for compact trolley)</b>	<b>8419911</b>
	<b>Transport coupling for transport supply unit (mounted to the trolley)</b>	<b>G93137</b>
	<b>Bed coupling for Savina® 300 trolley (not for compact trolley)</b>	<b>8419800</b>
<b>D</b>	<b>Bedside rail clamp for bed coupling</b>	<b>8419801</b>

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## HOSES FOR GAS SUPPLY

### 2. HOSES FOR GAS SUPPLY

<b>A Dräger supply hoses, selectable by the adjacent configuration number</b>	<b>8601697</b>
Design, color coding and connection thread of all Dräger supply hoses comply with the respective legal regulations in the indicated countries.	

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## ACCESSORIES FOR BASIC UNIT

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV

### 3. ACCESSORIES FOR BASIC UNIT

	<b>Monitor bracket Savina® 300</b>	<b>8418819</b>	■	■	■
	Mounting for patient monitors Delta/DeltaXL/Vista with GCX-VESA-Slide-Adapter.				
<b>A</b>	<b>CS connection hose holder</b>	<b>2M85446</b>	■	■	■
	<b>MEDIBUS cable</b>	<b>8306488</b>	■	■	■
	<b>Plug for nurse call connection</b>	<b>1846248</b>	■	■	■
	<b>Profile Definition, MEDIBUS.X®</b>	<b>9052608</b>	■	■	■
	<b>Rules and Standards for Implementation, MEDIBUS.X®</b>	<b>9052607</b>	■	■	■
	<b>Protocol Definition, Dräger RS 232 MEDIBUS</b>	<b>9028258</b>	■	■	■
	<b>Software Protocol MEDIBUS for Savina® 300</b>	<b>9052411</b>	■	■	■
<b>B</b>	<b>Earthing cable, 0.8 m (2.6 ft)</b>	<b>8301348</b>	■	■	■
	<b>Earthing cable, 3.2 m (10.5 ft)</b>	<b>8301349</b>	■	■	■



## ACCESSORIES FOR BASIC UNIT

### 3.1 HINGED ARMS

Keeping everything close at hand. To leave your hands free to care for the patient, all the joints of the Quickstop hinged arm can be fixed into place and released with just one quick hand – movement. All cables, lines and hose systems can be arranged in such a way as to more or less completely prevent accidental disconnection.

- |          |   |                |
|----------|---|----------------|
| <b>A</b> | <b>Hinged arm Quickstop 2, with universal ball holder system for optional accessories length 858 mm, max. load 0.4 kg, includes hose clamps: MP00493 (1×), 8404761 (2×), 8409841 (1×)</b> | <b>2M85706</b> |
| <b>B</b> | <b>Hinged arm 240°, 3 hinges, fixable, ball holder system, length 1,155 mm, max. load 0.3 kg, includes hose clamps: MP00493 (1×), 8409841 (1×)</b>  | <b>8409609</b> |

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## ACCESSORIES FOR BASIC UNIT

### 3.2 SENSORS

#### FLOW SENSORS BY DRÄGER

##### Hot-wire anemometry:

##### Accurate and intelligent flow measurement

Spirolog®, Infinity® ID, SpiroLife® and neonatal flow sensor ISO15 are based on proven hot-wire anemometry techniques. In this field we have gained experience and expertise through decades of development and manufacturing. Due to extremely fast and accurate measurement results, as well as extremely low inspiratory and expiratory resistances, hot-wire anemometry is clearly superior to other methods.

The neonatal flow sensor ISO15 offers a customized solution for the sole exception of the specific needs of neonates in intensive care.

The sterilizable neonatal flow sensor ISO15 and SpiroLife® combine the advantages of hot-wire technology with the possibility of sterilizability. Therefore, they integrate seamlessly into your hospital's sterilization processes.

<b>A</b>	<b>Flow sensor Spirolog®, PC/ABS, 5 pcs.</b>	<b>8403735</b>	■	■	■	■
<b>B</b>	<b>Flow sensor SpiroLife®, PPSU, sterilizable</b>	<b>MK01900</b>	■	■	■	■
<b>C</b>	<b>Flow sensor Infinity® ID, PC/ABS, 5 pcs.</b>	<b>6871980</b>	■	■	■	■

#### OXYGEN SENSORS BY DRÄGER

##### Oxygen for Life

Oxygen sensors have the important task of measuring and monitoring the inspiratory and/or expiratory O<sub>2</sub> concentration in anesthesia and respiratory devices. Oxygen sensors from Dräger satisfy the highest national and international requirements.

Furthermore, they have been extensively tested and fine-tuned for proper functional performance with Dräger devices and systems.

Dräger's years of experience in the field of sensors, and their high quality has even convinced NASA. Dräger O<sub>2</sub> sensors are electrochemical transducers used to measure O<sub>2</sub> partial pressure. They operate on the principle of a galvanic cell. Examples of locations where the sensors are used include the Fabius® family, Evita® family and the Primus®.

To take into account the high value of oxygen measurement, we only recommend using O<sub>2</sub> sensors from Dräger.

<b>D</b>	<b>O<sub>2</sub> sensor OxyTrace VE, for Savina® ventilator</b>	<b>MX01049</b>	■	■	■	■
----------	---	----------------	---	---	---	---

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



Not all articles are available worldwide.

## ACCESSORIES FOR BASIC UNIT

### 3.3 EXPIRATORY VALVES AND FILTERS

#### EXPIRATORY VALVES WITHOUT INFINITY® ID

	<b>Expiratory valve, for Savina®, reusable</b>	<b>8413660</b>				
<b>A</b>	<b>Expiratory valve, for Savina® 300, reusable</b>	<b>8417050</b>	■	■	■	■
<b>B</b>	<b>Replcement diaphragm for reusable expiratory valve</b>	<b>8413661</b>	■	■	■	■
<b>C</b>	<b>Expiratory valve, disposable, 10 pcs.</b>	<b>MP01061/ MP20109</b>	■	■	■	■

#### EXPIRATORY VALVES WITH INFINITY® ID

	<b>Expiratory valve, adult, reusable</b>	<b>8416750</b>	■	■	■	■
<b>D</b>	<b>Infinity® ID, expiratory valve, disposable, 10 pcs.</b>	<b>MP01060/ MP20108</b>	■	■	■	■
<b>E</b>	<b>Microfilter (Air intake filter)</b>	<b>6737545</b>	■	■	■	■
<b>F</b>	<b>Dust filter set</b>	<b>8417898</b>	■	■	■	■
<b>G</b>	<b>Expiratory filter, disposable, 20 pcs.</b>	<b>MP01781</b>	■	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



**ACCESSORIES FOR BASIC UNIT**

**3.4 CO<sub>2</sub> MONITORING**

CO<sub>2</sub> cuvettes are designed for the monitoring of CO<sub>2</sub> concentration during patient ventilation.

**DISPOSABLE**

<b>A</b>	<b>CO<sub>2</sub> cuvette disposable, adult, 10 pcs.</b>	<b>MP01062</b>	■ ■
<b>B</b>	<b>CO<sub>2</sub> cuvette disposable, paediatric, 10 pcs.</b>	<b>MP01063</b>	■ ■



**REUSABLE**

<b>C</b>	<b>CO<sub>2</sub> cuvette reusable, adult</b>	<b>6870279</b>	■ ■
<b>D</b>	<b>CO<sub>2</sub> cuvette reusable, paediatric</b>	<b>6870280</b>	■ ■

**ACCESSORIES**

<b>E</b>	<b>Parking holder for CO<sub>2</sub> sensor</b>	<b>8416352</b>	■ ■
	<b>CO<sub>2</sub> sensor Savina® 300</b>	<b>8420757</b>	■ ■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®





## ACCESSORIES FOR BASIC UNIT

### 3.5 TROLLEY AND ACCESSORIES FOR THE TROLLEY

<b>A</b>	<b>Trolley for Savina® 300</b>	<b>8417890</b>	■	■	■
	<b>Compact trolley for Savina® 300</b>	<b>8424046</b>	■	■	■
<b>ACCESSORIES FOR TROLLEYS</b>					
	<b>Cylinder holder set</b>	<b>8411970</b>	■	■	■
	<b>Kit of counterweights</b>	<b>8418211</b>	■	■	■
<b>B</b>	<b>Gas cylinder holder for Savina® 300 trolley</b>	<b>G93110</b>	■	■	■
	<b>Gas cylinder holder for Savina® 300 compact trolley</b>	<b>2601063</b>	■	■	■
	<b>Oxett head, configurable via OxyLine</b>	<b>MP00100</b>	■	■	■
	<b>Battery cable</b>	<b>8418810</b>	■	■	■
	<b>Ferrite for battery cable</b>	<b>1857053</b>	■	■	■
<b>C</b>	<b>External batteries set for Savina® 300 trolley, consisting of battery box and 2 batteries</b>	<b>8417191</b>	■	■	■
	<b>External batteries set for compact trolley, consisting of battery box and 2 batteries</b>	<b>8422736</b>	■	■	■
<b>D</b>	<b>Batteries 12 V/17 Ah (2 batteries included without box)</b>	<b>8422634</b>	■	■	■
	<b>DC battery cable kit for Savina®</b>	<b>8421430</b>	■	■	■
	<b>Transport coupling for transport supply unit (mounted to the trolley)</b>	<b>G93137</b>	■	■	■
	<b>Bed coupling for Savina® 300 trolley (not for compact trolley)</b>	<b>8419800</b>	■	■	■
<b>E</b>	<b>Bedside rail clamp for bed coupling</b>	<b>8419801</b>	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



**BREATHING GAS HUMIDIFIERS**

**4. BREATHING GAS HUMIDIFIERS**



**4.1 HUMIDIFIER AQUAPOR H300**

	<b>Aquapor H300 – Configuration</b>	<b>MP01000</b>				
<b>A</b>	<b>Aquapor H300 Basic unit 230V</b>	<b>MP07101</b>				
	<b>Aquapor H300 Basic unit 115V</b>	<b>MP07102</b>				
<b>ACCESSORIES FOR HUMIDIFIER AQUAPOR H300</b>						
	<b>VentStar® Aquapor, dual heated, A/P, 10 pcs.</b>	<b>MP07130</b>				
<b>B</b>	<b>VentStar® Aquapor, heated, A/P, 10 pcs.</b>	<b>MP07131</b>				
<b>C</b>	<b>VentStar® Aquapor, dual heated, N plus, 10 pcs.</b>	<b>MP07150</b>				
<b>D</b>	<b>VentStar® Aquapor, heated, N plus, 10 pcs.</b>	<b>MP07151</b>				
	<b>HI-Flow Star Kit Aquapor A, heated inspiratory breathing circuit for the H300 humidifier, heated inspiratory hose, humidifier chamber, breathing hose when used with a ventilator, O<sub>2</sub> hose when used with an O<sub>2</sub> blender, bed clip, disposable, 10 pcs.</b>	<b>MP17090</b>				
<b>E</b>	<b>HI-Flow Star Valve Kit w/ Pop-Up Valve and Connector 22 mm, disposable, 15 pcs.</b>	<b>MP05507</b>				
<b>F</b>	<b>H300 Heater wire adapter, insp</b>	<b>MP07106</b>				
<b>G</b>	<b>H300 Heater wire adapter, dual</b>	<b>MP07107</b>				
	<b>H300 Temp. sensor cable*, 160 cm</b>	<b>MP07108</b>				

\*Tradeware

Not all articles are available worldwide.

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## BREATHING GAS HUMIDIFIERS

### 4.2 FISHER & PAYKEL HUMIDIFIER AND RELATED ACCESSORIES

<b>F&amp;P configuration for MR810 and MR850</b>		<b>8414144</b>	■	■	■	■
<b>A</b>	<b>Humidifier chamber VentStar®, disposable, 10 pcs.</b>	<b>MP00234</b>	■	■	■	■
<b>B</b>	<b>Humidifier chamber, 10 pcs.</b>	<b>MP02605</b>	■	■	■	■
<b>C</b>	<b>Humidifier chamber* for F&amp;P MR 290, disposable, 10 pcs.</b>	<b>8418282</b>	■	■	■	■
<b>D</b>	<b>Humidifier chamber* for F&amp;P MR 340, reusable, N/P</b>	<b>8411047</b>	■	■	■	■
<b>E</b>	<b>Humidifier chamber* for F&amp;P MR 370, reusable, A</b>	<b>8412217</b>	■	■	■	■

\*Tradeware

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



**BREATHING GAS HUMIDIFIERS**

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®

<b>A</b>	<b>Temperature probe/flow sensor* for F&amp;P MR 850, for all breathing circuits, 1.1 m</b>	<b>8414989</b>	■	■	■	■
<b>B</b>	<b>Temperature probe/flow sensor* for F&amp;P MR 850, for all breathing circuits, 1.5 m</b>	<b>8414966</b>	■	■	■	■
<b>C</b>	<b>Hose heater adapter* for F&amp;P MR 850, for reusable breathing circuit</b>	<b>8414968</b>	■	■	■	■
<b>D</b>	<b>Hose heater adapter for F&amp;P MR 850, for single-heated disposable breathing circuit</b>	<b>8414992</b>	■	■	■	■
<b>E</b>	<b>Hose heater adapter for F&amp;P MR 850, for dual-heated disposable breathing circuit</b>	<b>MP01006</b>	■	■	■	■
<b>F</b>	<b>Draw wire* for hose heater 1.5 m, for reusable breathing circuit</b>	<b>8411050</b>	■	■	■	■
<b>G</b>	<b>Temperature probe adapter* 900 MR 532, for F&amp;P humidifiers</b>	<b>MP01007</b>	■	■	■	■

\*Tradeware



Not all articles are available worldwide.


## BREATHING GAS HUMIDIFIERS

### FISHER & PAYKEL HUMIDIFIER MOUNTING MATERIAL

<b>A</b>	<b>Rail clamp* for F&amp;P humidifier and Seattle PAP bubble box</b>	<b>8411074</b>	■	■	■	■
<b>B</b>	<b>Water bag holder* 900 MR 290</b>	<b>MP01008</b>	■	■	■	■
	<b>Universal holder* with DIN rail</b>	<b>G93140</b>	■	■	■	■
	<b>Humidifier holder* for chassis, swiveling</b>	<b>G93111</b>	■	■	■	■

\*Tradeware

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®




# MEDICATION NEBULIZER

## 5. MEDICATION NEBULIZER

### 5.1 PNEUMATIC MEDICATION NEBULIZER

<b>A</b>	<b>Pneum. Medication Nebulizer</b>	<b>8412935</b>
<b>B</b>	<b>Nebulizer hose</b>	<b>8412985</b>

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## MEDICATION NEBULIZER

### 5.2 AEROGEN NEBULIZER

**Aerogen** works with an electronic micropump. Aerogen is ideally suited to nebulizer therapy requirements for inline mechanical ventilation and it accommodates all patient age ranges from neonates to adults. It operates without changing patient ventilator parameters, without heating or degrading medications and it can be refilled without interrupting ventilation.

<b>A Configuration for Aerogen Pro and Solo Medication Nebulizer</b>	<b>MP01010</b>				
<b>AEROGEN PRO*</b>					
<b>B Nebulizer unit with filler cap for Aerogen Pro</b>	<b>MP01030</b>				
<b>C T-piece for Aerogen Pro, A, ISO 22 mm, 5 pcs.</b>	<b>MP01031</b>				
<b>D Controller cable for Aerogen Pro</b>	<b>MP01036</b>				
<b>E Filler cap for Aerogen Pro, silicone, 5 pcs.</b>	<b>MP01035</b>				
<b>AEROGEN SOLO*</b>					
<b>F Nebulizer head, 5 pcs.</b>	<b>MP01078</b>				
<b>G Nebulizer head, 10 pcs.</b>	<b>MP01079</b>				
<b>Aerogen S, T-Adapter A, 10 pcs.</b>	<b>MP01085</b>				
<b>Aerogen S, T-Adapter P, 10 pcs.</b>	<b>MP01086</b>				
<b>Aerogen T-Adapter**, 11.5 mm, N, 10 pcs.</b>	<b>MP02689</b>				
<b>Aerogen USB controller replacement w/o AC/DC adapter</b>	<b>MP01095</b>				
<b>Aerogen S, T-Piece, F&amp;P, N, 10 box</b>	<b>MP01077</b>				
<b>Aerogen S, T-Piece, N, 10 pcs.</b>	<b>MP01691</b>				
<b>Aerogen S, T-Piece Plug, 10 pcs.</b>	<b>MP13069</b>				

\*Tradeware

\*\*Made by Drägerwerk AG & Co. KGaA

Not all articles are available worldwide.

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

### 6. VENTILATION ACCESSORIES

Dual-limb systems can be used in all therapy types without changing circuits and humidification systems e.g., in case patients that may need to escalate from NIV to invasive ventilation without delay and vice versa. This does not only support the handling but may also decrease training efforts, costs, and usage of material.

#### 6.1 DISPOSABLE BREATHING CIRCUITS




##### INFINITY ID PRODUCTS FOR ADULTS

<b>A</b>	<b>Breathing circuit Infinity® ID, disposable, coaxial, 1.9 m, 20 pcs., not made with natural rubber latex</b>	<b>MP03375</b>				
<b>B</b>	<b>Breathing circuit Infinity® ID, disposable, coaxial, 2.3 m, 20 pcs., not made with natural rubber latex</b>	<b>MP03377</b>				
<b>C</b>	<b>Breathing circuit Infinity® ID, disposable, water trap, 2.1 m, 10 pcs., not made with natural rubber latex</b>	<b>MP00378</b>				

##### PRODUCTS WITHOUT INFINITY ID FOR ADULTS

<b>D</b>	<b>Breathing circuit Vent Set, disposable, basic, 1.5 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00300</b>				
<b>E</b>	<b>Breathing circuit VentStar®, disposable, basic, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00349</b>				
<b>F</b>	<b>Breathing circuit VentStar®, disposable, basic, 2.5 m, 10 pcs., not made with natural rubber latex</b>	<b>MP00350</b>				



SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

<b>A</b>	<b>Breathing circuit Vent Set, disposable, extendable, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00305</b>
<b>B</b>	<b>Breathing circuit VentStar®, disposable, extendable, 2.2 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00355</b>
<b>C</b>	<b>Breathing circuit VentStar®, disposable, coaxial, water trap, 2.1 m, 10 pcs., not made with natural rubber latex</b> <small>&lt;= Software 3.53</small>	<b>MP00373</b>
<b>D</b>	<b>Breathing circuit VentStar®, disposable, water trap, w/o Luer Lock, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00334</b>
<b>E</b>	<b>Breathing circuit VentStar®, disposable, water trap, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00337</b>
<b>F</b>	<b>Breathing circuit VentStar®, disposable, smoothbore, for MRI, 3 m, 10 pcs., not made with natural rubber latex</b>	<b>MP04800</b>



SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## VENTILATION ACCESSORIES

<b>A Breathing circuit VentStar®, disposable, coaxial, 1.7 m, 20 pcs., not made with natural rubber latex</b>	<b>MP03373</b>
<b>B Breathing circuit VentStar®, disposable, coaxial, 1.9 m, 20 pcs., not made with natural rubber latex</b>	<b>MP03374</b>
<b>C Breathing circuit VentStar®, disposable, coaxial, 2.3 m, 20 pcs., not made with natural rubber latex</b> > Software 3.53	<b>MP03376</b>


SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## VENTILATION ACCESSORIES

<b>A</b>	<b>VentStar® Helix, dual heated, 10 pcs.</b>	<b>MP02606</b>
<b>B</b>	<b>VentStar® Helix, heated, 10 pcs.</b>	<b>MP02607</b>
<b>C</b>	<b>Breathing circuit VentStar® Helix, disposable, N, dual heated, 1.7 m, 10 pcs., not made with natural rubber latex</b>	<b>MP02650</b>
<b>D</b>	<b>Breathing circuit VentStar® Helix, disposable, N, heated, 1.7 m, 10 pcs., not made with natural rubber latex</b>	<b>MP02608</b>
<b>E</b>	<b>Double connector 11ID/15OD, disposable, 10 pcs.</b>	<b>MP02655</b>

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®





## VENTILATION ACCESSORIES

### DISPOSABLE BREATHING CIRCUITS FOR PAEDIATRICS

<b>A</b>	<b>Breathing circuit VentStar®, disposable, P, basic, w/o Luer Lock, 1.5 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00338</b>
<b>B</b>	<b>Breathing circuit VentStar®, disposable, P, basic, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00351</b>
<b>C</b>	<b>Breathing circuit VentStar®, disposable, P, coaxial, 1.5 m, 20 pcs., not made with natural rubber latex</b>	<b>MP00379</b>
<b>D</b>	<b>Breathing circuit VentStar®, disposable, P, water trap, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00361</b>
<b>E</b>	<b>Breathing circuit VentStar®, disposable, P, water trap, w/o Luer Lock, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00362</b>

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®




## VENTILATION ACCESSORIES

### DISPOSABLE BREATHING CIRCUITS FOR NEONATES

<b>A Breathing circuit VentStar®, disposable, N, basic, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00353</b>
<b>B Breathing circuit VentStar®, disposable, N, water trap, 1.8 m, 25 pcs., not made with natural rubber latex</b>	<b>MP00363</b>

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®




**VENTILATION ACCESSORIES**

**6.2 REUSABLE BREATHING CIRCUITS/SETS AND RELATED ACCESSORIES**

**BLUESSET VENTILATION**

<b>A</b>	<b>BlueSet Ventilation, N</b>	<b>8412344</b>
	Consists of:	
	<b>1</b> Y-piece, reusable, N, 90°, with suction opening	<b>8403075</b>
	<b>2</b> Corrugated hose, reusable, N/P, 14 cm (5,5 in), connector 10 mm/10 mm, Ø 10 mm	<b>8403333</b>
	<b>3</b> Temperature sensor mount, reusable, N	<b>8403056</b>
	<b>4</b> Water trap with outlet, reusable N/P	<b>8404760</b>
	<b>5</b> Breathing hose, silicone, reusable, N/P, 60 cm, connector 10 mm/10 mm, Ø 10 mm	<b>2165848</b>
	<b>6</b> Water trap, reusable, N/P	<b>8409627</b>
	<b>7</b> Red cap, 7 mm	<b>8401645</b>
	<b>8</b> Breathing hose, silicone, reusable, N/P, 40 cm, connector 10 mm/22 mm, Ø 10 mm	<b>2165856</b>
	<b>9</b> Breathing hose, silicone, reusable, N/P, 60 cm, connector 10 mm/22 mm, Ø 10 mm	<b>2165821</b>
	<b>Dummy Plug</b>	<b>8410049</b>




SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®

## VENTILATION ACCESSORIES

<b>A BlueSet Ventilation, P/A</b>	<b>8413146</b>
Consists of:	
<b>1 Adapter, reusable, 22 mm/22 mm</b>	<b>M25647</b>
<b>2 Adapter, reusable, 22 mm/15 mm, 90°</b>	<b>8290274</b>
<b>3 Hose clamp, 10 pcs.</b>	<b>8403566</b>
<b>4 Y-piece with sensor mount, reusable, A</b>	<b>8405435</b>
<b>5 Breathing hose, silicone, reusable, A, 35 cm, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2165619</b>
<b>6 Breathing hose, silicone, reusable, A, 60 cm, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2165627</b>
<b>7 Water trap, reusable, A</b>	<b>8404985</b>
<b>Dummy Plug</b>	<b>8410049</b>


SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## VENTILATION ACCESSORIES

<b>A</b>	<b>BlueSet Ventilation HME, P/A</b>	<b>8412860</b>
	Consists of:	
<b>①</b>	<b>Breathing hose, silicone, reusable, A, 1.5 m, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2166046</b>
<b>②</b>	<b>Y-piece, polyamide, reusable, straight, A</b>	<b>M25650</b>



SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

	Silicone breathing hose (adult) with standard sleeves 22 mm/22 mm	
<b>A</b>	<b>Breathing hose, silicone, reusable, A, 35 cm, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2165619</b>
<b>B</b>	<b>Breathing hose, silicone, reusable, A, 60 cm, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2165627</b>
<b>C</b>	<b>Breathing hose, silicone, reusable, A, 1.1 m, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2166038</b>
<b>D</b>	<b>Breathing hose, silicone, reusable, A, 1.5 m, connector 22 mm/22 mm, Ø 22 mm</b>	<b>2166046</b>
	Silicone breathing hose (paediatric) with standard sleeves 10 mm on both sides and 10 mm inner diameter	
<b>E</b>	<b>Breathing hose, silicone, reusable, N/P, 60 cm, connector 10 mm/10 mm, Ø 10 mm</b>	<b>2165848</b>

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®





## VENTILATION ACCESSORIES

Silicone breathing hose (paediatric) with standard sleeves with one 10 mm and one 22 mm sleeve and 10 mm inner diameter

<b>A</b>	<b>Breathing hose, silicone, reusable, N/P, 40 cm, connector 10 mm/22 mm, Ø 10 mm</b>	<b>2165856</b>				
<b>B</b>	<b>Breathing hose, silicone, reusable, N/P, 60 cm, connector 10 mm/22 mm, Ø 10 mm</b>	<b>2165821</b>				
<b>C</b>	<b>Breathing hose, silicone, reusable, N/P, 1.1 m, connector 10 mm/22 mm, Ø 10 mm</b>	<b>2165651</b>				

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

<b>A</b>	<b>Y-piece, polyamide, reusable, straight, A</b>	<b>M25650</b>
<b>B</b>	<b>Y-piece for F&amp;P, reusable, A*</b>	<b>8412888</b>
<b>C</b>	<b>Y-piece with sensor mount, reusable, A</b>	<b>8405435</b>
<b>D</b>	<b>Y-piece, reusable, N, 90°, with suction opening</b>	<b>8403075</b>
<b>E</b>	<b>Temperature probe adapter 900 MR 532, for F&amp;P humidifiers*</b>	<b>MP01007</b>

\*Tradeware

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

<b>A</b>	<b>Adapter, reusable, 22 mm/22 mm</b>	<b>M25647</b>	■	■	■	■
	<b>Adapter straight, reusable, 18 mm/11 mm, A</b>	<b>8403076</b>	■	■	■	■
	<b>Adapter, reusable, 22 mm/22 mm, 90°</b>	<b>M25649</b>	■	■	■	■
<b>B</b>	<b>Adapter reusable, 11/11, N/P</b>	<b>8409897</b>	■	■	■	■
	<b>Temperature sensor mount, reusable, N</b>	<b>8403056</b>	■	■	■	■
<b>C</b>	<b>Water trap, reusable, A</b>	<b>8404985</b>	■	■	■	■
<b>D</b>	<b>Water trap, reusable, N/P</b>	<b>8409627</b>	■	■	■	■
	<b>Red cap, 7 mm</b>	<b>8401645</b>	■	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



Not all articles are available worldwide.

## VENTILATION ACCESSORIES

		SAVINA® 300	SAVINA® 300 CLASSIC	SAVINA® 300 SELECT	SAVINA® 300 NIV	SAVINA®
<b>Set of caps, set of 5 pcs.</b>	<b>8402918</b>	■	■	■	■	■
<b>Corrugated hose</b>	<b>8402041</b>	■	■	■	■	■
<b>Corrugated hose, reusable, N/P, 14 cm (5,5 in), connector 10 mm/10 mm, Ø 10 mm</b>	<b>8403333</b>	■	■	■	■	■
<b>Hose clamp, 10 pcs.</b>	<b>8403566</b>	■	■	■	■	■

**VENTILATION ACCESSORIES**

**6.3 BREATHING SYSTEM FILTERS AND HMEs**

**ELECTROSTATIC FILTER AND HME**

The **Filter/HME TwinStar® Plus** from Dräger is a combination of an efficient HME and a high-performance breathing system filter, thereby significantly contributing to infection prophylaxis in ventilation treatment.

- Combination of efficient HME and high-performance breathing system filter
- High retention rates for bacteria and viruses
- High return of moisture
- Sampling connector with tethered cap
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and blue color coding for quick identification

<b>A</b>	<b>Filter + HME TwinStar® 90 Plus, disposable, ~90 ml dead space, 100 pcs.</b>	<b>MP05800</b>	■	■	■	■
<b>B</b>	<b>Filter + HME TwinStar® 60A Plus, disposable, ~60 ml dead space, angled, 100 pcs.</b>	<b>MP05810</b>	■	■	■	■
<b>C</b>	<b>Filter + HME TwinStar® 55 Plus, disposable, ~55 ml dead space, 100 pcs.</b>	<b>MP05805</b>	■	■	■	■
<b>D</b>	<b>Filter + HME TwinStar® 25 Plus, disposable, ~25 ml dead space, 100 pcs.</b>	<b>MP05815</b>	■	■	■	■
<b>E</b>	<b>Filter + HME TwinStar® 9 Plus, disposable, ~9 ml dead space, 100 pcs.</b>	<b>MP05820</b>	■	■	■	■

**MECHANICAL FILTER AND HME**

<b>F</b>	<b>Mechanical filter + HME TwinStar® HEPA Plus, disposable, ~90 ml dead space, 100 pcs.</b>	<b>MP05801</b>	■	■	■	■
----------	---	----------------	---	---	---	---

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



Not all articles are available worldwide.

## VENTILATION ACCESSORIES

### ELECTROSTATIC FILTER

The **CareStar® Plus** breathing system filter significantly contributes to infection prophylaxis in anesthesia and ventilation treatment through high retention rates.

- High-performance electrostatic filter
- Patient-side or device-side use
- Sampling connector with tethered cap
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and yellow color coding for quick identification

<b>A</b>	<b>Electrostatic filter CareStar® 35 Plus, disposable, ~35 ml dead space, 100 pcs.</b>	<b>MP05755</b>
<b>B</b>	<b>Electrostatic filter CareStar® 20 Plus, disposable, ~20 ml dead space, 100 pcs.</b>	<b>MP05770</b>

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## VENTILATION ACCESSORIES

### MECHANICAL FILTER

The **SafeStar® Plus** mechanical high performance breathing system filter shows excellent retention rates, thereby significantly contributing to infection prophylaxis in anesthesia and ventilation treatment.

- High-performance mechanical HEPA filter
- High hydrophobicity
- Patient-side or device-side use
- Sampling connector with tethered cap
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and red color coding for quick identification

<b>A</b>	<b>Mechanical filter SafeStar® 90 Plus, disposable, ~90 ml dead space, 100 pcs.</b>	<b>MP05785</b>	■	■	■	■
<b>B</b>	<b>Mechanical filter SafeStar® 60A Plus, disposable, ~60 ml dead space, angled, 100 pcs.</b>	<b>MP05795</b>	■	■	■	■
<b>C</b>	<b>Mechanical filter SafeStar® 55 Plus, disposable, ~55 ml dead space, 100 pcs.</b>	<b>MP05790</b>	■	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

### HME

**HumidStar® Plus** moisture exchangers from Dräger partially take over the function of the upper airways and efficiently humidify the inspired air.

- High return of moisture
- Sampling connector with tethered cap
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and green color coding for quick identification

<b>A</b>	<b>HME HumidStar® 55 Plus, disposable, ~55 ml dead space, 100 pcs.</b>	<b>MP05730</b>	■ ■ ■ ■ ■
<b>B</b>	<b>HME HumidStar® 25 Plus, disposable, ~25 ml dead space, 100 pcs.</b>	<b>MP05735</b>	■ ■ ■ ■ ■
<b>C</b>	<b>HME HumidStar® Trach Plus, disposable, 6 ml dead space, 100 pcs.</b>	<b>MP05750</b>	■ ■ ■ ■ ■

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



**VENTILATION ACCESSORIES**


**6.4 CATHETER MOUNTS AND AIRWAY CONNECTORS**

- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages
- Extendable product provides a lightweight and flexible solution
- Double swivel connector increases ease of use
- Double cap to maintain PEEP during suctioning or bronchoscopy
- Each catheter mount and airway connector is equipped with a red safety cap so that the system is completely closed and protected until directly prior to use on the patient
- 22F/15M connector on device side as a two-in-one solution
- Standardized connectors for safe connection to other components
- Luer Lock port for gas sampling
- Disposable convenience enhances workflow
- PVC-free design of gas leading components eliminates risk of softeners such as DEHP

<b>A</b>	<b>Catheter mount ErgoStar® CM 40, disposable, 10 cm, double swivel elbow, 50 pcs.</b>	<b>MP01840</b>	■	■	■
<b>B</b>	<b>Catheter mount ErgoStar® CM 45, disposable, 15 cm, double swivel elbow, 50 pcs.</b>	<b>MP01845</b>	■	■	■
<b>C</b>	<b>Catheter mount ErgoStar® CM 50, disposable, 10 – 16 cm, extendable, 50 pcs.</b>	<b>MP01850</b>	■	■	■
<b>D</b>	<b>Catheter mount ErgoStar® CM 55, disposable, 15 cm, straight, 50 pcs.</b>	<b>MP01855</b>	■	■	■
<b>E</b>	<b>Catheter mount ErgoStar® CM 60, disposable, 10 cm, small, double swivel elbow, 50 pcs.</b>	<b>MP01860</b>	■	■	■
<b>F</b>	<b>Airway connector ErgoStar® AC 90, disposable, double-swivel elbow w/ double cap, 100 pcs.</b>	<b>MP01890</b>	■	■	■
<b>G</b>	<b>Airway connector ErgoStar® AC 95, disposable, w/ Luer Lock sampling port, 100 pcs.</b>	<b>MP01895</b>	■	■	■

Not all articles are available worldwide.

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®

### COMBINATIONS OF FILTER AND TUBE EXTENTION

<b>A CombiStar® Plus Filter + HME straight, combination of airway connector MP01845 and filter TwinStar® 55 Plus (MP05805), 25 pcs.</b>	<b>MP12060</b>
<b>B CombiStar® Plus Filter + HME flex, combination of flexible tube extention MP01850 and filter TwinStar® 55 Plus (MP05805), 25 pcs.</b>	<b>MP12061</b>
<b>C CombiStar® Plus Filter + HME HEPA flex, combination of flexible tube extention MP01850 and filter TwinStar® HEPA Plus (MP05801), 25 pcs.</b>	<b>MP12062</b>
<b>D CombiStar® Plus Filter flex, combination of flexible tube extention MP01850 and filter SafeStar® 55 Plus (MP05790), 25 pcs.</b>	<b>MP12063</b>



**VENTILATION ACCESSORIES**

**6.5 NON-INVASIVE PATIENT INTERFACES AND ACCESSORIES**

**NIV MASKS AND ACCESSORIES**

**ClassicStar® NIV oronasal masks** cover the nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing. The mask design ensures high patient comfort combined with anatomical fit and effective seal due to individually adjustable mask cushion. The 6-point headgear provides a stable and secure fixation. Adjustable forehead support and forehead pad allow for distributing the pressure from the bridge of nose towards the forehead for more patient comfort.

ClassicStar® NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option.

All ClassicStar® NIV masks are for single patient use. Sizing gauges for easy mask size identification as well as a pump ball are delivered together with the masks.

<b>A</b>	<b><u>NIV oronasal mask ClassicStar® SE, disposable, size S</u></b>	<b><u>MP01573</u></b>	■	■	■	■
	<b><u>NIV oronasal mask ClassicStar® SE, disposable, size M</u></b>	<b><u>MP01574</u></b>	■	■	■	■
	<b><u>NIV oronasal mask ClassicStar® SE, disposable, size L</u></b>	<b><u>MP01575</u></b>	■	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®




Not all articles are available worldwide.

## VENTILATION ACCESSORIES

Oronasal masks are often the first-line strategy in the initial management of acute respiratory failure as they cover nose and mouth. Nasal masks are often chosen if NIV therapy has to be prolonged due to improved comfort and reduced complications associated with oronasal masks. Similar to ClassicStar® oronasal masks, the nasal mask cushion can be individually adjusted to the patient's face providing high patient comfort and effective seal. The 4-point-headgear ensures a stable and secure fixation reducing mask displacement. All **ClassicStar® NIV nasal masks** are for single patient use. Sizing gauges for easy mask size identification as well as a pump ball are delivered together with the masks.

<b>A NIV nasal mask ClassicStar® non-vented, disposable, size S</b>	<b>MP01623</b>	■ ■ ■ ■
<b>NIV nasal mask ClassicStar® non-vented, disposable, size M</b>	<b>MP01624</b>	■ ■ ■ ■
<b>NIV nasal mask ClassicStar® non-vented, disposable, size L</b>	<b>MP01625</b>	■ ■ ■ ■

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



D-5990-2023

## VENTILATION ACCESSORIES

The **ClassicStar® NIV total-face mask** seals around the perimeter of the face, covering mouth, nose and eyes. These areas are less pressure sensitive and have smoother facial contours compared to the areas a oronasal masks covers (nasal bridge). The mask is suitable for patients with nasal bridge breakdown or facial irregularities. The air filled cushion can be adjusted individually to the patients face, providing high patient comfort and effective seal.

The 5-point headgear which includes a headstrap at the top of the mask provides a stable and secure attachment, thus preventing mask displacement and reducing pressure. Quick release clips are integrated for quick removal of the mask.

Sizing gauges for easy mask size identification as well as a pump ball are delivered together with the masks.

<b>A</b>	<b>NIV total-face mask ClassicStar® SE, disposable, size S</b>	<b>MP02910</b>	■	■	■	■
	<b>NIV total-face mask ClassicStar® SE, disposable, size M</b>	<b>MP02911</b>	■	■	■	■
	<b>NIV total-face mask ClassicStar® SE, disposable, size L</b>	<b>MP02912</b>	■	■	■	■

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®





D-6282-2018

## VENTILATION ACCESSORIES

The **ClassicStar® plus mask** series features a soft and anatomically shaped silicone cushion as sealing interface to the patients face. The mask has an easy to adjust forehead pad with silicone cushion which allows optimal adaption to the sensitive nasal bridge. The mask is available in four sizes (S, M, L, XL) as oronasal and nasal mask. ClassicStar® plus NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option.

<b>A</b>	<b>NIV oronasal mask ClassicStar® plus SE, disposable, size S, 5 pcs.</b>	<b>MP04711</b>	■	■	■	■
	<b>NIV oronasal mask ClassicStar® plus SE, disposable, size M, 5 pcs.</b>	<b>MP04712</b>	■	■	■	■
	<b>NIV oronasal mask ClassicStar® plus SE, disposable, size L, 5 pcs.</b>	<b>MP04713</b>	■	■	■	■
	<b>NIV oronasal mask ClassicStar® plus SE, disposable, size XL, 5 pcs.</b>	<b>MP04714</b>	■	■	■	■
<b>B</b>	<b>NIV nasal mask ClassicStar® plus non-vented, disposable, size S, 5 pcs.</b>	<b>MP04741</b>	■	■	■	■
	<b>NIV nasal mask ClassicStar® plus non-vented, disposable, size M, 5 pcs.</b>	<b>MP04742</b>	■	■	■	■
	<b>NIV nasal mask ClassicStar® plus non-vented, disposable, size L, 5 pcs.</b>	<b>MP04743</b>	■	■	■	■
	<b>NIV nasal mask ClassicStar® plus non-vented, disposable, size XL, 5 pcs.</b>	<b>MP04744</b>	■	■	■	■

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

The **NovaStar® plus NIV** mask series cover nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing. This range of products provides ultimate comfort and an effective seal due to the individual adaption of the silicone-gel cushion with a pliable embedded ring. NovaStar® plus NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option. All NovaStar® plus NIV masks are for multi patient use. Replacement headgears with magnetic clips are available individually (cf. accessories for NIV masks).

**A NIV oronasal mask NovaStar® plus, SE, reusable, size S, 3 pcs.**

**MP05051**

**NIV oronasal mask NovaStar® plus, SE, reusable, size M, 3 pcs.**

**MP05052**

**NIV oronasal mask NovaStar® plus, SE, reusable, size L, 3 pcs.**

**MP05053**

**NIV oronasal mask NovaStar® plus, SE, reusable, size XL, 3 pcs.**

**MP05054**

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



## VENTILATION ACCESSORIES

The **MiniMe2\*** nasal mask is designed for children ages 2 – 12. It features an ergonomic silicone-gel cushion and a shapeable shell mask to fit individual needs and thus ensuring an effective seal and a comfortable fit.

**A Paediatric nasal mask MiniMe2 non-vented, disposable, size S** **MP01555**

**Paediatric nasal mask MiniMe2 non-vented, disposable, size L** **MP01556**


\*Tradeware

### ACCESSORIES FOR NIV MASKS

**B Two-sided pump ball for air management** **MP01590**

**C NovaStar® plus Headgear inclusive magnetic clips, disposable 3 pcs.** **MP05080**

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®




**VENTILATION ACCESSORIES**

**6.6 HIGH-FLOW THERAPY NASAL CANNULA AND ACCESSORIES**

**High-flow therapy** allows to provide warm humidified air with enriched oxygen concentration at flow rates up to 60 l per minute to patients via nasal cannulas. Patient comfort is one of the advantages of High Flow therapy over other treatment forms for patient with compromised respiration and can in many cases successfully avoid using less comfortable or even invasive forms of respiratory support. The HI-Flow Star Nasal Cannula and the inspiratory breathing circuit offer advanced features designed for enhanced patient safety and comfort such as integrated clips in order to provide effective breathing circuit tension relief, a flexible, breathable hose minimizing internal condensation in the high flow nasal cannula as well as soft materials for the prongs and the headgear. Thanks to its symmetrical design, the nasal cannula can be connected from either the left or the right side of the bed. This helps caregivers save time and effort. The HI-Flow Star products are disposable.

<b>A</b>	<b>① HI-Flow Star Nasal Cannula, disposable, adult S, Flow Range &lt; 60 l/min, 25 pcs.</b>	<b>MP05511</b>				
<b>A</b>	<b>② HI-Flow Star Nasal Cannula, disposable, adult M, Flow Range &lt; 60 l/min, 25 pcs.</b>	<b>MP05512</b>				
<b>A</b>	<b>③ HI-Flow Star Nasal Cannula, disposable, adult L, Flow Range &lt; 60 l/min, 25 pcs.</b>	<b>MP05513</b>				
<b>B</b>	<b>HI-Flow Star System, heated inspiratory breathing circuit, consists of: heated inspiratory hose, humidifier chamber, breathing hose when used with a ventilator, O<sub>2</sub> hose when used with an O<sub>2</sub> blender, bed clip, disposable, 25 pcs.</b>	<b>MP05601</b>				
<b>C</b>	<b>HI-Flow Star Kit Aquapor A, heated inspiratory breathing circuit for the H300 humidifier, heated inspiratory hose, humidifier chamber, breathing hose when used with a ventilator, O<sub>2</sub> hose when used with an O<sub>2</sub> blender, bed clip, disposable, 10 pcs.</b>	<b>MP17090</b>				
<b>D</b>	<b>HI-Flow Star Valve Kit w/ Pop-Up Valve and Connector 22 mm, disposable, 15 pcs.</b>	<b>MP05507</b>				

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



Not all articles are available worldwide.

## VENTILATION ACCESSORIES

### 6.7 PACK2GO®

<b>A</b>	<b>Pack2Go® invasive ventilation pack with breathing circuit VentStar®, coaxial, 1.8 m, 10 pcs.</b>	<b>MP20375</b>
<b>B</b>	<b>Pack2Go® invasive ventilation pack with breathing circuit VentStar®, basic, 1.8 m, 10 pcs.</b>	<b>MP20376</b>

### 6.8 SET2GO

<b>C</b>	<b>With Dräger's Set2go solutions we offer you totally customised sets consisting of an individual selection of single-patient-use accessories. Dräger Set2go makes using and ordering consumables and accessories much simpler for you. You just need one article number to get all the items in a single package and accomplish your specific task. Please contact your local Dräger sales organisation to create your individual accessory solution.</b>
----------	---

SAVINA® 300  
 SAVINA® 300 CLASSIC  
 SAVINA® 300 SELECT  
 SAVINA® 300 NIV  
 SAVINA®



## VENTILATION ACCESSORIES

### 6.9 ACCESSORIES FOR MANUAL VENTILATION

<b>A</b>	<b>Test lung SelfTestLung, latex-free, reusable</b>	<b>MP20020</b>
<b>B</b>	<b>Test Lung, reusable, max. 1.5 L, not made with natural rubber latex</b>	<b>8403201</b>

SAVINA® 300  
SAVINA® 300 CLASSIC  
SAVINA® 300 SELECT  
SAVINA® 300 NIV  
SAVINA®



Not all products, features, or services are for sale in all countries. Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to [www.draeger.com/trademarks](http://www.draeger.com/trademarks) to find the current status.

**Corporate Headquarters**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**Manufacturer:**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23542 Lübeck, Germany

**Region Europe**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany  
☎ +49 451 882 0  
✉ +49 451 882 2080  
✉ [info@draeger.com](mailto:info@draeger.com)

**Region Middle East, Africa**

Drägerwerk AG & Co. KGaA  
Branch Office  
P.O. Box 505108  
Dubai, United Arab Emirates  
☎ +971 4 4294 600  
✉ +971 4 4294 699  
✉ [contactuae@draeger.com](mailto:contactuae@draeger.com)

**Region Asia Pacific**

Draeger Singapore Pte. Ltd.  
61 Science Park Road  
The Galen #04-01  
Singapore 117525  
☎ +65 6872 9288  
✉ +65 6259 0398

**Region Central  
and South America**

Dräger Indústria e Comércio Ltda.  
Al. Pucurui - 51 - Tamboré  
06460-100 - Barueri - São Paulo  
☎ +55 (11) 4689-4900  
✉ [relacionamento@draeger.com](mailto:relacionamento@draeger.com)

**Canada**

Draeger Medical Canada Inc.  
2425 Skymark Avenue, Unit 1  
Mississauga, Ontario, L4W 4Y6  
☎ +1 905 212 6600  
Toll-free +1 866 343 2273  
✉ +1 905 212 6601  
✉ [Canada.support@draeger.com](mailto:Canada.support@draeger.com)



Locate your Regional Sales  
Representative at:  
[www.draeger.com/contact](http://www.draeger.com/contact)

