

Dräger SelfTestLung™ Consumables and Accessories

To support your clinical workflow our Dräger SelfTestLung™ combines high-quality materials and an intelligent design. Tested and approved for compatibility with a wide range of our Dräger ventilators. To support you in reducing the risk of infections our SelfTestlung™ can be completely disassembled for easy cleaning.



Benefits

Approved for use with a wide range of our Dräger ventilators

Each patient needs individualised care in breathing therapy. So as to support you in giving the precise and necessary care for your patient, we offer you quality ventilation equipment that specifically tailors to your patient's individual needs. That is why our Dräger SelfTestLung is designed and tested for use with a wide range of our Dräger ventilators for any clinical setting.

Complete disassembly for easy cleaning

To reduce the risk of infections and adhere to those protocols, we have engineered a test lung that can be safely reused. With proper reprocessing, according to the instructions for use, you can safely reduce the spread of viruses and bacteria. To make cleaning easier, our Dräger SelfTestLung can be completely disassembled.

Uncompromising quality with responsibility in mind

Control quality, control emissions: As an international medical technology manufacturer, we stand for high quality. True to our commitment to high quality, our test lung is designed in Germany and manufactured in Europe. This enables us to lower emissions during the transport to our distribution warehouse.

Technical Data

Ambient conditions

In storage:

| | |
|-------------|-----------------------------|
| Temperature | -20 to +60°C (-4 to +140°F) |
|-------------|-----------------------------|

In operation:

| | |
|--------------------------------------|--|
| Temperature | +10 to +40°C (+50 to +104°F) |
| Air pressure | 500 to 1,100 hPa (7.25 to 15.95 psi) |
| Relative humidity | 10 to 95% (non condensing) |
| Max. volume | 1,000 mL |
| Material | Breathing bag: Silicone (VMQ) Connector: Polyamide (PA) Housing plate: Polyketone (PK) |
| Classification Medical Device Europe | Class I |

Ordering Information

| | |
|---|-----------|
| Test lung SelfTestLung™, latex-free, reusable | MP 20 020 |
|---|-----------|

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION DACH

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06460-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel: +65 6872 9288
Fax: +65 6259 0398
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969, USA
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

