

Dräger SSK 5000 Drug Testing Device

Ideal for surface drug testing applications, the Dräger SSK 5000 sampling system is quick and easy to use. When used with the Dräger DrugTest 5000, the sample can be transferred to a test cassette and screened for possible drug contamination. The sample can also be stored and transported for further testing. FOR LAW ENFORCEMENT USE ONLY.



Benefits

Surface Testing

The Dräger SSK5000 consists of a vial containing a buffer solution and an integrated swab in the cap for sampling. The swab is used to wipe surfaces and test for drugs via the Dräger DrugTest 5000 System.

Unknown Substances

Use the swab of the SSK5000 to collect small amounts of, for example, an unknown granular formulation, transfer this into the buffer solution and test this for drugs via the Dräger DrugTest 5000 System.

Confirmation Testing

The substance of surfaces and unknown formulation which are transferred into the buffer solution using the SSK5000 can be stored and transported safely in the sample container. For dosing of liquid, the sample container can be modified into a dropping bottle thanks to a special screw cap. The liquid to be tested can then be separated or reduced for subsequent analysis using droplet dosing.

Services



Dräger Service

When you rely on Dräger breath alcohol and drug screening equipment, you can rest assured you made the right choice. For more than 60 years, we've been the global market leader in advanced breath alcohol testing and employ cutting-edge technology to provide devices that are fast, reliable, and easy to use. That's the Dräger Service Advantage.

Related Products

ST-430-2006



Dräger Alcotest® 9510

The Dräger Alcotest® 9510 is an advanced breath alcohol measuring instrument for evidential applications. It is designed to comply with national and international requirements and regulations (OIML R 126). With dual sensor technology, intuitive color touch screen interface, and modern design the Alcotest® 9510 sets superior standards. FOR LAW ENFORCEMENT USE ONLY.

D-74695-2013



Dräger Alcotest® 6820

Dräger Alcotest® 6820 is the logical successor to the successful Alcotest® 6810. The compact hand-held measuring instrument Alcotest® is quickly ready for use and is adaptable to international guidelines. With an even more sturdy enclosure and the proven Dräger sensor, Alcotest® 6820 is a reliable partner for breath alcohol analyses. FOR LAW ENFORCEMENT PURPOSES OR USE IN DOT-REGULATED INDUSTRIES ONLY.

ST-15093-2008



Dräger Alcotest® 7510

The Alcotest® 7510 allows for breath alcohol analysis for any application, and is listed on the National Highway Traffic Safety Administration's Conforming Product List (NHTSA CPL) as an evidential breath tester. This instrument has received the OIML Certificate of Conformity. FOR LAW ENFORCEMENT USE ONLY.

D-54720-2012



Dräger DrugTest 5000

The Dräger DrugTest 5000 system is a fast, accurate means of testing oral fluid samples for drugs of abuse, such as amphetamines, designer amphetamines, opiates, cocaine and metabolites, benzodiazepines, cannabinoids, and methadone. The analyzer offers easy data management with the Dräger Diagnostics Software. FOR LAW ENFORCEMENT USE ONLY.

Related Products



ST-189-2005

Dräger Interlock® XT

The Dräger Interlock® XT is a breath-alcohol measuring instrument that can prevent a vehicle from being started if the driver's breath alcohol concentration is higher than his or her state's legal limit. It meets worldwide device requirements, offers the greatest possible convenience for the user, and sets a new standard for alcohol interlocks. FOR USE WITH AN IGNITION INTERLOCK RESTRICTION ONLY.

Ordering Information

Item	Order No.
Draeger DrugTest 5000 Analyzer	83 19 900
Draeger DrugTest 5000 Test Kit 6-drug panel (20 per box)	83 19 965
Accessories 2	
Draeger SSK 5000 Device for surface samples (20 per box)	83 20 490

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

USA

Draeger, Inc.
7256 S. Sam Houston Parkway W.,
Suite 100
Houston, TX 77085
1 800 4DRAGER
(1 800 437 2437)

Locate your Regional
Sales Representative at:
www.draeger.com/contact

