

Dräger SSK 5000 Drug Testing Device

Ideal for surface drug testing applications, the Dräger SSK 5000 sampling system is quick and easy to use. When used with the Dräger DrugTest 5000, the sample can be transferred to a test cassette and screened for possible drug contamination. The sample can also be stored and transported for further testing.



Benefits

Surface Testing

The Dräger SSK 5000 consists of a vial containing a buffer solution and an integrated swab in the cap for sampling. The swab is used to wipe surfaces and test for drugs via the Dräger DrugTest 5000 System.

Confirmation Testing

The substance of surfaces and unknown formulation which are transferred into the buffer solution using the SSK 5000 can be stored and transported safely in the sample container. For dosing of liquid, the sample container can be modified into a dropping bottle thanks to a special screw cap. The liquid to be tested can then be separated or reduced for subsequent analysis using droplet dosing.

Unknown Substances

Use the swab of the SSK 5000 to collect small amounts of, for example, an unknown granular formulation, transfer this into the buffer solution and test this for drugs via the Dräger DrugTest 5000 System.

Related Products



Dräger DrugTest 5000

No pipetting, no drips, no timing: With the Dräger DrugTest 5000, a drug test is carried out simply and quickly. The collected saliva sample can be analyzed immediately – for accurate results on the spot.

Related Products



Dräger DCD 5000

Simple, fast and secure detection of intoxicants and drugs:
The Dräger DCD 5000 does not only enable an easy collection of saliva (oral fluid), you can also store the samples and send them to the laboratory for analysis – without the sample being falsified along the way.



Dräger DrugCheck® 1200 STK

The reliable method to prevent drug abuse in a cost-effective way.
Current drug influence can quickly and easily be detected with the saliva test.

Ordering Information

Dräger SSK 5000

Dräger SSK 5000, 20 Pieces	83 20 490
Pack sizes: 20 Pieces buffer-filled collection devices incl. swab-tap	

Dräger DrugTest 5000

Dräger DrugTest 5000 Analyzer, 1 Piece	83 19 900
Dräger DrugTest 5000 Test-Kits, 20 Pieces	83 19 970

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

REGION EUROPE

Dräger Safety AG & Co. KGaA
Revalstraße 1
23560 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Safety AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel: +65 6872 9288
Fax: +65 6259 0398
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06406-100 - Barueri - SP
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

