


Anesthetic Agent Consumption

Total	Isofluran	Desfluran	Sevofluran
7,836 EUR	n/a EUR	2,047 EUR	5,789 EUR
	Consump. n/a ml	Consump. 4,651 ml	Consump. 9,649 ml
	Price/btl. €50.00	Price/btl. €110.00	Price/btl. €150.00



Efficiency

Consumption vs. Uptake



Avg. FG-Flow



Costs

Cost / Case

11.36€

↓1.78%

Cost / Minute

0.13€

↑1.42%

Applied FG-Flows

- Minimal Flow 0 - 0.5 l/min
- Low Flow 0.5 - 1 l/min
- Medium Flow 1 - 8 l/min
- High Flow 8 - 20 l/min

