

Declaration of Compatibility <i>Kompatibilitätserklärung</i>	Date / Datum 2021-06-30
European Directive 93/42/EEC, Annex I, 9.1 <i>Europäische Richtlinie 93/42/EWG, Annex I, 9.1</i>	Document ID / Dokument Nr. HD108-045-2106-008-0

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck
Germany

EC Certificate: G10105780037Rev.01
Valid until: 2024-05-26

hereby declares under its sole responsibility, that the CE marked products
erklärt hiermit in alleiniger Verantwortung, dass die CE- gekennzeichneten Produkte

Product name / Produktbezeichnung	Manufacturer / Hersteller
VentStar Coax	Drägerwerk AG & Co. KGaA
Anaesthesia Set Coax	Drägerwerk AG & Co. KGaA
Venstar Coax GSL	Drägerwerk AG & Co. KGaA

and

Primus, Primus IE, Apollo, Perseus A500, Zeus, Zeus IE, Fabius	Drägerwerk AG & Co. KGaA
Evita 2dura / 4 / XL, Savina 300, Savina Color, Evita V300, Infinity Acute Care System Workstation Critical Care (Evita V500)	Drägerwerk AG & Co. KGaA

may be operated in combination as a system. Combination of the listed products has been tested and assessed. The combination is done according to the instructions for use of the respective devices in therein indicated configurations only.

Any deviation of the indicated configurations not authorized by Drägerwerk AG & Co. KGaA will invalidate this declaration.

in Kombination als ein System betrieben werden dürfen. Der Anschluss der oben genannten Geräte wurde getestet und bewertet. Die Kombination erfolgt gemäß der Gebrauchsanweisungen der entsprechenden Geräte in den darin kenntlich gemachten Konfigurationen.

Jede nicht durch Drägerwerk AG & Co. KGaA autorisierte Abweichung von den aufgeführten Konfigurationen führt zur Ungültigkeit dieser Erklärung.

Director
Quality & Regulatory
Business Unit Hospital Consumables & Accessories



Timo Harms



Director
Research & Development
Business Unit Hospital Consumables & Accessories



Kalle Heckmann

Attachment: Schedule 1 to HD108-045-2106-008-0

Schedule 1 to HD108-045-2106-008-0

		Primus IE *	Primus	Apollo	Perseus A500	Zeus	Zeus IE	Fabius**	Evita XL	Evita 4	Evita 2 dura	Savina 300 (≤ SW3.53)	Savina 300 (> SW3.53)	Savina	Evita V300	Infinity Acute Care System Workstation Critical Care	Evita V600/V800
MP00373	VentStar Coax Watertrap 210	x	x	x	x	x	x					x					
MP00378	Infinity ID Breathing Circuit Coax Watertrap 210	x	x	x	x	x	x					x					
MP00379	VentStar Coax (P) 150	x	x	x	x	x	x	x	x	x	x	x					
MP03373	VentStar Coax	x	x	x					x	x	x	x	x	x			x
MP03374	VentStar Coax 180	x	x	x					x	x	x	x	x	x	x	x	x
MP03375	ID Coax 180	x	x	x					x	x	x	x	x	x	x	x	
MP03376	VentStar Coax 230	x	x	x					x	x	x		x		x	x	x
MP03377	ID Coax 230	x	x	x					x	x	x		x		x	x	
MP03378	Ventstar Coax GSL 230		x		x			x									
MP03379	Ventstar Coax WT GSL 210		x		x			x									
MP00358	Anesthesia Set Coax GSL 230		x		x			x									
MP00359	Anesthesia Set Coax WT GSL 210		x		x			x									
MP00357	Ventstar ID Coax GSL 230		x		x			x									
MP00367	Ventstar ID Coax WT GSL 210		x		x			x									
MP00368	Anesthesia ID Set Coax GSL 230		x		x			x									
MP00370	Anesthesia ID Set Coax WT GSL 210		x		x			x									

* Device supports RFID Functionality

** Fabius Family includes Fabius Tiro/ Fabius GS/ Fabius GS Premium/ Fabius MRI/ Fabius OS/ Fabius Plus

Drägerwerk AG & Co. KGaA
 Moisinger Allee 53-55
 23558 Lübeck, Germany
 Postal address:
 23542 Lübeck, Germany
 Tel +49 451 882-0
 Fax +49 451 882-2080
 info@draeger.com
 www.draeger.com
 VAT no. DE135082211

Bank details:
 Commerzbank AG, Lübeck
 IBAN: DE95 2304 0022 0014 6795 00
 Swift-Code: COBA DE FF 230
 Sparkasse zu Lübeck
 IBAN: DE15 2305 0101 0001 0711 17
 Swift-Code: NOLADE21SPL

Registered office: Lübeck
 Commercial register:
 Local court Lübeck HRB 7903 HL
 General partner: Drägerwerk Verwaltungs AG
 Registered office: Lübeck
 Commercial register:
 Local court Lübeck HRB 7395 HL

Chairman of the Supervisory Board
 for Drägerwerk AG & Co. KGaA
 and Drägerwerk Verwaltungs AG:
 Prof. Dr. Nikolaus Schweickart
 Executive Board:
 Stefan Dräger (chairman)
 Rainer Klug
 Gert-Hartwig Lescow
 Dr. Reiner Piske
 Anton Schrofner