



Dräger Jaundice Meter JM-105

SKILLS CHECKLIST

Date:

Participant Name:

Hospital:



PASS	NEEDS REVIEW	SKILL	REVIEWER
		GENERAL	
		Locate the Jaundice Meter JM-105 within the unit.	
		Perform the Light Checker procedure.	
		State whether the Checker range is acceptable for use.	
		Perform the proper cleaning procedure prior to use.	
		Locate the ON/OFF switch and turn the unit ON.	
		Note the unit of measure: mg/dL or $\mu\text{mol/L}$.	
		Enter the patient's information.	
		Demonstrate how to take a reading.	
		State where on the infant's body the reading will be taken.	
		Flag a patient with the priority flag.	
		Look up a patient file in the HISTORY .	
		Turn the unit OFF.	
		Explain what the device displays when the measurement is out of range.	
		Replace the Jaundice Meter JM-105 on its charging cradle and explain what the orange light means.	

This Skills Checklist is provided for information purposes only and focuses on the Dräger Medical product referenced. The Skills Checklist assists in training *how to use or operate* this product. It does not address issues of medical practice or clinical judgment related to use of the product.

The checklist is intended only as a suggestion for your consideration; it may or may not be consistent with your institutional expectations for staff competency. It is the responsibility of each institution or provider to make their own determinations regarding matters of training and clinical protocol for their employees and the products that they use.