


BabyFlow® – Nasal CPAP Consumables and Accessories

Patient comfort, highly effective therapy and cost-efficiency – these are just some of the requirements fulfilled by BabyFlow®.



Benefits

Babyflow – nasal cpap

The Dräger nCPAP BabyFlow interface is at the cutting edge of medicine. It not only allows therapy that is beneficial to the patient but, as a non-invasive mode of ventilation, also protects against complications that may arise from intubation. The BabyFlow can also save money, as the high investment costs for additional ventilation devices do not apply when the BabyFlow is used on existing Babylog or Evita ventilators*.

* Prerequisite here is a software that conforms with the operation of the BabyFlow, as well as the NIV/NeoFlow option for the Evita range.

Prongs, masks and bonnets

The carefully produced silicone prongs and masks are available in different sizes. Their shapes are optimised anatomically, and ensure both safe use and an outstanding fit, despite specific facial features. The practical measuring tape makes it simple to determine the correct size. The bonnets are all colour-coded to optimize processes and simplify handling. They are made of high-quality fabric and can be easily and reliably fastened, while still allowing fast access to the head.

Gentle care

The BabyFlow not only supports the necessary flexibility in everyday clinic operations - its technical solution is also gentle on the patient. For example, the selected CPAP pressure level is as independent of the flow as possible, thus resulting in quiet respiratory support. This means the BabyFlow is fully in accord with the philosophy of gentle care, which benefits the development of the infant.

Ordering Information

OVERVIEW OF SIZES FOR DISPOSABLE BONNETS

Size	Colour	Head circumference (max.)		Weight		Order number
		cm	inches	grams	Lbs	
S	Yellow	23	9.1	700 – 1,000	1.5 – 2.2	8418534
M	Light orange	26	10.2	1,000 – 1,300	2.2 – 2.9	8418535
L	Light salmon colour	28	11.0	1,300 – 1,600	2.9 – 3.5	8418536
XL	Lavender	30	11.9	1,600 – 2,400	3.5 – 5.3	8418537
XXL	Blue	33	13.0	2,400 – 3,000	5.3 – 6.6	8418538
XXL plus	Green	36	14.2	3,000 >	6.6 >	8418539

Headgear

XXS 5 pieces	MP03880
XS 5 pieces	MP03881
S 5 pieces	MP03882
M 5 pieces	MP03883
L 5 pieces	MP03884
XL 5 pieces	MP03885

BabyFlow, disposable (pack of 20 pcs.)

8418583

Prongs, XS (pack of 10 pcs.)	8418415
Prongs, S (pack of 10 pcs.)	8418605
Prongs, M (pack of 10 pcs.)	8418416
Prongs, L (pack of 10 pcs.)	8418531
Prongs, XL (pack of 10 pcs.)	8418417
NeoMask, S (pack of 10 pcs.)	8418491
NeoMask, M (pack of 10 pcs.)	8418490
NeoMask, L (pack of 10 pcs.)	8418619

Standard prongs

Prongs XXL 10 pieces	84 18 617
----------------------	-----------

Wide prongs

Prongs 3, 10 pieces	84 18 603
Prongs 4, 10 pieces	84 18 607
Prongs 5, 10 pieces	84 18 609
Prongs 6, 10 pieces	84 18 611
Soft connector N, 20 pieces	MP03826

Demo package (disposable)

8418550

Consisting of:	
BabyFlow, disposable (5 pcs.)	
1 x Prong, medium (pack of 10 pcs.)	
1 x Prong, medium plus (pack of 10 pcs.)	
1 x NeoMask, M (pack of 10 pcs.)	
1 x Cap, L (pack of 5 pcs.)	
1 x Cap, XL (pack of 5 pcs.)	

Each component of the demo package is also available separately

Notes

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck, Germany

REGION DACH

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama S. de R.L.
59 East Street, Nuevo Paitilla
House 30, San Francisco Town
Panama City, Panama
Tel +507 377 9100
Fax +507 377 9130
servicioalcliente@draeger.com

REGION ASIA PACIFIC

Draeger Singapore Pte. Ltd.
25 International Business Park
#04-20/21 German Centre
Singapore 609916
Tel +65 6308 9400
Fax +65 6308 9401
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

