

Dräger BabyFlow® plus Consumables and Accessories

The superior non-invasive respiratory support system, designed to maximize performance and patient comfort.



D-488842-2015

Benefits

Performance

The BabyFlow plus respiratory support system was designed based on the results of state-of-the-art technology (Computational Flow Analysis) to minimize the imposed work of breathing for the patient. The streamlined adapter head optimizes the flow distribution.

Enhanced comfort for patient and caregiver

We want your patients to feel comfortable. The less discomfort a patient experiences, the greater the opportunity for healthy growth and development during this critical time in the patient's recovery.

The adjustable angle allows you to fit the prongs and masks properly on the patient, thus reducing leaks and providing more stable CPAP levels. More stable CPAP levels reduce noise levels, making BabyFlow plus a low noise level respiratory support system, therefore, supporting Developmental Care practices and offering the best possible therapy for tiny patients.*

Ease of use

The ease of use of BabyFlow plus in combination with innovative components offer effective and gentle therapy where it's needed.

The BabyFlow plus respiratory support system is designed to address the current issues with non-invasive respiratory support therapy. Every single component, including individually sized masks and nasal prongs, has been engineered to maximize performance and patient comfort.

The flexible, kink resistant hoses enable optimal positioning of the patient. Swivel connectors reduce the necessary disconnections of the breathing circuit to a minimum, thus ensuring stable CPAP pressures.

In combination with the Dräger Babylog® VN500, the BabyFlow plus supports a comprehensive selection of modalities for non-invasive respiratory support. Babylog VN500 automatically adjusts the flow according to the existing leaks to ensure stable CPAP pressure and prevent nuisance alarms.

* L. Kirchner, M. Wald, V. Jeitler, A. Pollack, „In vitro Comparison of Noise Levels produced by different CPAP generators“, in: Neonatology, 2012; 101: 95-100

* Karam et al., „Noise levels during nCPAP are flow-dependent but not device dependent“, in: Archives of Disease in Childhood Neonatal Edition, 2008: 93: F 132-F134

Accessories

D-33872-2015-D-37401-2015



Prongs

The nasal prongs are available in 11 sizes and two different widths. In combination with the BabyFlow plus, the prongs deliver an effective seal and reduce pressure on the infant's septum, therefore, minimizing skin breakdown.

D-23992-2015-D-34333-2009



Masks

The nasal masks are available in six different sizes. Sizes S, M and L are well-known and accepted silicone masks for smaller patients up to approximately 2 kg (4.41 lbs). The new masks for patients from 2 kg (4.41 lbs) up to 5 kg (11 lbs) are designed for maximum therapeutic effectiveness and patient comfort using 3D infant facial scan technology. The flexible bellows allow the mask to find a natural position over the infant's nose, which helps minimize leaks and maintain a consistent nCPAP level.

D-23991-2015



Headgear

The headgear helps secure the BabyFlow plus quickly and easily while minimizing pressure points. It can be gently wrapped around the infant's head and secured with adjustable closures. It is available in six sizes to fit head circumferences from 17 cm to 42 cm (6.69 inches to 16.54 inches).

D-48861-2015



Microfiber Caps

The flexible and breathable microfiber caps are available in seven sizes. The disposable caps include side straps and top straps to help secure and stabilize the BabyFlow plus.

Ordering Information

BabyFlow® plus disposable

S	20 pieces	MP03700
M	20 pieces	MP03701
L	20 pieces	MP03702

Microfiber Caps

XS	white	5 pieces	84 18 533
S	yellow seam	5 pieces	84 18 534
M	red seam	5 pieces	84 18 535
L	green seam	5 pieces	84 18 536
XL	orange seam	5 pieces	84 18 537
XXL	blue seam	5 pieces	84 18 538
XXL plus	white/black seam	5 pieces	84 18 539

Headgear

XXS	5 pieces	MP03880
XS	5 pieces	MP03881
S	5 pieces	MP03882
M	5 pieces	MP03883
L	5 pieces	MP03884
XL	5 pieces	MP03885

Neo mask

S	10 pieces	84 18 491
M	10 pieces	84 18 490
L	10 pieces	84 18 619
4	10 pieces	MP01494
5	10 pieces	MP01495
6	10 pieces	MP01496

Standard prongs

Prongs XS	10 pieces	84 18 415
Prongs S	10 pieces	84 18 605
Prongs M	10 pieces	84 18 416
Prongs L	10 pieces	84 18 531
Prongs XL	10 pieces	84 18 417
Prongs XXL	10 pieces	84 18 617

Wide prongs

Prongs 3	10 pieces	84 18 603
Prongs 4	10 pieces	84 18 607
Prongs 5	10 pieces	84 18 609
Prongs 6	10 pieces	84 18 611
Soft connector N	20 pieces	MP03826

Demo pack BabyFlow® plus, size S

MP03703

BabyFlow® plus, size S, 10 pieces

2 x Prongs XS

2 x Prongs S

2 x Prongs M

2 x Prongs 3

2 x Neo mask S

2 x Neo mask M

2 x Neo mask L

2 x Headgear XXS
 2 x Headgear XS
 2 x Headgear S
 Tape measure and sizing gauge
 Soft connector N, 20 pieces

Demo pack BabyFlow® plus, size M

MP03704

BabyFlow® plus, size M, 10 pieces

2 x Prongs M
 2 x Prongs L
 2 x Prongs XL
 2 x Prongs 4
 2 x Neo mask M
 2 x Neo mask L
 2 x Neo mask 4
 2 x Headgear S
 2 x Headgear M
 2 x Headgear L

Tape measure and sizing gauge
 Soft connector N, 20 pieces

Demo pack BabyFlow® plus, size L

MP03705

BabyFlow® plus, size L, 10 pieces

2 x Prongs XL
 2 x Prongs XXL
 2 x Prongs 5
 2 x Prongs 6
 2 x Neo mask 4
 2 x Neo mask 5
 2 x Neo mask 6
 2 x Headgear M
 2 x Headgear L
 2 x Headgear XL

Tape measure and sizing gauge
 Soft connector N, 20 pieces

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany

Locate your Regional Sales
Representative at:
www.draeger.com/contact



REGION EUROPE CENTRAL AND EUROPE NORTH

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com